

## Platinum Care Foam Mattress – 519 Series

- This premium support surface offers outstanding comfort, pressure reduction and pressure redistribution
- The heel section uses smaller cells to address the sensitive pressure reduction needs of the heels and ankles
- Aids in the prevention and treatment of stage I pressure ulcers
- Zippered, breathable top cover moves with the patient to reduce shearing
- Soft and durable top covers are bacteria resistant, flame retardant, fluid resistant, anti-fungal, anti-static, non-allergenic, have near 0% vapor transmission, and are easy to clean
- Soft and durable bottom covers are bacteria resistant, flame retardant, fluid resistant, anti-static, non-allergenic, stain resistant, and easy to clean
- Meets 16 CFR 1633 requirements
- Designed to meet the requirements of HCPCS code: E0184
- Maximum Weight Capacity 425 lb - EVENLY DISTRIBUTED
- Five year limited warranty
- Latex Free



	w/ zipper	w/ zipper + DPM*
75" x 35" x 6"	51975-1633N	51975P-1633N
80" x 35" x 6"	51980-1633N	51980P-1633N
80" x 42" x 6"	5198042-1633N	5198042P-1633N
84" x 35" x 6"	51984-1633N	51984P-1633N
84" x 42" x 6"	5198442-1633N	5198442P-1633N

w/ heel slope		
	w/ zipper	w/ zipper + DPM*
75" x 35" x 6"	51975H-1633N	51975PH-1633N
80" x 35" x 6"	51980H-1633N	51980PH-1633N
80" x 42" x 6"	5198042H-1633N	5198042PH-1633N
84" x 35" x 6"	51984H-1633N	51984PH-1633N
84" x 42" x 6"	5198442H-1633N	5198442PH-1633N

\*defined perimeter



### Specifications

Weight .....	24.5 lb
Warranty .....	5 years
Maximum Weight Capacity .....	425 lb
Foam Core.....	Antimicrobial, 100 % High Density Polyurethane
Cover .....	Nylon
Pressure Reduction .....	Best

