


You May Also Need



Invacare® Heavy-Duty Bed Rail
Model# BAR6640IVC



Invacare® Reliant™ 600 Heavy-Duty Lift
Model# RPL600-1



Invacare® Softform Premier™
Mattress
Model# IPM1080B42

The **Invacare® BAR600IVC Bariatric Bed** is a heavy-duty full-electric bed frame designed for bariatric individuals. It is capable of supporting patients who weigh up to 600 lb. The bed frame extends to the edges of the bed deck for better support than most other split-spring bariatric beds. Using the hand pendant, the patient or caregiver is able to change the positioning of the head and foot sections, as well as the bed height. Despite being substantial in size, bed setup within the home is easy – no special tools are needed.

Features

Heavy-duty frame design ensures added strength for client support

Heavy-duty bed sleep surface is 17% larger than conventional beds to ensure comfortable positioning

Split-spring design with universal bed ends for ease in setup

No special tools necessary for in-home setup

Easy-release motor mounts allow quick motor removal

Robotically welded, steel frame construction

Pendant allows for dynamic client positioning

Quiet, smooth operation

Use with Invacare bariatric mattress (BARMATT42) and bed rails (BAR6640IVC):

- BARMATT42 meets fire safety regulations 16CFR1633 and 16CFR1632

600 lb. client weight capacity

W
S
T
Q
D

Invacare® BAR600IVC Model # BAR600IVC




Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
BAR600IVC	88" L 14"-22" H 42" W	80" L 42" W	600 lb.	N/A	N/A	UL962	5 Years Welds/Frame/ Structure 1 Years Mechanical/ Electrical
BAR5000IVC Headspring				114 lb.	43" L 10" W 43" D		
BAR5490IVC Footspring				97 lb.	42.5" L 9.5" W 43" D		
BAR5301LOW Bed End				96 lb.	45" L 7" W 22.75" D		

Invacare® Bariatric Foam Mattress and Invacare® Reduced Gap Heavy-Duty Half-Length Bed Rail

Model # BARMATT(39/42/48/EXP)

Model # BAR6640IVC



Dimensions
LxWxD



Gauge



Regulation
Compliance



Limited
Warranty

Model #	Dimensions LxWxD	Gauge	Regulation Compliance	Limited Warranty
BARMATT39	80"x39"x6.5"	N/A	16CFR1632/33	1 Year
BARMATT42	80"x42"x6.5"	N/A	16CFR1632/33	1 Year
BARMATT48	80"x48"x6.5"	N/A	16CFR1632/33	1 Year
BARMATTEXP	80"x39"-48"x6.5"	N/A	16CFR1632/33	1 Year



Rail Length*



Rail Height
Above Deck
Raised



Rail Height
Above Deck
Lowered



Width
Between
Rails



Shipping
Weight



Limited
Warranty

Model #	Rail Length*	Rail Height Above Deck Raised	Rail Height Above Deck Lowered	Width Between Rails	Shipping Weight	Limited Warranty
BAR6640IVC	28.75"	15"	3.5"	42"	25 lb.	1 Year

W
S
T
Q
D

Invacare® Bariatric Foam Mattress

Durable 11-ounce vinyl cover is anti-bacterial and waterproof for easy use and care (BARMATT42 only)

Durable nylon cover (BARMATT39/48/EXP only)

Solid inner core

Mattress is reversible and flippable

Constructed of latex-free, high resilience foam for extended life and service for the homecare provider

Invacare® Reduced Gap Heavy-Duty Half-Length Bed Rail

Heavy-duty 1" chrome-plated, welded steel construction
Bolt-on hardware for added strength

For optimal safety and convenience, four half-length rails are recommended (one set for head area and one set for foot area)

Use with BAR600IVC bed only

You May Also Need



Invacare® Reliant™ 600 Heavy-Duty Lift
Model# RPL600-1



Invacare® Bariatric Floor Stand with Trapeze
Model# BARTRAP



Invacare® BAR600IVC Bariatric Bed
Model# BAR600IVC

*Measurement on flat bed with crossbraces installed as per instruction sheet.

You May Also Need



Invacare® Heavy-Duty Paddle Walker
Model# 6291-HDA



Invacare® Reliant™ 600 Heavy-Duty Lift
Model# RPL600-1



Invacare® BAR750 Bariatric Bed
Model# BAR750

The **Invacare® Trapeze** unit is an important patient room accessory designed to help the users change positions while in bed and to aid in transferring from bed to chair with easy attendant assistance. Trapeze units have a wide range of height adjustments and handbar positions to maximize patient accessibility.

W
S
T
Q
D

Features

Three locking overhead boom angle positions for comfortable, secure use

Tool-free design simplifies assembly

Four piece sub-assembly significantly reduces setup time

Lightweight design for easy delivery

Free-standing structure permits use with Bariatric, Home Care or non-medical bed

Assisted weight capacity of 1000 lb.

Invacare® Bariatric Floor Stand with Trapeze

Model # **BARTRAP**



Model #	Height	Base	Reach	Diameter	Weight	Shipping Weight	Limited Warranty
BARTRAP	73"	50.5" L 40" W	35.5"	2.25"	72 lb.	80 lb.	3 Years