

Nasogastric Feeding Tubes — Pediatric



8888260802

Argyle™ PVC Feeding Tubes

- Thermosensitive material with Slide-Tex finish
- Colored radiopaque Sentinel Line™ tubing
- Two eyes
- Proximal end accepts male Luer connector
- Funnel connector
- Attached cap to cover proximal end
- Sterile

CODE	MARKINGS	FRENCH SIZE	LENGTH	SHIP CASE
8888261206	From 5cm to 25cm in 1cm increments	3.5 (1.2mm)	12" (31cm)	10
8888260802		5 (1.7mm)	16" (41cm)	50
8888261008		5 (1.7mm)	36" (91cm)	50
8888260711		6.5 (2.2mm)	16" (41cm)	50
8888260604		8 (2.7mm)	16" (41cm)	50
8888260406		8 (2.7mm)	42" (107cm)	50
8888260208		10 (3.3mm)	42" (107cm)	50

Curity™ PVC Feeding Tubes

155720*	8" (20cm)	5 (1.7mm)	15" (38cm)	50
155721*	8" (20cm)	5 (1.7mm)	36" (91cm)	50
155722*	8" (20cm)	8 (2.7mm)	15" (38cm)	50
155723*	8" (20cm)	8 (2.7mm)	42" (107cm)	50

* Does not contain radiopaque Sentinel Line or cm markings.