

Deluxe Preferred Care® Recliner Series-Wide (Model FR565WG)

- The Lumex FR565WG Patented recliner mechanism allows for almost infinite positioning of recliner
- **Proudly made in the USA** from domestic and globally sourced components!
- Unique Position Lock can stabilize recliner position for patient comfort
- A unique suspension system self-adjusts to each person, cradling the resident from head to heel and virtually eliminating the risk of “bottoming out”
- Extra-wide padded armrests and side panels maximize patient comfort
- Durable epoxy coated welded steel frame is easy to clean and maintain
- Quick-release upholstery allows for easy cleaning, removal and replacement
- Vinyl upholstered back and seat are sulfide stain, mildew and oil resistant
- A variety of accessories are available to customize the chair for specific individual needs
- Reduced armrest and seat height enables the recliner to fit under dining/ activity tables
- Four TENTE® swivel casters facilitate transport and maneuverability in confined areas
 - Three total locking casters, with easy to locate black locking tab, lock wheel and swivel with a single motion
 - Directional wheel lock, with high visibility green locking tab on right rear caster, enables the recliner to track straight while moving
- 41 Colors and Fabrics. Meets CA TB-133



Specs

Overall Height.....	48"	Seat Height.....	20"
Overall Width.....	32"	Seat Width.....	24"
Overall Depth.....	37"	Seat Depth.....	19.5"
Width Between Arms	24"	Shipping Weight	91 lb
Arm Height From Floor.....	28"	Maximum Weight Capacity	
Arm Height From Seat.....	8.5"	EVENLY DISTRIBUTED:	400 lb



Position Stabilizing System



Manufactured by:
 GF Health Products, Inc.
 Atlanta, GA 30360
Proudly made in the USA

www.grahamfield.com

© 2014 GF Health Products, Inc.
 Atlanta, Georgia 30360 GF0900047RevC14



Colors may vary slightly from actual vinyl

FR565WG401	Vintage Gold
FR565WG409	Warm Taupe
FR565WG427	Blue Ridge
FR565WG432	Imperial Blue
FR565WG454	Royal Blue
FR565WG6705	Dolce Sand
FR565WG6714	Dolce Burgundy
FR565WG6717	Dolce Moss
FR565WG6724	Dolce Jet
FR565WG6725	Dolce Indigo
FR565WG6727	Dolce Meteor
FR565WG8504	Sand
FR565WG8516	Chestnut
FR565WG8533	Port

FR565WG8561	Mulberry
FR565WG8567	Dove
FR565WG8584	Taupe
FR565WG8602	Oatmeal
FR565WG8604	Concrete
FR565WG8611	Terra Cotta
FR565WG8615	Berry
FR565WG8624	Granny Smith
FR565WG8626	Eucalyptus
FR565WG8628	Sea
FR565WG8801	Marine Blue
FR565WG8805	Jet
FR565WG8808	Champagne
FR565WG8809	Taupe

FR565WG8818	Sage
FR565WG9201	Cobblestone
FR565WG9202	Periwinkle
FR565WG9205	Gypsum
FR565WG9207	Ice Blue
FR565WG9208	Wineberry
FR565WG9210	Cypress
FR565WG9212	Blue Jade
FR565WG9214	Steel Blue
FR565WG9715	Eucalyptus
FR565WG9722	Desert Sage
FR565WG9764	Steel Grey



Manufactured by:
GF Health Products, Inc.
 Atlanta, GA 30360
Proudly made in the USA

www.grahamfield.com

© 2014 GF Health Products, Inc.
 Atlanta, Georgia 30360 GF0900047RevC14