




Texium™ IV Bag / Line Access Devices

REF 2300E-0500

Packaging	Quantity
Shelf pack	N/A
Case	100



BD SmartSite™ Bag Spike. Connection between IV bag and infusion set. Built-in BD SmartSite™ Needle-free Connector for closed transfer to IV bag. Does not contain natural rubber latex. BPA free.

REF 2309E


Packaging	Quantity
Shelf pack	N/A
Case	100



BD SmartSite™ Bag Spike with check valve. Connection between IV bag and infusion set. Built-in BD SmartSite™ Needle-free Connector for closed transfer to IV bag while mating with a Texium™ Closed Male Luer. Check valve allows no backflow of fluid into the IV bag. Does not contain natural rubber latex. BPA free.

REF 2000EA

Packaging	Quantity
Shelf pack	N/A
Case	100



BD SmartSite™ Needle-free Connector. Permits a closed injection and aspiration of fluids without the use of needles. Labeled for use with low-pressure power injectors up to 325 psi and maximum flow rate of 10 mL / second. Does not contain natural rubber latex. BPA free.

