

DonJoy® Armor™

The Armor FourcePoint™ is DonJoy's strongest prefabricated knee ligament brace. It is ideally suited for extreme and contact sports. The strong, aircraft-grade aluminum frame helps provide strength and confidence to aid for athletes who demand a knee brace that provides protection. Swooping thigh cuff helps allow inner thigh clearance and bilateral use. Indications include support for mild to severe ACL, PCL, MCL and LCL instabilities.

FEATURES & BENEFITS

- Clinically-proven 4-Points-of-Leverage System™
- Backed by DonJoy's Knee Protection Program, an industry first for prefabricated bracing (FourcePoint™ version only)
- Swooping Thigh Cuff for inner thigh clearance and bilateral brace use
- Strong aircraft grade aluminum frame (0.125" thick)
- 1 year warranty on frame and hinges

OPTIONS

- ACL, PCL, or Combined Instability (CI) strap configurations
- Standard or short frames available
- Logo included

ARMOR ACTION

PART NO.	LIGAMENT	HINGE	CALF LGTH.	LEG
11-1440-x	ACL	FourcePoint™	Short	Right
11-1441-x	ACL	FourcePoint™	Short	Left
11-1442-x	ACL	FourcePoint™	Standard	Right
11-1443-x	ACL	FourcePoint™	Standard	Left
11-1004-x-26000	ACL	Standard	Short	Right
11-1005-x-26000	ACL	Standard	Short	Left
11-1444-x	ACL	Standard	Standard	Right
11-1445-x	ACL	Standard	Standard	Left
11-1029-x-26000	CI	Standard	Short	Right
11-1030-x-26000	CI	Standard	Short	Left
11-1447-x	CI	Standard	Standard	Right
11-1448-x	CI	Standard	Standard	Left

ACCESSORIES

PART NO.	DESCRIPTION	SIZE
11-0075-x-00000	Lycra® Suspension Sleeve	S-XXL
11-0122-x	Neoprene Suspension Sleeve	XS-XXL
11-1006-x-06000	Impact Guard	S-XXL
11-1095-9-06000	Calf Pinch Guard	Universal
11-0956-0-06000	Condyle Pad	Universal
11-0398-9-00000	Full Extension Stop	Universal
25-0422-0-00000	Flexion / Extension Stop Kit	Universal
25-2428	Mesh Brace Bag	Universal
11-0016-x-06000	Sports Brace Cover (Standard)	S-XXL
11-1525-x	Tru-Pull® Advance Attachment	XS-XXXL
11-1415	Sili-Grip Strap Pad Kit	Universal
11-0543-x-06000	Neoprene Suspension Strap Kit	XS-XXL

Sizes available – page 97.

