

**MEDICARE
REIMBURSABLE!**

JOBST® FarrowWrap® 4000

The JOBST® FarrowWrap® 4000 is an easy-to-apply legpiece with just four bands, making this wrap a simple and effective alternative to bandaging and stockings. Contains Elastic Short-Stretch compression technology for comfortable resting and effective working compression. Also, includes a FarrowHybrid® sock which provides compression in the foot and ankle while acting as a non-compressive liner above the ankle.



Product Features:

- Easy to don with four bands
- Anatomically contoured to prevent gapping
- Inner sleeve assists with donning
- Cost-effective and includes six month warranty
- 30-40 mmHg* compression
- Ready-to-wear sizes XS-L, black or tan, and regular or tall
- Easy maintenance - machine wash and dry (no heat)


HCPCS Code: A6545 - Gradient Compression Wrap, Non-Elastic, Below Knee, 30-50 mmHg*, Each

Products are covered with the diagnosis of an open venous leg ulceration. The modifier AW must be used in conjunction with A6545 to signify the presence of an open wound. The A6545-AW code does not cover other conditions, including venous insufficiency without an ulcer present or lymphedema without an ulcer present. LT (left) and RT (right) must be used as additional modifiers to designate which limb will be treated with the compression garment. There is a separate modifier used if there are bilateral ulcers. The A6545 code will not reimburse both legs unless both legs have an active venous ulceration. Coverage is limited to 1 garment per leg every 6 months.

For additional information on coverage, coding and documentation requirements for billing A6545, suppliers should refer to the Surgical Dressings LCD and related Policy Article, which can be found on the PDAC website – www.dmepdac.com.

JOBST® FarrowWrap® 4000 Legpiece and FarrowHybrid® Liner Sizing

FarrowWrap® 4000 Leg	XS	S	M	L	Length (B-D)
Ankle (cB)	17-23 cm	19-28 cm	25-37 cm	29-43 cm	Reg 30-33 cm
Calf (cC)	25-38 cm	32-48 cm	37-56 cm	43-67 cm	Tall 34-37 cm



JOBST® FARROWWRAP® 4000

Legpiece	BSN Code	Size	Length	Color
	7666100	Extra Small	Regular	Tan
	7666101	Small	Regular	Tan
	7666102	Medium	Regular	Tan
	7666103	Large	Regular	Tan
	7666108	Extra Small	Regular	Black
	7666109	Small	Regular	Black
	7666110	Medium	Regular	Black
	7666111	Large	Regular	Black
	7666104	Extra Small	Tall	Tan
	7666105	Small	Tall	Tan
	7666106	Medium	Tall	Tan
	7666107	Large	Tall	Tan
	7666112	Extra Small	Tall	Black
	7666113	Small	Tall	Black
	7666114	Medium	Tall	Black
	7666115	Large	Tall	Black



FARROWWRAP® 4000

JOBST® FARROWHYBRID®

BSN Code	Size	Description
7666400	Hybrid sock CCL 1 size Small taupe	Mens shoe size 3-6 Womens shoe size 4-7
7666401	Hybrid sock CCL 1 size Medium taupe	Mens shoe size 7-9 Womens shoe size 8-10
7666402	Hybrid sock CCL 1 size Large taupe	Mens shoe size 10-12 Womens shoe size 11-13
7666403	Hybrid sock CCL 1 size Extra Large taupe	Mens shoe size 13-15 Womens shoe size 14-16
7666404	Hybrid sock CCL 2 size Small taupe	Mens shoe size 3-6 Womens shoe size 4-7
7666405	Hybrid sock CCL 2 size Medium taupe	Mens shoe size 7-9 Womens shoe size 8-10
7666406	Hybrid sock CCL 2 size Large taupe	Mens shoe size 10-12 Womens shoe size 11-13
7666407	Hybrid sock CCL 2 size Extra Large taupe	Mens shoe size 13-15 Womens shoe size 14-16
7666412	Hybrid sock CCL 1 size Small wide taupe	Mens shoe size 3-6 Womens shoe size 4-7
7666413	Hybrid sock CCL 1 size Medium wide taupe	Mens shoe size 7-9 wide Womens shoe size 8-10 wide
7666414	Hybrid sock CCL 1 size Large wide taupe	Mens shoe size 10-12 wide Womens shoe size 11-13 wide
7666415	Hybrid sock CCL 1 size Extra Large wide taupe	Mens shoe size 13-15 wide Womens shoe size 14-16 wide
7666408	Hybrid sock CCL 2 size Small wide taupe	Mens shoe size 3-6 wide Womens shoe size 4-7 wide
7666409	Hybrid sock CCL 2 size Medium wide taupe	Mens shoe size 7-9 wide Womens shoe size 8-10 wide
7666410	Hybrid sock CCL 2 size Large wide taupe	Mens shoe size 10-12 wide Womens shoe size 11-13 wide
7666411	Hybrid sock CCL 2 size Extra Large wide taupe	Mens shoe size 13-15 wide Womens shoe size 14-16 wide
7666416	Hybrid sock CCL 1 size Small wide black	Mens shoe size 3-6 wide Womens shoe size 4-7 wide
7666417	Hybrid sock CCL 1 size Medium wide black	Mens shoe size 7-9 wide Womens shoe size 8-10 wide
7666418	Hybrid sock CCL 1 size Large wide black	Mens shoe size 10-12 wide Womens shoe size 11-13 wide
7666419	Hybrid sock CCL 1 size Extra Large wide black	Mens shoe size 13-15 wide Womens shoe size 14-16 wide



FARROWHYBRID®

cB = circumference at the B
cC = circumference at the C

* The mean compression for an average size ankle.
Please refer to the product label and / or package insert for full instructions on the safe use of these products.

BSN medical Inc., an Essity company
5825 Carnegie Blvd. Charlotte, NC 28209-4633
Tel. (+1) 704 554 9933 Fax (+1) 800 835 4325
To order toll-free: BSN medical (+1) 800 552 1157
63150 RN ©2017 BSN medical Inc. REV 12/17

