



**ARMREST** 

**ADJUSTABLE** 

**FOOTREST** 

♦ ADJUSTABLE BASE PLATE ♦ FLIP UP ARMREST ♦ ALL STAINLESS STEEL CONSTRUCTION ♦ DURABLE POWDER COAT FINISH ♦ RECHARGEABLE BATTERY & WALL MOUNT CHARGER ◆ FIVE YEAR STRUCTURAL WARRANTY ♦ FIVE YEAR PRO-RATED ELECTRONICS WARRANTY

> CLICK HERE TO VIEW THIS PRODUCT: HTTPS://AQUACREEKPRODUCTS.COM/

> > **MATERIAL: 304 STAINLESS STEEL** FINISH: POWDER COAT - PER ORDER APPROX. WEIGHT: 117.38 LBS.

**Aqua Creek Products** 

**BASE PLATE** 

**DIMENSIONS ARE IN INCHES** 

**UNLESS OTHERWISE SHOWN** FRACTIONAL: ..... ±1/16

ANGULAR: ..... ±1° TWO PLACE DECIMAL: .... ±.01

THREE PLACE DECIMAL: .. ±.005

**TOLERANCES** 

**9889 GARRYMORE LANE** MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602 www.AquaCreek.com © 2020 AQUA CREEK PRODUCTS. ALL RIGHTS RESERVED. [45.72cm] 18"

# **PROPRIETY AND CONFIDENTIAL INFORMATION**

**RANGER 2 POOL LIFT** 

**CATALOG NO. F-RNGR2** 

**PART NO. INV01538** 

THIS DRAWING AND ALL INFORMATION CONTAINED HEREIN IS THE SOLE PROPERTY OF AQUA CREEK PRODUCTS LLC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION FROM AQUA CREEK IS PROHIBITED.

\*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

REVISIONS					
REV.	DESCRIPTION	DATE	APPROVED		
P1	INITIAL RELEASE	4/1/2021	J.MAXWELL		

P1 - Wednesday, August 11, 2021 11:17:28 AM

PRODUCT NAME: Ranger2™ Pool/Spa Access Lift CSI SECTION: 13 11 46 Swimming Pool Accessories

**APPLICATIONS:** 

• In Ground Pools & Spas with flush pool decks having up to 2-inch (5.08cm) high curbs.

### SUBMITTALS:

Alternate make or model lifts will not be considered unless equal to the specified product in every respect. Alternative products must be submitted for approval no less than ten (10) calendar days prior to the bid due date. Submittal data must include complete documentation relating to all the specified features and include all manufacturer's sales literature, specification sheets, installation/operation manuals, engineered drawings and documents that support the strength, service capacity, and the ability to meet the appropriate Federal Regulatory Standards associated with the product.

### **QUALITY ASSURANCE:**

- The product must be a standard product of a company registered with the United States Food and Drug Administration (FDA) as a class 1 medical device manufacturer and that is regularly engaged in the manufacture of commercial grade swimming pool and spa accessibility equipment.
- The product shall be structurally capable of supporting a 350 pound (158.76kg) load capacity and static load tested to 1.5 times the rated load.
- The product shall be certified to be compliant with UL test standard 60335-2-1000 Standard of Household and Similar Appliances: Particularly Requirements for Electrically Powered Pool Lifts.
- The product shall be tested and confirmed to comply with design standards outlined in the United States ADA Accessibility Standards for Accessible Designs.

#### SPECIFICATIONS:

- The products main structure shall be constructed of type 304 stainless steel; no dissimilar metallic components shall be allowed. The product shall be finished with a minimum 4-mil thick, UV resistant powder coat finish. All fastening hardware shall be of type 316 stainless steel. Product shall be supplied with a Linak Actuator suitable for medical use, a separate wall mount battery charger, a hand-held controller, and a 24-volt rechargeable battery.
- The chair of the product shall be of a rigid one-piece, polyethylene, molded design. The minimum seat width shall be 18-inches (45.72cm) and the seat back shall be a minimum height of 17.5-inches (44.45cm). The seat shall feature dual flip-up arms and supplied with a 2-inch (5.08cm) wide adjustable safety belt with hook and loop fasteners.
- The product shall be capable of accommodating a deck to water distance (water draft) of up to 14-inches (35.56cm) and capable of an anchor set back distance of up to 19-inches (48.26cm).
- Product can be field reversible and includes an adjustable and removable footrest that moves with the seat.

# WARRANTY:

Product shall have a minimum 5 year structural and 5 year pro-rated electronic warranty.

## **ANCHORING OPTIONS:**

INV01766

INV01910

F-RNGR2-W

F-RNGR2-C

Item No	Order No	Description
F-04CAJP	INV01411	Ranger 2 Lift Anchor for Standard Concrete Decks
F-08CAJP	INV01414	Ranger 2 Lift Anchor for Paver Applications
F-046QAB-AT1	INV01406	Ranger 2 Lift Anchor for Standard Concrete Decks with Quick Attach Knobs
F-40220-AT1	INV01444	Ranger 2 Lift Anchor for Wood Deck Applications

Ranger 2 Lift White Powder Coat w/White Seat

Ranger 2 Lift (Specify Seat and Powder Coat Colors)

## PRODUCT OPTIONS:

TROBUGI GI HORO.				
Order No	Description			
INV01990	Ranger 2 Seat Cover in Blue			
INV00525	Ranger 2 Lift Protective Lift Cover in Blue			
INV00528	Ranger 2 Lift Protective Lift Cover in Tan			
INV00527	Ranger 2 Lift Protective Lift Cover in Gray			
INV01509	Ranger 2 Lift Adjustable Height Seat Arm			
INV01693	Ranger 2 Lift Solar Charger			
INV01455	Ranger 2 Lift Adjustable Head Rest			
INV01457	Ranger 2 Lift Chest Strap			
INV01458	Ranger 2 Lift Transport Cart			
INV01950	Ranger 2 Lift Pull Out Leg Rest Blue			
INV01952	Ranger 2 Lift Pull Out Leg Rest Tan			
INV01954	Ranger 2 Lift Pull Out Leg Rest Gray			
INV01953	Ranger 2 Lift Pull Out Leg Rest White			
INV01379	Ranger 2 Lift Cycle Attachment			
INV01509	Ranger 2 Lift Adjustable Height Seat Arm			
	INV01990 INV00525 INV00528 INV00527 INV01509 INV01693 INV01455 INV01457 INV01458 INV01950 INV01952 INV01954 INV01953 INV01379			

Aqua Creek Products LLC 9889 Garrymore Lane Missoula, MT 59808

www.aquacreekproducts.com

Tel: 888-687-3552 / 406-549-0769

Fax: 406-541-8958



9889 GARRYMORE LANE MISSOULA, MT 59808 (406) 549-0769 FAX (406) 549-2602 www.AguaCreek.com © 2020 AQUA CREEK PRODUCTS. ALL RIGHTS RESERVED. **PROPRIETY AND CONFIDENTIAL INFORMATION** 

THIS DRAWING AND ALL INFORMATION CONTAINED HEREIN IS THE SOLE PROPERTY OF AQUA CREEK PRODUCTS LLC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION FROM AQUA CREEK IS PROHIBITED.

\*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

P1 - Wednesday, August 11, 2021 11:17:28 AM