



# Independent Total Body Trainer

## Independent Arms & Legs Removable 20" Swivel Seat Commercial Grade

PhysioMax	
Technical Details	
Swivel Seat	20" Oversized Seat Rotates 360° for Easy Access
Seat Position	Fully Adjustable Forward and Backward 14 Positions
Removable Seat	Seat Removes for Wheelchair Access (ADA)
Resistance System	Electromagnetic
Resistance Levels	1 to 16
Programs	18 Set Programs (See Computer Details)
Handles	Adjustable Comfort Grip Handles
Pedals	Othopedic Rehab Pedals with Heel Support
Stride	Independent Cycle Motion
Max User Weight	500 lbs.
Power	AC 110V-240V   50/60Hz
Product Type	Commercial Grade
Model #	TBT-1000
Computer Details	
Display	High Definition 9" LCD Display
Backlight	Easy to Read High Contrast Blue
USB Summary	Export Workout Data Including Peaks and Averages
Feedback	Time, Distance, Pulse, Steps Per Minute (SPM), Calories, WATTS, METs, & Resistance Level "1-16"
Programs	<b>18 Programs:</b> Manual, Warm Up, Interval, Ramps, Mountain, Rolling, Climb, Random, Hill, Fat Burn, Cardio, Strength, HRC 65%, HRC 85%, Constant SPM, Constant Watts, Goal Calories, Goal Distance
Heart Rate	Built in Receiver Compatible with Polar and other Manufactures *Chest Belt <b>Not</b> Included
Warranty	
Frame	Lifetime
Mechanical Parts	3 Years
Electronic Parts	3 Years
Labor	1 Year
Unit Dimensions	
Compact Dimensions	60" L X 25" W X 59" H
Unit Weight	350 lbs.



9" LCD Display



Adjustable Arms



Removable ADA Seat



Ergonomic Design



**HEALTHCARE**  
INTERNATIONAL

www.hcfitness.com

