Pink Nitrile® Exam Gloves

McKesson Pink Nitrile ® exam gloves feature a soft, stretchy material designed to make donning easy by quickly conforming to your hand and while also providing a wide range of motion. The gloves also feature improved tactile sensitivity, grip and fit. Fingertips are textured for decreased slippage during a procedure. Gloves are ambidextrous and non-sterile. They feature a beaded cuff and are available in four sizes.

Features and Benefits

- Powder-free
- · Beaded cuff
- Improved conformability provides a high quality fit and extended wear comfort
- Textured fingertips provide tactile sensitivity and dexterity
- · Non-sterile
- Not made with natural rubber latex



ORDERING INFORMATION

| DESCRIPTION | MFR # | QTY |
|---|-----------|-------------------|
| Pink Nitrile® Exam Gloves, Powder-Free, Textured Fingertips, Non-Sterile, Small | 14-6NPNK2 | 250ea/bx, 10bx/cs |
| Pink Nitrile® Exam Gloves, Powder-Free, Textured Fingertips, Non-Sterile, Medium | 14-6NPNK4 | 250ea/bx, 10bx/cs |
| Pink Nitrile® Exam Gloves, Powder-Free, Textured Fingertips, Non-Sterile, Large | 14-6NPNK6 | 250ea/bx, 10bx/cs |
| Pink Nitrile® Exam Gloves, Powder-Free, Textured Fingertips, Non-Sterile, X-Large | 14-6NPNK8 | 230ea/bx, 10bx/cs |

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on Care SM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2017 McKesson Medical-Surgical Inc.