

OxyHOME[®]

Stationary Concentrator

Keep Going.sm

NEW

BENEFITS

- Quiet Operation and Small Footprint
- Purity Check by Percentage
- Alarm Indicators: No Flow, Low/High Pressure, Purity, Power
- Twin Blower Cooling Design
- Contemporary Design
- 3-Year Replacement Warranty

PRODUCT SPECS **Part Number 1400-5100**

Flow Specifications	0 - 5 LPM
Oxygen Concentration	93% (+3%, -3%)
Oxygen Output Pressure	8 (+2, -2) psi
Rated Voltage/Frequency	120 VAC ±10% 60Hz
Power Consumption	350 Watts
Sound	≤ 45 dB(A)
Dimensions	13.4" W x 11.8" D x 25.6" H 34 cm x 30 cm x 65 cm
Weight	41.8 lbs (19 kg)



1400-5100

CONTACT US TODAY TO LEARN MORE

888-327-7301 • orders@oxygo.life • oxygo.life