

3M

Tegaderm™

Wound Care Solutions

3M™ Tegaderm™ Ag Mesh Dressing with Silver

The 3M™ Tegaderm™ Ag Mesh Dressing with Silver* provides a fast-acting, long-lasting antimicrobial barrier in an affordable, versatile and easy-to-use dressing.

*Caution: Federal law (U.S.A.) restricts the device to sale by or on the order of a licensed health care professional.



3M™ Tegaderm™ Ag Mesh Dressing with Silver*

Key benefits



Can be cut, folded or fluffed



May be used as primary dressing or with absorbent wound fillers



Effective against a wide range of microbes, including yeast, fungi** and antibiotic-resistant bacteria (MRSA and VRE)




Effective up to seven days



Can be moistened with wound exudate, sterile normal saline, sterile water or liquid hydrogel

Suggested applications

- ▶ Pressure, arterial, neuropathic (diabetic) and venous leg ulcers
- ▶ Surgical wounds
- ▶ Trauma wounds
- ▶ Superficial and partial thickness burns
- ▶ Abrasions
- ▶ Donor sites



Ordering information

Cat. No.	Size	Items/Box	Boxes/Case	HCPCS Code
90500	2 in x 2 in (5 cm x 5 cm)	5	4	A6402
90501	4 in x 5 in (10 cm x 12.7 cm)	5	4	A6403
90502	4 in x 8 in (10 cm x 20 cm)	5	4	A6403
90503	8 in x 8 in (20 cm x 20 cm)	3	4	A6404



For more information, contact your 3M Health Care Sales Representative, call the 3M Health Care Customer Helpline at **1-800-228-3957** or visit **3M.com/Medical**.



3M Medical Solutions Division
2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-228-3957
Web 3M.com/medical

*Caution: Federal law (U.S.A.) restricts the device to sale by or on the order of a licensed health care professional.
**Based upon *in vitro* studies. 3M data on file.

HCPCS codes have been provided to assist you in the preparation of insurance claims. Please note, however, that the reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coverage and coding for 3M products. Insurers' reimbursement policies can vary and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level. Health care providers should exercise independent clinical judgment in choosing the codes which most accurately describe the products provided.

3M and Tegaderm are registered trademarks of 3M. Please recycle.
Printed in USA. © 3M 2019. All rights reserved. 70-2011-7814-5