

BD® Blunt Fill and BD® Blunt Filter Needles

BD® Blunt Fill and BD® Blunt Filter Needles help reduce the risk of needlestick injuries during medication preparation due to a higher skin penetration force.*†

BD® Blunt Fill Needle

- Cost-effective, safe alternative for medication preparation that helps reduce needlestick risk.*
- Maintains easy vial access.*



Easily recognizable color



Beveled blunt needle
(enlarged to show detail)

BD® Blunt Filter Needle

- Includes a 5-micron filter to remove foreign matter, such as glass particles from ampoules.



Easily recognizable color



Beveled blunt needle
(enlarged to show detail)

- BD® Blunt Fill and BD® Blunt Filter Needles require at least 10x the force of an 18 G BD® Conventional Hypodermic Needle to penetrate simulated skin material but require no more than 2.5x the penetration force on vial stoppers.†
- Packaged in a red shield alerting users that the product is different from other needles.

The BD® Blunt Fill Needle and BD® Blunt Filter Needle are more preferred than competitive brands globally.†

Scan this code to watch a video on how to use the **BD® Blunt Fill** and **BD® Blunt Filter Needles**:



* As compared to a BD® 18 G Conventional Hypodermic Needle.

† Bench test results may not necessarily be indicative of clinical performance.

‡ Based on a 2025 global double-blinded survey conducted by BD among n=575 pharmacists (including pharmacy assistants and technicians) and nurses across the US, Canada, UK, France, Australia, Brazil, and India (no pharmacists were included in the India sample). All respondents primarily worked in hospital or acute care settings, had 3–35 years of experience, and regularly used hypodermic needles and syringes as part of their professional role. Each participant had experience with at least one BD needle and syringe and a competitive product within the past year. Product categories evaluated included conventional needles, blunt fill/filter needles, safety pivoting needles, safety shielding needles, and syringes.



BD® Blunt Fill and BD® Blunt Filter Needles

BD® Blunt Fill Needles

Product number	Description	Packaging
305180	18 G x 1 1/2 in.	100/box, 1000/case
305181	18 G x 1 in.	100/box, 1000/case

BD® Blunt Fill Needle (filter 5 micron)

Product number	Description	Packaging
305211	18 G x 1 1/2 in.	100/box, 1000/case

3 mL syringe with BD® Blunt Fill Needle

Product number	Description	Packaging
305060	BD Luer-Lok™ 3 mL with 18 G x 1 1/2 in.	100/box, 800/case

5 mL syringe with BD® Blunt Fill Needle

Product number	Description	Packaging
305062	BD Luer-Lok™ 5 mL with 18 G x 1 1/2 in.	100/box, 400/case

10 mL syringe with BD® Blunt Fill Needle

Product number	Description	Packaging
305064	BD Luer-Lok™ 10 mL with 18 G x 1 1/2 in.	100/box, 400/case

For more information on our products, please contact Medical Support Services at **800.555.7422** or medical.services@bd.com

To speak with your BD sales representative, call **888.237.2762**

[bd.com](https://www.bd.com)

