

TRANSITIONS®



Threshold Ramps, Mats, and Plates

EZ-ACCESS® knows that a gentle slope can make a world of difference. TRANSITIONS® threshold solutions allow for a smooth ground-to-sill transition over doorway thresholds throughout your home. Reduce preventable trip hazards and make your home more comfortable and safe for mobility users.



The TRANSITIONS Advantage

- Sturdy, seam-free surface makes TRANSITIONS ramps, mats, and plates the safe solution for wheelchair, scooter, walker, and cane users.
- Attractive designs are available in multiple styles and materials to suit your needs and taste.
- Slip-resistant tread and surfaces.
- Durably designed for years of continued use.
- Quick and easy setup and removal – installs in minutes, so there's little to no disruption for you and your family.
- Ideal for doorways, sliding glass doors, and raised landings.
- Portable, convenient, and durable.
- The entire line is backed by our lifetime warranty.



EZ-ACCESS®

MADE IN AMERICA

Visit ezaccess.com for product specifications

Safety starts at your front door.

TRANSITIONS Modular Entry Ramp



34" wide
700 lb. weight capacity

- Modular design with interchangeable extensions that can accommodate thresholds from 3/4" to 6-1/4" high.
- Installation is simple; add and remove extensions to accommodate your needs.
- Pivoting flap attaches to the top of the ramp to provide 1/4" of additional adjustability.
- Made of aluminum and features an extruded, slip-resistant surface for superior traction.

TRANSITIONS Angled Entry Ramp



36" wide
700 lb. weight capacity

- Swivel feet accommodate varying heights, angles, and surface types.
- Three lengths are available (12", 24", and 36") for thresholds ranging from 1-3/8" to 5-7/8" high.
- Manufactured using high-strength aluminum that will not rust or corrode.
- Features an extruded, slip-resistant surface for superior traction.

TRANSITIONS Modular Entry Mat



48" wide
850 lb. weight capacity

- Made of 100% recycled rubber.
- Provides a 2-1/2" height, or can be combined with optional risers to reach a maximum height of 4-3/4".
- Can be trimmed or notched for a custom fit.
- Ideal for both indoor and outdoor use.

TRANSITIONS Angled Entry Mat



Usable width of 36"
(40"-43" overall)
850 lb. weight capacity

- This "welcome" mat, made of recycled rubber, is eco-friendly and can be used safely on nearly any surface.
- Stands up well in harsh outdoor conditions.
- Features beveled edges for cross-traffic access.
- Available in heights ranging from 1/2" to 2-1/2" in quarter-inch increments, and can be trimmed or notched if needed.
- Offered in black or grey, with or without "Welcome".

TRANSITIONS Angled Entry Plate



32" wide
700 lb. weight capacity

- Formed aluminum construction for simple, one-step threshold navigation.
- Three lengths are available (8", 10", and 12") to accommodate thresholds up to 2" high.
- Applied, slip-resistant surface offers great traction.



Lifetime Warranty
www.ezaccess.com

customerservice@ezaccess.com | 1.800.451.1903