

Sorbalgon® and Sorbalgon® Ag


A Calcium Alginate, that is stronger than Aquacel®

Sorbalgon® and **Sorbalgon® Ag** calcium alginate dressings are designed to quickly absorb wound exudate. As wound exudate is absorbed, the alginate forms a gel which assists in maintaining a moist healing environment and keeps the dressing from adhering to the wound bed for intact removal.

Sorbalgon® Ag contains silver particles which release silver ions in the dressing for an antimicrobial effect against a wide range of pathogens including *Staphylococcus aureus*, *Pseudomonas aeruginosa* and others.




Dressing Absorbency¹



The absorbency capacity of **Sorbalgon®** is 20% greater than **Aquacel® Extra**. The absorbency capacity of **Sorbalgon® Ag** is similar to **Aquacel® Extra Ag**. The dressings are designed to be fast and highly absorbent to help stop bleeding.

Minimal Wet Tensile Strength - any direction¹



The tensile strength of **Sorbalgon®** is 5% greater than **Aquacel® Extra**. The tensile strength of **Sorbalgon® Ag** is 30% greater than **Aquacel® Extra Ag**. The dressing removal is easy, without tissue irritation and without leaving gel behind in the wound bed.

HARTMANN provides healthcare professionals with clinical solutions for acute and chronic wounds. We offer a comprehensive advanced wound care portfolio that makes it easy to choose the right product to address healing needs **in the wound, on the wound, and around the wound.** Sorbalgon® and Sorbalgon® Ag are used in the wound.

 In Wound

 On Wound

 Around Wound

Unlike synthetic hydrofiber, **Sorbalgon® and Sorbalgon® Ag** are made from calcium alginate which is derived from natural materials (seaweed). After use, the dressing is biodegradable to help reduce the impact on the planet.

Calcium alginate fibers are hydrophilic, non-adherent and highly absorbent.



Once inside the wound, the fibers turn into a gel, which provides optimal wound moisture balance.



Order Number	Dressing Size	Dressings per Box	HCPCS Code
Sorbalgon®			
49200001	2" x 2" / 5cm x 5cm	10	A6196
49210001	4" x 4" / 10cm x 10cm	10	A6196
49230001	4" x 8" / 10cm x 20cm	5	A6197
49220001	1" x 12" / 2.5cm x 30cm, Ribbon	5	A6199
Sorbalgon® Ag			
999612	2" x 2" / 5cm x 5cm	10	A6196
999611	4" x 4" / 10cm x 10cm	10	A6196
999610	6" x 6" / 15cm x 15cm	5	A6197
999609	1" x 12" / 2.5cm x 30cm, Ribbon	5	A6199



Contact your local HARTMANN representative, or visit us at www.hartmannusa.com for more information.

HARTMANN USA, Inc.
Rock Hill, SC 29730
1-800-243-2294

www.hartmannusa.com

XLIT 2835 Rev. 3 (0820)

