

Confiderm 3.0 Nitrile exam gloves

McKesson Confiderm® 3.0 nitrile exam gloves feature a soft, stretchy material designed to make donning easy by quickly conforming to your hand and facilitating range-of-motion. The gloves also feature improved tactile sensitivity, grip, and fit. Fingertips are lightly textured for decreased slippage during a procedure. Gloves are ambidextrous and non-sterile. They feature a beaded cuff and are available in four sizes.

Features and benefits

- Thin with high tactile sensitivity
- Improved conformability provides superior fit and extended wear comfort
- Textured fingertips for an improved grip
- Dark blue color for easy identification
- Tested per ASTM D6319
- Not made with natural rubber latex



Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

MFR #	DESCRIPTION	QTY
14-6N32	3.0 Nitrile Exam Gloves, Small, Powder-Free	250 ea/bx, 10 bx/cs
14-6N34	3.0 Nitrile Exam Gloves, Medium, Powder-Free	250 ea/bx, 10 bx/cs
14-6N36	3.0 Nitrile Exam Gloves, Large, Powder-Free	250 ea/bx, 10 bx/cs
14-6N38	3.0 Nitrile Exam Gloves, X-Large, Powder-Free	230 ea/bx, 10 bx/cs
14-6N32EC	3.0 Nitrile Exam Gloves, Small, Powder-Free, Textured	100/bx, 10 bx/cs
14-6N34EC	3.0 Nitrile Exam Gloves, Medium, Powder-Free, Textured	100/bx, 10 bx/cs
14-6N36EC	3.0 Nitrile Exam Gloves, Large, Powder-Free, Textured	100/bx, 10 bx/cs
14-6N38EC	3.0 Nitrile Exam Gloves, X-Large, Powder-Free, Textured	100/bx, 10 bx/cs