

# STATLOCK<sup>®</sup>


FOLEY STABILIZATION DEVICE

## “Tape-Free” Stabilization for Foley Catheters

Addresses CDC’s guideline for securing Foley catheters<sup>1</sup>

### Designed to:

- Reduce Foley catheter movement
- Minimize accidental catheter dislodgements
- Maximize patient comfort by:
  - Eliminating circumferential compression
  - Alleviating traction of urethral catheters



# STATLock® Foley Stabilization Device

## STATLock® Device Retainer Design

- Lock-tight design stabilizes catheter to prevent pistoning and accidental dislodgement.
- Swivel design allows catheter movement with the patient without exerting a pull force on the catheter.
- Releasable design allows for easy cleaning of patient.
- Accommodates both 2-way and 3-way latex 8-22 Fr. and silicone 6-26 Fr. catheters for the ultimate in versatility.

## STATLock® Device Anchor Pad Design

- Available in breathable latex-free tricot fabric for added patient comfort.
- Available in both pediatric and adult sizes, giving you the versatility needed for effective Foley catheter stabilization.
- Designed to be easier to use than tape and circumferential leg straps.
- Unique alcohol-soluble adhesive is gentle on the skin.



## Ordering Information

Reference Number	Description	QTY/Box
<b>Adult</b>		
FOL0100	STATLock® Foley Device, Foam Anchor Pad w/Perspiration Holes, for Silicone Catheters	25
FOL0101	STATLock® Foley Device, Foam Anchor Pad w/Perspiration Holes, for Latex Catheters	25
FOL0102	STATLock® Foley Device, Tricot Anchor Pad, for Silicone Catheters	25
FOL0105	STATLock® Foley Device, Tricot Anchor Pad for 3-Way Foley Catheters	25
<b>Pediatric</b>		
VFDSSP	STATLock® Foley Device, Pediatric, Foam Anchor Pad, for Silicone Catheters	25
VFDSP	STATLock® Foley Device, Pediatric, Foam Anchor Pad, for Latex Catheters	25

STATLock® devices are sterile and latex-free.

**BARD | MEDICAL**

1. Centers for Disease Control and Prevention. Guidelines for Prevention of Catheter-associated Urinary Tract Infections. February 1981.

C. R. Bard, Inc.  
8195 Industrial Blvd.  
Covington, GA 30014

Customer Service 800.526.4455  
Fax 800.852.1339

Please consult product inserts and labels for indications, contraindications, hazards, warnings, cautions, and directions for use.  
Bard and STATLock are registered trademarks of C. R. Bard, Inc.  
© Copyright 2010 C. R. Bard, Inc. All rights reserved.  
1401-27 R07/10 THP P04/14 2.5M

**StatLock®**  
Stabilization Devices