



Easy•Touch® U-100 Insulin Syringes

Comfort you can feel!

EasyTouch® Insulin Syringes are designed and manufactured to the highest quality levels possible. Consumers are quick to recognize this high level of quality on a personal level as they experience a more comfortable injection.

The quality level found in our Insulin Syringes can also be found throughout the entire EasyTouch® line of diabetic products.

Gauges

- 27G – 1cc, 5/8" Needle
- 27G – .5cc & 1cc, 1/2" Needle
- 28G – .5cc & 1cc, 1/2" Needle
- 29G – .5cc & 1cc, 1/2" Needle
- 30G – .3cc, .5cc, & 1cc, 1/2" Needle
- 30G – .3cc, .5cc, & 1cc, 5/16" Needle
- 31G – .3cc, .5cc, & 1cc, 5/16" Needle
- 31G – .3cc, .5cc, & 1cc, 1/4" Needle **NEW**




Needle Technology

Our EasyTouch® needles are specifically engineered to provide patients with the maximum level of comfort!

- 1 Triple-Bevel Cut for Superior Sharpness
- 2 Anti-Coring Technology
- 3 Electropolished to Remove Burrs
- 4 Thin-Walled for Maximum Flow
- 5 Film-Coated for a Comfortable Injection

**Contains No Natural Rubber Latex*



Easy•Touch®**U-100 Insulin Syringes – 10 x 10ct. bags****" ● " New Product!**

Item #	Gauge	mL/CC	Length	NRC# 08496-	Box Color	Box Qty	Case Qty	
827158	27 G	1 cc	5/8 in. (16mm)	7020-01	Orange	100/box	5 retail boxes	
827155	27 G	1 cc	1/2 in. (12.7mm)	2715-01	Orange	100/box	5 retail boxes	
827555	27 G	.5 cc	1/2 in. (12.7mm)	2755-01	Orange	100/box	5 retail boxes	
828155	28 G	1 cc	1/2 in. (12.7mm)	2815-01	Purple	100/box	5 retail boxes	
828555	28 G	.5 cc	1/2 in. (12.7mm)	2855-01	Purple	100/box	5 retail boxes	
829155	29 G	1 cc	1/2 in. (12.7mm)	2915-01	Green	100/box	5 retail boxes	
829555	29 G	.5 cc	1/2 in. (12.7mm)	2955-01	Green	100/box	5 retail boxes	
830155	30 G	1 cc	1/2 in. (12.7mm)	3015-01	Blue	100/box	5 retail boxes	
830555	30 G	.5 cc	1/2 in. (12.7mm)	3055-01	Blue	100/box	5 retail boxes	
830355	30 G	.3 cc	1/2 in. (12.7mm)	3035-01	Blue	100/box	5 retail boxes	
830165	30 G	1 cc	5/16 in. (8mm)	3016-01	Blue	100/box	5 retail boxes	
830565	30 G	.5 cc	5/16 in. (8mm)	3056-01	Blue	100/box	5 retail boxes	
830365	30 G	.3 cc	5/16 in. (8mm)	3036-01	Blue	100/box	5 retail boxes	
831165	31 G	1 cc	5/16 in. (8mm)	3116-01	Yellow	100/box	5 retail boxes	
831565	31 G	.5 cc	5/16 in. (8mm)	3156-01	Yellow	100/box	5 retail boxes	
831365	31 G	.3 cc	5/16 in. (8mm)	3136-01	Yellow	100/box	5 retail boxes	
● 831114	31 G	1 cc	1/4 in. (6mm)	7067-01	Yellow	100/box	5 retail boxes	
● 831514	31 G	.5cc	1/4 in. (6mm)	7068-01	Yellow	100/box	5 retail boxes	
● 831314	31 G	.3cc	1/4 in. (6mm)	7069-01	Yellow	100/box	5 retail boxes	

Easy•Touch®**U-100 Insulin Syringes – Individually Wrapped (100ct.)****" ● " New Product!**

Item #	Gauge	mL/CC	Length	NRC# 08496-	Box Color	Box Qty	Case Qty	
627155	27 G	1 cc	1/2 in. (12.7mm)	6510-01	Orange	100/box	5 retail boxes	
627555	27 G	.5 cc	1/2 in. (12.7mm)	6520-01	Orange	100/box	5 retail boxes	
628155	28 G	1 cc	1/2 in. (12.7mm)	6511-01	Purple	100/box	5 retail boxes	
628555	28 G	.5 cc	1/2 in. (12.7mm)	6521-01	Purple	100/box	5 retail boxes	
629155	29 G	1 cc	1/2 in. (12.7mm)	6512-01	Green	100/box	5 retail boxes	
629555	29 G	.5 cc	1/2 in. (12.7mm)	6522-01	Green	100/box	5 retail boxes	
630165	30 G	1 cc	5/16 in. (8mm)	7040-01	Blue	100/box	5 retail boxes	
630565	30 G	.5 cc	5/16 in. (8mm)	7038-01	Blue	100/box	5 retail boxes	
630365	30 G	.3cc	5/16 in. (8mm)	7042-01	Blue	100/box	5 retail boxes	
631165	31 G	1 cc	5/16 in. (8mm)	7044-01	Yellow	100/box	5 retail boxes	
631565	31 G	.5 cc	5/16 in. (8mm)	7037-01	Yellow	100/box	5 retail boxes	
631365	31 G	.3cc	5/16 in. (8mm)	7045-01	Yellow	100/box	5 retail boxes	
● 631114	31 G	1cc	1/4 in. (6mm)	7050-01	Yellow	100/box	5 retail boxes	
● 631514	31 G	.5cc	1/4 in. (6mm)	7065-01	Yellow	100/box	5 retail boxes	
● 631314	31 G	.3cc	1/4 in. (6mm)	7066-01	Yellow	100/box	5 retail boxes	