


Adhesive Bandages



Curity™ Fabric Bandages

- Flexible fabric stretches and conforms with movement
- Ideal for fingers, toes, elbows and other bending areas
- Not made with natural rubber latex


Cat. No.	Style	Bandage Dim.	Pad Dim.	Sterility	Pkg. string	Case qty.
44100	Flexible	¾ in. x 3 in.	¾ in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 24 bx/cs	1,200
44101-	Flexible	1 in. x 3 in.	1 in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 24 bx/cs	1,200
44102-	Flexible XL	2 in. x 3 ¾ in.	1 ¼ in. x 1 7/8 in.	Sterile	1 ea/po, 50 po/bx, 12 bx/cs	600
44106	Knuckle	1 ½ in. x 3 in.	1 in. x 7/8 in.	Sterile	1 ea/po, 30 po/bx, 40 bx/cs	1,200
44107	Spot	7/8 in.	3/8 in. x 3/8 in.	Sterile	1 ea/po, 50 po/bx, 24 bx/cs	1,200
44109	Oval	1 ¼ in. x 1 in.	½ in. x ½ in.	Sterile	1 ea/po, 50 po/bx, 24 bx/cs	1,200



Curity™ Plastic Bandages

- Perforated vinyl film and porous adhesive allow skin to breathe
- Not made with natural rubber latex

Cat. No.	Type	Bandage Dim.	Pad Dim.	Sterility	Pkg. string	Case qty.
44112	Plastic	¾ in. x 3 in.	¾ in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 24 bx/cs	1,200
44113-	Plastic	¾ in. x 3 in.	¾ in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 72 bx/cs	3,600
44114	Plastic	1 in. x 3 in.	1 in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 24 bx/cs	1,200
44115	Plastic	1 in. x 3 in.	1 in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 72 bx/cs	3,600
44116	Patch	1 ½ in.	5/8 in. x 5/8 in.	Sterile	1 ea/po, 50 po/bx, 24 bx/cs	1,200



Curity™ Sheer Bandages

- Thin plastic material blends with skin to make the bandage less noticeable
- Not made with natural rubber latex

Cat. No.	Type	Bandage Dim.	Pad Dim.	Sterility	Pkg. string	Case qty.
44118	Sheer	¾ in. x 3 in.	¾ in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 72 bx/cs	3,600
44119	Sheer	1 in. x 3 in.	1 in. x 7/8 in.	Sterile	1 ea/po, 50 po/bx, 72 bx/cs	3,600
44120	Sheer	7/8 in.	3/8 in. x 3/8 in.	Sterile	1 ea/po, 100 po/bx, 36 bx/cs	3,600