

MD by Drive Air-Cane

HCPCS:




FIGURE A



FIGURE B



FIGURE C

Most canes on the market are commodities the Air-Cane breaks the mold with Airgonomic technology that directly addresses the most common user complaint: hand and wrist discomfort. Its patented Air-Handle cushions the hand to reduce strain, while the FlexGrip Air-Tip enhances stability and confidence with every step.

- **Airgonomic Comfort:** Adaptive Air-Handle features an integrated air bladder that cushions the grip.
- **Improved Traction & Stability:** Patented FlexGrip Air Tip grips the earth and compresses to soften every step.
- **Compact Packing & Travel:** Collapsible 3-part design folds in seconds.
- **Fits Most Users:** Easy to adjust for individuals 4'9" - 6'7" up to 350 lb.

PRODUCTS	Description
MD20100BK	Air Cane, Onyx
MD20100BL	Air Cane, Sapphire
MD20100GN	Air Cane, Forest
MD20100GR	Air Cane, Slate
MD20100PR	Air Cane, Lavender
MD20100RD	Air Cane, Crimson

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	1 lb
Product Dimensions	5.5 (W) x 32 - 37 (H)
Material	Anodized Aluminum
Handle Height	32 to 37
Weight Capacity	350 lb
Height	Adjustable for users 49 - 67
Retail Box Dimensions	6 (L) x 2.5 (W) x 16 (H)