

MEDNEB+ Compressor Nebulizer

HCPCS: E0570+A7003-1,E0570+A7003-1



FIGURE A

- Includes Airial compressor, disposable and optionable reusable nebulizer kit, 7' kink resistant tubing, mouthpiece, and extra filters
- Compact, Desk-Top Unit
- Piston-Driven System
- Optional carry bag with MQ5700B model

PRODUCTS	Description
MQ5600	Compressor Nebulizer - MEDNEB 1/ea, 8/cv
MQ5700	MEDNEB+ w/Reusable Neb Cup 1/ea, 8/cv
MQ5700B	MEDNEB+, Reusable Neb Cup 1/ea, 8/cv

SPECIFICATIONS	Description
Warranty	Limited 5 Year
Compressor Free Air Flow	5-8 LPM
Operating Humidity	15-90% non-condensing
Operating Voltage	120VAC, single phase, 60Hz
Power Consumption	180 VA
Size	6.7" x 5.4" x 3.5"
Storage/Transport Humidity	10-90%
Weight	2.6 lbs

ACCESSORIES	Description
18254	Reusable Neb Kit 10/case, CS
MASK 001A	Adult Aerosol mask CS/50, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS