

Assura® Flat MIDI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 350 mL

Units: 30

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
12130	12140	Cut-to-fit	$1\frac{3}{16}$ - $2\frac{1}{2}$ " (20-65 mm)	A4416
12134	12144	Pre-cut	1" (25 mm)	A4416
12135	12145	Pre-cut	$1\frac{3}{16}$ " (30 mm)	A4416
12136	12146	Pre-cut	$1\frac{3}{8}$ " (35 mm)	A4416
NA	12147	Pre-cut	$1\frac{9}{16}$ " (40 mm)	A4416
NA	12149	Pre-cut	2" (50 mm)	A4416

Assura® Flat MINI Closed Pouch

- Standard wear
- **Filter**
- Soft cloth on front and back
- **Not made with natural rubber latex**


Opaque 5" (12 cm) / 170 mL

Units: 30

Code	Cut	Stoma Size	HCPCS*
12110	Cut-to-fit	$1\frac{3}{16}$ - $2\frac{1}{2}$ " (20-65 mm)	A4416

Assura® Convex Light MAXI Closed Pouch



- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



8 ½" (22 cm) / 560 mL

Units: 10

Code	Code	Cut	Stoma Size	HCPCS*
Transparent	Opaque			
14441	14444	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A4417
14442	14445	Cut-to-fit	$\frac{5}{8}$ - $1\frac{11}{16}$ " (15-43 mm)	A4417

Indicates convexity

Assura® Convex Light MIDI Closed Pouch

- Standard wear
- **Filter**
- Belt tabs for use with 4215, 4220
- Soft cloth on back (and front of opaque)
- **Not made with natural rubber latex**



7" (18 cm) / 350 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
------	-----	------------	--------

Transparent Opaque

14431	14434	Cut-to-fit	$\frac{5}{8}$ - $1\frac{5}{16}$ " (15-33 mm)	A4417
14432	14435	Cut-to-fit	$\frac{5}{8}$ - $1\frac{11}{16}$ " (15-43 mm)	A4417
NA	14451	Pre-cut	$\frac{7}{8}$ " (22 mm)	A4417
NA	14452	Pre-cut	1" (25 mm)	A4417
NA	14453	Pre-cut	$1\frac{1}{8}$ " (28 mm)	A4417
NA	14454	Pre-cut	$1\frac{1}{4}$ " (31 mm)	A4417

Assura® Flat Stoma Cap

- Standard Wear
- **Filter**
- Belt tabs
- Soft cloth on back
- Small pouch size for intimate moments or colostomates who irrigate
- **Not made with natural rubber latex**

**Opaque**

Units: 30

Code	Cut	Stoma Size	HCPCS*
2501	Cut-to-fit	$1\frac{3}{16}$ - $2\frac{1}{4}$ " (20-55 mm)	A5055

Urostomy**Assura® Flat MAXI Urostomy Pouch**

- Standard wear
- Multi-chambers minimize sloshing
- Soft outlet can be connected to urostomy night bag #21365 or any standard night bag
- Soft cloth on back
- **Not made with natural rubber latex**

**Transparent** 10 $\frac{3}{4}$ " (27 $\frac{1}{2}$ cm) / 470 mL

Units: 10

Code	Cut	Stoma Size	HCPCS*
14222	Cut-to-fit	$\frac{3}{8}$ - $2\frac{1}{4}$ " (10-55 mm)	A5071

 Indicates convexity

*For important information regarding the use of HCPCS codes listed and code definitions see page 180