

Delta Pro Homecare Bed System

HCPCS:



- Patent pending Universal Frame can configure with various component to create up to SIX unique bed configurations offering more homecare bed options without the need to store more beds!
- Assists with inventory constraints by offering 1 homecare bed system that can transform to suit multiple needs using color coded component options
- Product Safety & Performance Standards: Meets IEC 60601-2-52; 60601-1; 60601-1-2; 60601-1-11
- Safe working load capacity of 450 lbs
- Choose between semi-electric and full electric motor options to provide head, foot, and height adjustability of bed frame as needed (excludes 15003PN, manual crank homecare bed system)
- Ultra lightweight head and foot frames easy to set up by 1 individual with no links to connect and no tools required
- Square tube steel frame slat deck provides durable strength, no-sag mattress support, and is easy to clean
- HD Polyethylene molded bed ends are attractive, labeled for proper assembly, and easy to clean and maintain
- Bed End safety lock slides over hanger pins to prevent disengaging from bed ends
- UL approved motor is self contained for reduced weight and noise and can be easily installed and removed with patient in bed (excludes 15003PN)
- Reduce service calls with large glow in the dark hand pendant that is easy-to-use for adjusting head and foot deck sections; built in battery backup and LED indicators for power and lockout included (excludes 15003PN)
- Power failure protection available on dual style motors via a 9-volt battery (not included) located in the hand pendant which can return the frame to a horizontal position. Manual cranks included for all styles (excludes 15003PN, no electrical components means independent of electrical power availability)
- Frame includes mounted retractable full rail receiver brackets to eliminate cross braces ensuring proper placement and reducing equipment loss
- Half (15561HRAIL) and full (15561FRAIL) rails are available, sold separately, and meet FDA entrapment guidelines as well as ASTM standards
-
-

PRODUCTS

Description

15003PN	Delta Pro Bed - Manual
15004PN	Delta Pro Bed-Standard Linear Semi Elect
15005LPN	Delta Pro Low Linear Full Electric
15030N	Delta Pro Bed -Standard Semi Electric
15033NC	Delta Pro Bed - Standard Full Electric
15235N	Delta Pro Bed -Low Full Electric

SPECIFICATIONS

Description

Battery	15030N, 15033NC, 15235N: Backup available but not included
Bed Height	15030N: 13.5" - 23" 15033NC: 13.5" - 23" 15235N: 10" - 23" 15003PN: 13.5" - 23" 15004PN: 13.5" - 23" 15005LPN: 10" - 23"

SPECIFICATIONS	Description
Bed Height w/3 Casters	15030N, 15033NC, 15003PN and 15004PN: 28" - 38" 15235N and 15005LPN: 27" - 40"
Carton Shipping Weight	15030N: 160 lbs 15033NC: 168 lbs 15235N: 171 lbs 15003PN: 166 lbs 15004PN: 173 lbs 15005LPN: 181 lbs
Casters	3" casters (2 locking and 2 unlocking)
Compressor Free Air Flow	Brown
Motor	15030N: Dual Semi- Electric 15033NC and 15235N: Dual Full Electric w/Hi Lo 15003PN: No Motors 15004PN: Linear Semi- Electric 15005LPN: Linear Full Electric w/Hi lo
Sleep Surface	36" W x 80" L
Warranty	Motor and Hand Control: 5 years Welds and Frame: Limited Lifetime All other parts and components: 1 Year
Weight Capacity	Safe working load 450lbs (to include occupant and any accessories)
Material Composition	Bed Frame: Steel Bed Ends: HD Polyethylene Motors: Electrical Components
ACCESSORIES	Description
15550MOHC	Delta Pro Dual Motor w/pendant, EA
15561FRAIL	Compet II Full Rail, 1/PR, PR
15561HRAIL	Compet II Half Rail, 1/PR, PR