

The First Responder solution for active patient warming in any casualty situation:



Emergency Rescue




Disaster Preparedness



Air Transport


Ready-Heat II

- The #1 selling active warming blanket
- Quick warming blanket, with 8+ hours of active warmth
- Ideal for all patient transport and emergency rescue
- Water resistant top layer made of coated medical grade non-woven material and a foam core




Ready-Heat II Torso

- Compact size and smaller footprint for greater portability
- Deploys in seconds and heats in minutes
- Supplies consistent heat for 8 hours
- 5 year shelf life
- Water resistant top layer made of coated medical grade non-woven material and a foam core



Ready-Heat Infant Cocoon

- Designed to actively warm the most precious cargo
- Consists of 6 thermal pads strategically placed to warm infant's body
- Hood assists in the retention of body warmth and the padded foam core and dual hook and loop design provide maximum comfort and fit



For the First Responder

Stay warm on the job with proprietary Ready-Heat technology in a disposable heated vest.

Include Ready-Heat in Your:

- First Aid Kits
- Wilderness Survival Kits
- Mountain Rescue Kits
- Pediatric Transport Kits
- Water Rescue Kits
- Trauma Kits



MANUFACTURED BY TECHTRADE LLC IN THE USA 

Visit our website for more detailed information and IFUs - www.ready-heat.com/IFU
www.ready-heat.com • info@ready-heat.com • 1-866-491-TECH • F: +1 201 706 8130


0459



ISO 13485:2016



4 PANEL BLANKET

NSN #: 6532-01-525-4062

Part Number / Per Case: S4RHMD / Case of 8



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 48" / 86 cm x 122 cm


Weight: 1.75 lbs. / 794 gm

Description: Torso Sized Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250 lbs./115 kg). Military Medical Protocol.

6 PANEL BLANKET

NSN #: 6532-01-525-4063

Part Number / Per Case: S6RHLG / Case of 8



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 60" / 86 cm x 152 cm

Weight: 2.5 lbs. / 1.13 kg

Description: Full Body Blanket. Used to Treat Shock and Hypothermia. Suitable as an evacuation litter (up to 250lbs./115 kg). Number 1 choice of US Air Force transports.

1 PANEL BLANKET

NSN #: 6532-01-538-1525

Part Number / Per Case: S1RHSM / Case of 12



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 16" x 18" / 41 cm x 46 cm


Weight: 7 oz. / 198 gm

Description: Mini Blanket. Suitable for injured extremities. Fully adjustable with 2 hook and loop closures. Perfect for First Aid and Survival Kits.

READY-HEAT II BLANKET

NSN #: 6532-01-538-0759

Part Number / Per Case: G12RHZ / Case of 8



Warming Temperature:
100°F / 38°C in 8-10 Minutes

Maintains Heat: Up to 8 Hours


Size: 34" x 48" / 86 cm x 122 cm

Weight: 1.5 lbs. / 794 gm

Description: First Responder Blanket. Used to Treat Shock, Trauma and Pre-Hypothermia. Has a foam core for water resistance and is lightweight and comfortable. Quick to reach operating temperature. Works well at high altitudes.

READY-HEAT II TORSO BLANKET

Part Number / Per Case: G6RH2T / Case of 12



Warming Temperature:
100°F / 38°C in 8-10 Minutes

Maintains Heat: Up to 8 Hours

Size: 34" x 28" / 86 cm x 61 cm

Weight: 14 oz. / 397 gm

Description: First Responder Blanket. Used to treat shock, trauma and pre-hypothermia. Smaller overall footprint both in bag and on patient. All the benefits of the full-sized blanket, G12RHZ.

READY-HEAT PERSONAL USE VEST

Part Number / Per Case: GRHV06A / Case of 12



Warming Temperature:
100°F / 38°C in 10-15 Minutes

Maintains Heat: Up to 8 Hours


Size: Overhead Style with Adjustable Hook and Loop Straps. One Size Fits All.

Weight: 10 oz. / 284 gm

Description: Designed to keep the First Responder warm during long missions. Also great for outdoor activities (Hiking, Skiing, Camping, Hunting, Fishing, Marathons and Mountain Climbing).

INFANT WARMING COCOON

Part Number / Per Case: GIW6C / Case of 12



Warming Temperature:
100°F / 38°C in 15 Minutes

Maintains Heat: Up to 8 Hours


Size: 31.5" x 28" / 80 cm x 71 cm

Weight: 12 oz. / 340 gm

Description: Award Winning Infant Warming Cocoon is designed for the most precious cargo. Six warming pads to cover the vital organs and two adjustable hook and loop fasteners for size adjustment and attached hood for added heat retention.

TEMPERATURE MANAGEMENT FULL BODY

Part Number / Per Case: SB9RH9120 / Case of 10



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours


Size: 34" x 60" / 86 cm x 152 cm

Weight: 18 oz. / 510 gm

Description: Patient Warming Blanket. Fully mobile, can move with the patient from pre-op to intra-op and post-op. Prevents hypothermia, increases patient comfort for higher quality of care, fully disposable, self-contained system.

TEMPERATURE MANAGEMENT HALF BODY

Part Number / Per Case: SB6RH9120 / Case of 20



Warming Temperature:
104°F / 40°C in 15-20 Minutes

Maintains Heat: Up to 10 Hours

Size: 34" x 24" / 86 cm x 61 cm

Weight: 12 oz. / 350 gm

Description: Patient Warming Blanket. Fully mobile, can move with the patient from pre-op to intra-op and post-op. Half-sized blanket can be moved and adjusted for different types of surgeries or access. All the benefits of the full-sized blanket, SB9RH9120.

MANUFACTURED BY TECHTRADE LLC IN THE USA 

Visit our website for more detailed information and IFUs - www.ready-heat.com/IFU
www.ready-heat.com • info@ready-heat.com • 1-866-491-TECH • F: +1 201 706 8130


0459


ISO 13485:2016

