

Coverflex®Grip



Spec.-No.: D 9.4050

Department: -

Date: 2018-08-22

1. General Product Description

Coverflex®Grip is an elastic tubular bandage.

Coverflex®Grip carries the CE mark according to EU directive 93/42/EEC for medical devices.

The product is classified as a class I medical device.

A conformity assessment has been performed for Coverflex®Grip and it has been shown to be in compliance with all applicable requirements of the above-mentioned directive.

The safe use of Coverflex®Grip, therefore, is ensured if the product is used in line with the intended purpose (see 2. Application/Indication).

2. Application / Indication

The tubular compression bandages are intended to be used by medical professionals and patients for secondary fixation,
as skin protection under cast
and support to joints and ligaments.

The Compression effect of the product is because of material and product design but not a product which is intended for compression therapy.

3. Presentations

Coverflex®Grip is supplied in rolls of 10 m in dispenser boxes.

	Size	Width in cm
Infant legs & arms	A	4,00
Small hands	B	6,25
Hands, feet & arms	C	6,75
Arms & legs	D	7,50
Legs & knees	E	8,75
Thighs & knees	F	10,00
Large thighs	G	12,00
Small trunks	J	18,00
Medium trunks	K	22,00
Large trunks	L	32,50
X-Large trunks	M	36,00

Coverflex®Grip



Spec.-No.: D 9.4050

Department: -

Date: 2018-08-22

4. Product Characteristics

Material composition:

82 % Unbleached cotton

9 % Elastane

9 % Polyamide

5. Product Requirements

Sterilizing capability: Steam sterilization according to DIN EN ISO 17665 with validated processes
ETO sterilization according to EN ISO 11135 with validated processes

6. Labelling

Lot-No. with 9-digit code:

Lot.-No.	3	999	01	00	1
Explanation	[3 = 2013]	[999 – 000]	Internal Code		
Definition	Last figure	Serial production order			

Expiry date:.....Year – Month or Year – Month – Day

Example:..... 2020-03.....or..2020 – 03 - 01

Shelf Life.....5 years

7. Packaging

Coverflex®Grip 10 m rolls are packed individually in dispensers. The dispensers are packed in transport cartons, sealed with adhesives tapes and packed onto euro pallets.

Date: 22nd August 2018



-Hospidana ApS-