




SB2t

Dimensions = Inches



* Footrests rotated inwards




Product features:

- Weight capacity - 350lb
- Tilt angle 0 - 30 degrees, infinite positioning
- Rolls over most standard toilets
- Low maintenance
- Arms - removable / lockable / flip-back
- Height adjustable
- Safety mechanisms & side stabilizers on bath base
- Footrests - removable and swing-away

- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap & chest belts plus calf strap
- Caster Wheels - 5" x 4 locking
- Commode system
- Neck rest, fully adjustable
- Tool-free assembly

Bird's eye view



T = 3.50 inch

A = 6.75 inch

S = 4.25 inch

Accessory items

- Refer to the accessory pages or on our website