




SB2t

Dimensions = Millimetres



* Footrests rotated inwards




Product features:

- Weight capacity - 160kg
- Tilt angle 0 - 30 degrees, infinite positioning
- Rolls over most standard toilets
- Low maintenance
- Arms - removable / lockable / flip-back
- Height adjustable
- Safety mechanisms & side stabilizers on bath base
- Footrests - removable and swing-away

- No rust aluminum frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts - lap & chest belts plus calf strap
- Caster Wheels - 125mm x 4 locking
- Commode system
- Neck rest, fully adjustable
- Tool-free assembly

Bird's eye view



T = 90mm
A = 170mm
S = 110mm

Accessory items

- Refer to the accessory pages or on our website