

Hoyer® Loop Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

Lift Types		Spreader Bar	Sling SWL
Advance Elara HPL450 HPL500	HPL700 Presence Stature	Loop 	500 lbs.

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
Access  Material: Padded	XS	Brown	NA1600									
	S	Red	NA1601									
	M	Yellow	NA1603									
	L	Green	NA1606									
	XL	Blue	NA1608									
Access, Head Support  Material: Padded	S	Red	NA1626									
	M	Yellow	NA1627									
	L	Green	NA1628									
	XL	Blue	NA1629									
Quickfit  Material: Padded	XS	Brown	NA1002P									
	S	Red	NA1000P									
	M	Yellow	NA1003P									
	L	Green	NA1006P									
	XL	Blue	NA1008P									
Quickfit  Material: Mesh	S	Red	NA1001									
	M	Yellow	NA1004									
	L	Green	NA1007									

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
Full Back  Material: Padded	XS	Brown	NA1067									
	S	Red	NA1068									
	M	Yellow	NA1069									
	L	Green	NA1070									
	XL	Blue	NA1071									
Full Back  Material: Mesh	XS	Brown	NA1080									
	S	Red	NA1081									
	M	Yellow	NA1082									
	L	Green	NA1083									
	XL	Blue	NA1084									
Full Back Disposable  Material: Disposable	XS	Brown	NA25620									
	S	Red	NA25621									
	M	Yellow	NA25622									
	L	Green	NA25623									
	XL	Blue	NA25624									
Full Body Hammock  Material: Polyester	S	Red	NA3171									
	M	Yellow	NA3172									
	L	Green	NA3173									

Hoyer® Loop Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

Lift Types		Spreader Bar	Sling SWL
Advance Elara HPL450 HPL500	HPL700 Presence Stature	Loop 	500 lbs.

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
Quickfit Deluxe  Material: Polyester	XS	Brown	NA1048									
	S	Red	NA1050									
	M	Yellow	NA1053									
	L	Green	NA1056									
Quickfit Deluxe  Material: Mesh	XL	Blue	NA1058									
	S	Red	NA1051									
	M	Yellow	NA1054									
	L	Green	NA1057									
Quickfit Deluxe  Material: Padded	XL	Blue	NA1059									
	S	Red	NA1050P									
	M	Yellow	NA1053P									
	L	Green	NA1056P									

*All slings on this page may be used for amputees

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
Long Seat  Material: Polyester	S	Red	NA1100									
	M	Yellow	NA1103									
	L	Green	NA1106									
	XL	Blue	NA1109									
Long Seat  Material: Mesh	S	Red	NA1101									
	M	Yellow	NA1104									
	L	Green	NA1107									
Long Seat Commode  Material: Polyester	S	Red	NA1125									
	M	Yellow	NA1128									
	L	Green	NA1131									
Long Seat Commode  Material: Mesh	S	Red	NA1126									
	M	Yellow	NA1129									
	L	Green	NA1132									

NA2650 Extension Straps available add length to the existing sling straps

Hoyer® Clip Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

Lift Types	Cradle	Sling SWL
Presence Stature	Clip 	500 lbs.

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
Comfort Standard  Material: Polyester	XS	Brown	NA25501									
	S	Red	NA25502									
	M	Yellow	NA25503									
	L	Green	NA25504									
	XL	Blue	NA25505									
Comfort  Material: Mesh	S	Red	NA25512									
	M	Yellow	NA25513									
	L	Green	NA25514									
	XL	Blue	NA25515									
Comfort  Material: Padded	XS	Brown	NA25506									
	S	Red	NA25507									
	M	Yellow	NA25508									
	L	Green	NA25509									
Comfort Disposable  Material: Disposable	XS	Brown	NA25610									
	S	Red	NA25611									
	M	Yellow	NA25612									
	L	Green	NA25613									
	XL	Blue	NA25614									

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
Comfort Access  Material: Padded	S	Red	NA25571									
	M	Yellow	NA25572									
	L	Green	NA25573									
	XL	Blue	NA25574									
Comfort Access Head Support  Material: Padded	M	Yellow	NA25577									
	L	Green	NA25578									
	XL	Blue	NA25579									
Comfort Amputee*  Material: Polyester	S	Red	NA1201									
	M	Yellow	NA1202									
	L	Green	NA1203									
Comfort Sling Accessory			Style		Size		Part Number					
Comfort Pad 			Padded		One Size		NA25517					
			Standard		One Size		NA25516					

*Must be used with HOY_AMPUTEE_SB

Hoyer® Loop and Clip Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

Lift Types	Cradle	Sling SWL
Ascend** Elevate Journey	Loop 	440 lbs.

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	440
Deluxe Stand-Aid  Material: Polyester	S	Red										
	M	Yellow										
	L	Green										
Deluxe Transport  Material: Polyester	S	Red										
	M	Yellow										
	L	Green										

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	440
Standing Disposable  Material: Disposable	S	Red										
	M	Yellow										
	L	Green										

**Ascend is discontinued - for reference only

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	440
Stand-Aid Clip*  Material: Polyester	M	Yellow										
	L	Green										


*Elevate only. Must have 0Y0488 type spreader bar for clip style sling.

Lift Types	Cradle	Sling SWL
Elevate	Clip 	440 lbs.

Hoyer® Loop Style Sling Color/Size Chart - Bariatric Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

Lift Types	Spreader Bar	Sling SWL
Calibre HPL700	Bariatric Loop 	850 lbs.

Sling	Size	Color	Color Example	Part Number
Quickfit  Material: Polyester 	M	Yellow with White Stripes		NA1009P
	L	Green with White Stripes		NA1010P
	XL	Blue with White Stripes		NA1011P
	XXL	Black with White Stripes		NA1012P
Full Back  Material: Polyester 	M	Yellow with White Stripes		NA1072
	L	Green with White Stripes		NA1073
	XL	Blue with White Stripes		NA1074
	XXL	Black with White Stripes		NA1075
Hammock*  Material: Spacer 	XL	Red with White Stripes		NA1030
	XXL	Yellow with White Stripes		NA1031
	3XL	Green with White Stripes		NA1032
	4XL	Blue with White Stripes		NA1033
	5XL	Black with White Stripes		NA1034


*Only smaller sizes should be used with the HPL700 due to clearance over beds.

** Due to the varying nature of bariatric patient shapes and sizes, it becomes increasingly difficult to provide a guide on sling size based on weight alone. Patients should be individually assessed to ascertain the correct size sling.

Hoyer® Loop Style Sling Color/Size Chart

Recommended weight range should only be used as a guideline as other factors can affect sizing.

Lift Types		Spreader Bar
Advance Elara HPL450	HPL500 Presence Stature	Loop 
Calibre	HPL700	Bariatric Loop 

Sling	Item Number	Sling SWL	Sling	Item Number	Sling SWL	Sling	Item Number	Sling SWL
Repositioning, Standard 	NA25710	1000 lbs.	Twin Turner, Regular 	NA25705	625 lbs.	Limb, Multi Loop 	NA25700	400 lbs.

Sling	Sling SWL	Size	Color	Recommended Weight Range (lbs)									
				35	75	100	125	150	175	200	275	350	500
Standing Harness  Material: Polyester 	500 lbs	XS	Brown	NA26600									
		S	Red		NA26601								
		M	Yellow				NA26602						
		L	Green						NA26603				

Sling	Size	Part Number
Slide Sheet  Material: Polyester 	S	NA25525
	M	NA25526
	L	NA25527
	XL	NA25528
	XXL	NA25529

Hoyer® Loop Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

Lift Types		Spreader Bar	Sling SWL
HML400 HPL402	CHLA-2 CHLA-2T	Loop 	600 lbs.

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
U Sling Head Support  Material: Padded	S	Red			70013							
	M	Yellow				70012						
	L	Green						70011				
	XL	Blue								70010		
U Sling Head Support  Material: Mesh	S	Red			50027							
	M	Yellow				50026						
	L	Green						50025				
U Sling  Material: Padded	S	Red			70003							
	M	Yellow				70002						
	L	Green						70001				
	XL	Blue								70000		

Sling	Size	Color	Recommended Weight Range (lbs)									
			35	75	100	125	150	175	200	275	350	500
U Sling  Material: Mesh	S	Red			50023							
	M	Yellow				50022						
	L	Green						50021				
	XL	Blue								50020		
One Piece  Material: Polyester	S	Red			70053							
	M	Yellow				70052						
	L	Green						70051				
One Piece Commode  Material: Polyester	S	Red			70058							
	M	Yellow				70057						
	L	Green						70056				
	XL	Blue								70055		

*One Piece slings may be used for amputees