



HEALTHCARE

IMPROVING QUALITY OF LIFE
– BECAUSE WE CARE





Shaped Pads

Maximum comfort, protection, and discretion

ABENA Light Premium and ABENA San Premium are a range of thin, anatomically shaped pads suitable for people living with light to moderate incontinence. All pads have peel-off tape on the back to ensure the product stays in place. Combined with our odor control system and leakage barriers, ABENA Light and ABENA San allow users to feel protected all day and night.

- High-performance core and Top-Dry Technology ensure a dry feeling
- Invisible leakage barriers for extra protection
- Anatomically shaped for a comfortable fit
- Fully breathable materials support a healthy skin balance
- Odor Control system for maximum discretion



Learn more and find

Incontinence

Shaped Pads



Shaped pad, ABENA Light



Item number	Size	Dimensions	Absorption Rothwell	Quality level	Pcs per pk / pk per cll
1999905382	Ultra Mini 0	19 × 7,5 cm	80 ml	Premium	24 / 12

ABENA Light is our Premium range of anatomically shaped pads with self-adhesive tape developed for light to moderate incontinence. The efficient core with Top-Dry Technology provides fast absorption and a dry surface, even after repeated urination, and ensures optimal skin balance and increases comfort. The fully breathable products feature the Odor control system that minimizes the risk of unpleasant smells. The dermatologically tested products carry the Nordic Swan Ecolabel and are made from PEFC-certified materials.



ABENA Light



Item number	Size	Properties	Dimensions	Absorption Rothwell	Quality level	Pcs per pk / pk per cll
1000021337	Mini 1	Single packed	22 × 10 cm	180 ml	Premium	20 / 10
1000021338	Mini Plus 1A	Single packed	28 × 10 cm	200 ml	Premium	16 / 16
1000021339	Normal 2	Single packed	26 × 11 cm	350 ml	Premium	12 / 12
1000021340	Extra 3	Single packed	33 × 11 cm	500 ml	Premium	10 / 20
1000021341	Extra Plus 3A	Single packed	33 × 11 cm	650 ml	Premium	10 / 20
1000021342	Super 4		42 × 20 cm	850 ml	Premium	30 / 9
1999905383	Maxi 4A	Single packed	41 × 15 cm	1000 ml	Premium	8 / 6

ABENA Light is our Premium range of anatomically shaped bladder protection pads for light to moderate incontinence with adhesive release tape. The high performance core with Top-Dry Technology ensures skin integrity and optimal skin care by keeping the exposed area comfortable and dry even after multiple voids. The fully breathable pads are fitted with invisible lying leakage barriers or has a bowl shape due to elasticated foam edges offering extra protection against leakage. When combining these 3 functionalities- breathability, protection and Top-Dry- it keeps the skin in balance and maintains a natural pH skin balance. The Odor control system minimizes unpleasant smells. The dermatologically tested products have been given the Nordic Swan Ecolabel and are made from FSC™ certified materials.



Shaped pad, ABENA San



Item number	Size	Properties	Dimensions	Absorption Rothwell	Quality level	Pcs per pk / pk per cll
1000021301	1	Single packed	22 × 10 cm	200 ml	Premium	30 / 10
1000021302	1A	Single packed	28 × 10 cm	250 ml	Premium	28 / 12
1000021303	2		26 × 10 cm	350 ml	Premium	30 / 8
1000021304	3		33 × 11 cm	500 ml	Premium	28 / 7
1000021305	3A		33 × 11 cm	600 ml	Premium	28 / 7
1000021306	4		42 × 20 cm	800 ml	Premium	30 / 9

ABENA San 1-4 is our Premium range of anatomically shaped bladder protection pads for light to moderate incontinence with adhesive release tape. The high performance core with Top-Dry Technology ensures fast absorption and a dry surface even after multiple voids. The fully breathable pads are fitted with invisible lying leakage barriers for extra protection. When combined, these functionalities help protect the skin and maintain a natural pH skin balance. The Odor control system minimizes





ABENA San

Optimum leakage protection

ABENA San Premium in sizes 5–11 is a range of anatomically shaped pads that are used together with special fixation pants. This two-piece system ensures that the pad stays in place while reducing the risk of leakages. Further, ABENA's life cycle assessment has shown that the two-piece system uses fewer raw materials and has less environmental impact as the fixation pants are re-useable.

- High-performance core with Top-Dry Technology ensures a dry feeling
- 360° leakage barriers reduce workload and improve user comfort
- Wetness indicator with graduation scale to help guide product changes
- Fully breathable materials support a healthy skin balance
- Two-piece system minimizes the usage of raw materials



Learn more and find

Incontinence

Shaped Pads



Shaped pad, ABENA San



Item number	Size	Dimensions	Absorption Rothwell	Quality level	Pcs per pk / pk per cll
1000021307	5	54 × 28 cm	1200 ml	Premium	36 / 4
1000021308	6	63 × 30 cm	1600 ml	Premium	34 / 3
1000021309	7	63 × 36 cm	2100 ml	Premium	30 / 4
1000021310	8	63 × 36 cm	2500 ml	Premium	22 / 4
1000021311	9	73 × 37 cm	2400 ml	Premium	28 / 4
1000021312	10	73 × 37 cm	2800 ml	Premium	25 / 4
1000021313	11	73 × 37 cm	3400 ml	Premium	21 / 4

ABENA San 5-11 is our Premium range of anatomically shaped pads for moderate to heavy incontinence. The high performance core with Top-Dry Technology ensures skin integrity and optimal skin care by keeping the exposed area comfortable and dry even after multiple voids. The fully breathable pads are fitted with standing leakage barriers combined with front and rear built-in cross barriers offering 360° leakage protection. When combined, these functionalities help protect the skin and maintain a natural pH skin balance. The Odor control system minimizes unpleasant smells. The wetness indicator with graduation scale guides the user and caretaker to the product with the right absorption level and product changes. The dermatologically tested products have been given the Nordic Swan Ecolabel. ABENA San in sizes 5-11 should always be used with special fixation underwear, like our fixation pants, to ensure a tight fit and reduce the risk of leakage. ABENA's two-piece systems, consisting of a pair of fixation pants and an incontinence pad, are a cost-effective solution, as they use less raw materials, thereby having less impact on natural resources and the environment.



Shaped pad, ABENA San



Item number	Size	Dimensions	Absorption Rothwell	Quality level	Pcs per pk / pk per cll
1000021315	7/XXL	84 × 46 cm	2000 ml	Premium	26 / 4
1000021316	11/XXL	84 × 46 cm	3400 ml	Premium	16 / 4

ABENA San 7/XXL and 11/XXL are ABENA's range of anatomically shaped pads designed for bariatric use. The longer pad provides greater security for all body shapes. The high performance core with Top-Dry Technology ensures optimal skin care by keeping the exposed area comfortable and dry even after multiple voids. The fully breathable pads are fitted with standing leakage barriers combined with built-in cross barriers at the rear. When combined, these functionalities help protect the skin and maintain a natural pH skin balance. The Odor control system minimizes unpleasant smells. The wetness indicator with graduation scale guides the user and caretaker to determine the right absorption level and product changes. The dermatologically tested products have been given the Nordic Swan Ecolabel. These bladder protection pads should always be used with special fixation underwear, like our ABENA Fix Pants Super 4-6XL products, both to ensure a tight fit and to reduce the risk of leakage. ABENA's two-piece systems, consisting of a pair of fixation pants and an incontinence pad, are a cost-effective solution, as they use less raw materials, thereby having less impact on natural resources and the environment.



Shaped pad, ABENA San



Item number	Size	Dimensions	Absorption Rothwell	Quality level	Pcs per pk / pk per cll
1000021314	Special	73 × 37 cm	2000 ml	Premium	30 / 4

ABENA San Special is our anatomically shaped pad specifically designed for fecal incontinence. The pad can be used for both fecal and urinary incontinence and features a light green textile-like breathable backsheet foil for optimum comfort and discretion. The pad has effective leakage barriers with built-in pockets that keep the fecal contents inside. The high performance core with Top-Dry Technology ensures skin integrity and optimal skin care by keeping the exposed area comfortable and dry even after multiple voids. The fully breathable pads are fitted with standing leakage barriers combined with front and rear built-in cross barriers offering a 360° leakage protection. When combined, these functionalities help protect the skin and maintain a natural pH skin balance. The Odor control system minimizes unpleasant smells. The wetness indicator guides the user and caretaker to the product with the adequate absorption level and product changes. The dermatologically tested products have been given the Nordic Swan Ecolabel. ABENA San Special should always be used with special fixation underwear, like our ABENA fixation pants, to ensure a tight fit and reduce the risk of leakage.

