

Self-Cath® Plus Coudé Tip

This coudé tip intermittent catheter has fire-polished eyelets designed for comfort and a siliconized surface created for smooth insertion. The catheter has a hydrophilic coating that activates immediately upon exposure to water. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. **This catheter is not made with natural rubber latex.**



Self-Cath® Plus Tapered Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
4608	● 8 FR	30	16" length, funnel end	A4352
4610	● 10 FR	30	16" length, funnel end	A4352
4612	○ 12 FR	30	16" length, funnel end	A4352
4614	● 14 FR	30	16" length, funnel end	A4352



Self-Cath® Plus Olive Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
4808	● 8 FR	30	16" length, funnel end	A4352
4810	● 10 FR	30	16" length, funnel end	A4352
4812	○ 12 FR	30	16" length, funnel end	A4352
4814	● 14 FR	30	16" length, funnel end	A4352
4816	● 16 FR	30	16" length, funnel end	A4352

Self-Cath® Coudé Tip

Self-Cath coudé tip intermittent catheters have fire-polished eyelets designed for comfort, a siliconized surface created for smooth insertion. The blue Guide Stripe® along the top of the catheter allows monitoring of the coudé tip throughout the catheterization process. They are available in a variety of options. **This catheter is not made with natural rubber latex.**



Self-Cath® Olive Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
806	● 6 FR	30	16" length, funnel end	A4352
808	● 8 FR	30	16" length, funnel end	A4352
810	● 10 FR	30	16" length, funnel end	A4352
812	○ 12 FR	30	16" length, funnel end	A4352
814	● 14 FR	30	16" length, funnel end	A4352
816	● 16 FR	30	16" length, funnel end	A4352
818	● 18 FR	30	16" length, funnel end	A4352



Self-Cath® Tapered Coudé Tip with Guide Stripe®

Code	Size	Units	Description	HCPCS*
608	● 8 FR	30	16" length, funnel end	A4352
610	● 10 FR	30	16" length, funnel end	A4352
612	○ 12 FR	30	16" length, funnel end	A4352
614	● 14 FR	30	16" length, funnel end	A4352