

SpeediCath® Standard

The world's first ready-to-use hydrophilic-coated intermittent catheter

Easy to carry, store and use. The catheter is pre-hydrated and ready-to-use, eliminating wait time. The hydrophilic coating is smooth and uniform, designed to minimize urethral friction compared to uncoated catheters. The eyelets are polished for comfort. The uncoated Grip Zone™ area is created to offer control in handling. **This catheter is not made with natural rubber latex, DEHP or other phthalates.**

SpeediCath® Standard

SpeediCath catheter only. No insertion supplies.

**Straight Tip**

Code	Size	Units	Description	HCPCS*
28408	● 8 FR	30	14" length	A4351
28410	● 10 FR	30	14" length	A4351
28412	○ 12 FR	30	14" length	A4351
28414	● 14 FR	30	14" length	A4351
28416	● 16 FR	30	14" length	A4351
28418	● 18 FR	30	14" length	A4351

Female Straight Tip

28506	● 6 FR	30	6" length	A4351
28508	● 8 FR	30	6" length	A4351
28510	● 10 FR	30	6" length	A4351
28512	○ 12 FR	30	6" length	A4351
28514	● 14 FR	30	6" length	A4351
28516	● 16 FR	30	6" length	A4351

Pediatric Straight Tip

28706	● 6 FR	30	6" length	A4351
28708	● 8 FR	30	6" length	A4351
28710	● 10 FR	30	6" length	A4351

Boy Straight Tip

28606	● 6 FR	30	10" length	A4351
28608	● 8 FR	30	10" length	A4351
28610	● 10 FR	30	10" length	A4351
28612	○ 12 FR	30	10" length	A4351

*For important information regarding the use of HCPCS codes listed and code definitions see page 180 | For important safety information see page 8