



Cardinal Health™

Open Suction



Table of contents

Open Suction Products

Tips, valves and packaging	3
----------------------------	---

Single Suction Catheters

With chimney valve	4
With directional valve	5

Suction Catheter Dry Kits and Mini Trays

With chimney valve	5
Mini tray with chimney valve	6
Mini tray with directional valve	6

Suction Catheter Wet Trays

With sterile saline	7
With sterile water	7

Irrigating Solutions

Sterile saline	7
Sterile water	7

Mucus Traps

With vacuum breaker and filter	8
With vacuum breaker	8
Suction catheter with mucus trap	8

Specimen Containers

Lukens specimen container	9
Specimen trap	9
Luki disposable aspirating tubes	9

Tracheostomy Care

Tracheostomy tube holder	10
--------------------------	----

Open Suction Products

These products are available with different tip, valve and packaging options to meet each patient's needs. Every suction catheter has two staggered eyes to allow for efficient suctioning. Both valve options deliver easy fingertip control with virtually no residual vacuum to minimize trauma.



Tips

DeLee

- Standard tip for all pediatric and directional valve suction catheters
- Smooth tip to minimize trauma

Whistle

- Round, angled tip to provide efficient nontraumatic flow rate



Valves

Chimney valve

- Minimal residual vacuum in the open position

Directional valve

- Minimizes risk of cross-contamination



Suction catheter packaging

Straight

Coiled

Single Suction Catheters



31420



30888

31400

Cardinal Health™ Argyle™ Single Suction Catheter with Chimney Valve

- Available with the DeLee tip for pediatric sizes (5-8 Fr) and whistle tip for adult sizes (10-18 Fr). Also available with no control valve.
- 1-cm graduation marks on all catheters with control valves from 4-40 cm
- Sterile, single use
- Not made with natural rubber latex

Code	French size	Tip	Ship case
Straight-packed			
30500-	5 Fr, pediatric	DeLee	50
30688	6 Fr, pediatric	DeLee	50
30888	8 Fr, pediatric	DeLee	50
31000-	10 Fr	Whistle	50
31200-	12 Fr	Whistle	50
31400	14 Fr	Whistle	50
31445-	14 Fr, with Coudé tip	Whistle	50
31600	16 Fr	Whistle	50
31800	18 Fr	Whistle	50
Coil-packed			
30520	5 Fr, pediatric	DeLee	50
30620	6 Fr, pediatric	DeLee	50
30820	8 Fr, pediatric	DeLee	50
31020	10 Fr	Whistle	50
31220-	12 Fr	Whistle	50
31420	14 Fr	Whistle	50
31620-	16 Fr	Whistle	50
31820	18 Fr	Whistle	50
Straight-packed with straight connector with control valve			
33425	14 Fr	Whistle	50
Coil-packed with straight connector without control valve			
33000	10 Fr	Whistle	50
33400-	14 Fr	Whistle	50
33800	18 Fr	Whistle	50

Single Suction Catheters



Cardinal Health™ Argyle™ Single Suction Catheter with Directional Valve

- Directional valve with ISO standard colors to easily identify French sizes
- Sterile, single use
- Not made with natural rubber latex

Code	French size	Tip	Color	Ship case
Straight-packed				
141906	6 Fr, pediatric, graduated	DeLee	Light green	100
141908	8 Fr, pediatric, graduated	DeLee	Light blue	100
141900	10 Fr	DeLee	Black	100
141901	12 Fr	DeLee	White	100
141902	14 Fr	DeLee	Green	100
141903	16 Fr	DeLee	Orange	100
141904	18 Fr	DeLee	Red	100

Suction Catheter Dry Kits and Mini Trays



Cardinal Health™ Argyle™ Suction Catheter Kit with Chimney Valve

- Sterile, coil-packed suction catheter with chimney valve
- Two (2) walled blue nitrile exam gloves
- 1-cm graduation marks on all catheters with control valves from 4-40 cm
- Pop-up solution cup in sterile inner wrap
- Not made with natural rubber latex

Code	French size	Tip	Ship case
30579	5 Fr	DeLee	50
30677	6 Fr	DeLee	50
30877	8 Fr	DeLee	50
31077	10 Fr	DeLee	50
31079	10 Fr	Whistle	50
31279	12 Fr	Whistle	50
31479	14 Fr	Whistle	50
31679	16 Fr	Whistle	50
31879	18 Fr	Whistle	50

Suction Catheter Dry Kits and Mini Trays

37424



Cardinal Health™ Argyle™ Suction Catheter Mini Tray with Chimney Valve

- Sterile, coil-packed suction catheter with chimney valve
- 1-cm graduation marks on all catheters with control valves from 4-40 cm
- Two (2) walled blue nitrile exam gloves
- Mini tray for maximum efficiency
- Not made with natural rubber latex

Code	French size	Tip	Ship case
36524	5 Fr	DeLee	50
36626	6 Fr	DeLee	50
36826	8 Fr	DeLee	50
37024-	10 Fr	Whistle	50
37224	12 Fr	Whistle	50
37424	14 Fr	Whistle	50
37524	16 Fr	Whistle	50
37724	18 Fr	Whistle	50

14092



Cardinal Health™ Argyle™ Suction Catheter Mini Tray with Directional Valve

- Sterile, coil-packed suction catheter with directional valve
- Two (2) walled blue nitrile exam gloves
- Mini tray for maximum efficiency
- Not made with natural rubber latex

Code	French size	Tip	Ship case
140988	8 Fr, pediatric, graduated	DeLee	100
140980	10 Fr	DeLee	100
140982	14 Fr	DeLee	100

Suction Catheter Wet Trays




Cardinal Health™ Argyle™ Suction Catheter Tray with Solution

- Sterile, coil-packed suction catheter tray with chimney valve
- Two (2) walleded blue nitrile exam gloves
- 1-cm graduation marks on all catheters with control valves from 4-40 cm
- Not made with natural rubber latex

Code	French size	Tip	Ship case
Sterile saline — 100-mL bottles			
12062-	6 Fr	DeLee	24
12082	8 Fr	DeLee	24
12102-	10 Fr	Whistle	24
12142	14 Fr	Whistle	24
12153	14 Fr with drape	Whistle	24
12182-	18 Fr	Whistle	24
Sterile Water — 100-mL bottles			
10062	6 Fr	DeLee	24
10082	8 Fr	DeLee	24
10102	10 Fr	Whistle	24
10122	12 Fr	Whistle	24
10142	14 Fr	Whistle	24
12191	14 Fr with drape	Whistle	24
10182	18 Fr	Whistle	24

Irrigating Solutions



Cardinal Health™ Argyle™ Sterile Saline and Sterile Water

- Sterile water and saline for suction and irrigation
- Safety seal

Code	Description	Ship case
Sterile saline		
1020	3.38 oz (100 mL) Saline, (8) six packs	48
1022-	3.38 oz (100 mL) Saline, bulk	48
Sterile water		
10000-	3.38 oz (100 mL) Water, (8) six packs	48
1024-	3.38 oz (100 mL) Water, bulk	48

Mucus Traps



8888257469

Cardinal Health™ Argyle™ Mucus Traps with Vacuum Breaker and Filter

- For oral or wall/mechanical suction
- Advanced hydrophobic filter
- Whistle tip suction catheter with Slide-Tex™ finish
- Transparent 20-cc graduated mucus trap with removable cap
- Additional cap for transport to lab
- Sterile



Caution:
This product contains natural rubber latex which may cause allergic reactions.

Code	Size	Catheter French size	Ship case
8888257485	20 cc	6.5 Fr (2.2 mm)	50
8888257477	20 cc	8 Fr (2.7 mm)	50
8888257469	20 cc	10 Fr (3.3 mm)	50



8888257378

Cardinal Health™ Argyle™ Mucus Trap with Vacuum Breaker

- For wall/mechanical suction
- Directional valve
- Whistle tip suction catheter with Slide-Tex™ finish
- Transparent 20-cc graduated mucus trap with removable cap
- Additional cap for transport to lab
- Sterile
- Not made with natural rubber latex

Code	Size	Catheter French size	Ship case
8888257386	20 cc	6.5 Fr (2.2 mm)	50
8888257378	20 cc	8 Fr (2.7 mm)	50
8888257360	20 cc	10 Fr (3.3 mm)	50



8888257550

Cardinal Health™ Argyle™ Suction Catheter with Mucus Trap

- For oral suction
- Whistle tip catheter with Slide-Tex™ finish
- Transparent 20-cc graduated mucus trap with removable cap
- Suction tube with molded mouthpiece
- Sterile
- Not made with natural rubber latex

Code	Size	Catheter French size	Ship case
8888257550	20 cc	6.5 Fr (2.2 mm)	50
8888257535	20 cc	8 Fr (2.7 mm)	50
8888257527	20 cc	10 Fr (3.3 mm)	50

Specimen Containers



8888258608

Cardinal Health™ Argyle™ Lukens Specimen Container

- 14 Fr whistle-tip suction catheter with Slide-Tex™ finish
- 20-cc specimen trap graduated at 2-cc intervals
- Extra trap cap
- Suction tube with integral tapered connector
- Patient label
- Sterile
- Not made with natural rubber latex

Code	French size	Ship case
8888258608	14 Fr (4.7 mm)	50



8884724500

Cardinal Health™ Argyle™ Specimen Trap

- Screw-on cap
- Removable male connector
- Full 40- or 70-cc capacity with 1-cc graduation marks
- Patient label
- Sterile
- Not made with natural rubber latex
- DEHP-free

Code	Capacity	Ship case
8884724500	40 cc	50
8884724504	70 cc	50



C30200A

Cardinal Health™ Medi-Vac® Mucus Specimen Trap

- Removable male connector
- 40-cc capacity with 1-cc graduation marks
- Sterile
- Not made with natural rubber latex

Code	Capacity	Ship case
C30200A	40 cc	24



8886864304

Cardinal Health™ Argyle™ Luki Disposable Aspirating Tube

- Sterile
- Disposable



Caution:
This product contains natural rubber latex which may cause allergic reactions.

Code	Capacity	Tube size	Ship case
8886864304	8 cc	3.25 in.	48
8886864604	20 cc	6.25 in.	48

Tracheostomy Care

Cardinal Health™ Argyle™ Tracheostomy Tube Holder

- Easy to thread through tracheostomy plate
- Velcro® fasteners
- Not made with natural rubber latex
- Single use
- Non-sterile

Code	Description	Ship case
49000	Tracheostomy tube holder, Adult	100



