


Sports Medicine product catalog



actimovesport.com


Actimove

COACH
by Actimove
ZONAS
by Actimove
ELASTIKON
by Actimove

Delta-Cast®

Leukoplast®

essity



Forward

It's just not a direction.

It's an attitude.

It's about not resting on our laurels,

but leading the way.



Forward-thinking is what makes us a leader in health, hygiene, and medical solutions.

Our extensive portfolio of leading therapeutic brands reflects the spirit of innovation we bring to everything we do.

We have the courage to lead, to break barriers, and to always push for better.

Better Products.

Better Service.

Better Value.

Better Outcomes.

Our wide range of therapeutic solutions, clinical expertise and world-class brands **mean one partner to support your needs.**

Orthopedics: **Actimove®**

COACH
by Actimove

ZONAS
by Actimove

ELASTIKON
by Actimove

Delta-Cast®

Wound Care: **Leukoplast®**

Your advantage

Convenience

- One partner for a range of therapeutic solutions
- The right assortment tailored to your needs

Medical expertise

- On-going educational tools and support
- Quality products based on validated evidence for better outcomes
- Confidently recommend the right product

Our goal

Being your ally in better outcomes for you and those in your care. Now and for the future.



Contents

Sports Taping & Bandaging

Introduction	6
Portfolio overview	7
1.1 ELASTIC ADHESIVE TAPES	
ELASTIKON® by Actimove®	8
TensoSport®	8
Tensoplast®	9
Tensoplast® LF	9
1.2 LIGHTWEIGHT ELASTIC ADHESIVE TAPES	
Lightplast® Pro	10
1.3 NON-ELASTIC ADHESIVE STRAPPING TAPES	
Leukotape® P	10
COACH® by Actimove®	11
ZONAS® by Actimove®	11
1.4 ADHESIVE PROTECTION SPRAY	
Tensospray®	11
1.5 KINESIOLOGY TAPES	
Leukotape® K	12
1.6 NON-ADHESIVE SUPPORT BANDAGES	
BSN Elastic Bandage	12
1.7 COHESIVE LIGHT SUPPORT BANDAGES	
Co-Plus® LF	13
1.8 TUBULAR SUPPORT BANDAGES	
Tensogrip®	13

Fracture Management

Introduction	14
Portfolio overview	15
2.1 CAST TAPE	
Delta-Lite® Plus	16
Delta-Lite® Conformable	16
2.2 PADDINGS	
Delta-Rol® Synthetic Cast	17
Specialist® Micro Pleat Cotton Blend Cast	17
Sof-Rol® Absorbent Rayon Cast	18
Delta-Dry® Water Resistant Padding	18
AquaCast® Liner	19
2.3 STOCKINETTE	
Delta-Net® Orthopedic Synthetic	20
Specialist® Orthopedic Cotton	20
Terry-Net™ S Synthetic Terry Cloth	20
Terry-Net™ Thumb Spica	21
2.4 CAST ACCESSORIES	
Delta-Net® Adhesive Felt Padding	21
Sifoam™ Interface Padding	21

Fracture Management (cont.)

Delta Terry-Net™ Adhesive Fleece Edger	22
Delta Terry-Net™ Adhesive Hook & Stretch Loop Closures	22
Zip Stick	22
American Orthopaedic™	23
Essity Pail	23
Saw Stop® DE-FLEX™	23
Stainless Steel Cast Spreader	24
FCT Scissors	24
Esmarch Scissors	24
2.5 SPLINTING	
ORTHO-GLASS®	25

Orthopedic Soft Goods

Introduction	26
Portfolio overview	27
3.1 Actimove® UPPER EXTREMITY SUPPORTS	
Actimove® Shoulder Support	28
Actimove® Mitella Comfort Arm Sling	28
Actimove® Elbow Strap	29
Actimove® Elbow Support	29
Actimove® Episport® Elbow Clasp	
Adjustable	29
Actimove® Wrist Stabilizer	30
Actimove® Thumb Stabilizer	30
Actimove® Alufoam	31
Actimove® Rhizo Forte Thumb	31
3.2 Actimove® LOWER EXTREMITY SUPPORTS	
Actimove® Knee Brace	32
Actimove® Knee Brace	32
Actimove® Knee Brace	33
Actimove® Knee Stabilizer	33
Actimove® Knee Support	34
Actimove® Post-Op ROM	34
Actimove® Actimove® Thigh Support	35
Actimove® Thigh Support	36
Actimove® PowerMotion Thigh Support	36
Actimove® PowerMotion Calf Support	37
Actimove® Talocast Ankle Brace	37
Actimove® Standard Walker Low / High	38
Actimove® Standard Air Walker	
Low / High	38
Actimove® Closed Shell Air Walker	
Low / High	38
Actimove® Ankle Stabilizer	39
Actimove® Patella Strap	40

Wound Care

Introduction	42
Portfolio overview	43
4.1 STRETCHABLE LIGHTWEIGHT FABRIC DRESSING	
Leukoplast® elastic	44
4.2 TRANSPARENT DRESSING WITH ABSORBENT PAD	
Leukomed® T Plus	45
4.3 WIDE AREA FIXATION	
Leukoplast® Hypafix®	45
4.4 ELASTIC GAUZE BANDAGE	
Elastomull®	46
4.5 HYDROPHOBIC MICROBE BINDING FILM DRESSING	
Leukomed® Sorbact®	46
4.6 CONFORMING, ADHESIVE FIXATION BANDAGE	
Cover-Roll® Stretch	47
4.7 WOUND & EYE IRRIGATION	
SplashCap®	47
EyeCap®	47



Athletic Taping & Bandaging

Essity has a full line of quality sports tapes to help support your athletes.

Portfolio overview

ELASTIC ADHESIVE TAPES



ELASTIKON®
by Actimove®



TensoSport®



Tensoplast®



Tensoplast® LF

LIGHTWEIGHT ELASTIC ADHESIVE TAPES



Lightplast® Pro

NON-ELASTIC ADHESIVE STRAPPING TAPES



Leukotape® P



COACH®
by Actimove®



ZONAS®
by Actimove®

ADHESIVE PROTECTION SPRAY



Tensospray®

KINESIOLOGY TAPES



Leukotape® K

NON-ADHESIVE SUPPORT BANDAGES



BSN Elastic Bandage

COHESIVE BANDAGES



Co-Plus® LF

TUBULAR SUPPORT BANDAGES



Tensogrip®

ELASTIKON® by Actimove®

Soft cotton elastic tape, centerline for ease of application

ELASTIKON
by Actimove

- Reliable compression for support of sprains, strains and muscle injuries
- Soft cotton elastic tape conforms to joints such as elbows, knees and shoulders
- Allows skin to breathe and moisture to pass through
- Treatment of conditions which require limitation of movement, e.g. sprains, sports injuries
- Prophylactic use in active sport
- Retention of dressings and devices, e.g. post-orthopedic and general surgery or over a muscular wound



REF/BNR	Measurements	Color	UOM
5170	2 in x 2.5 yds	Cream	24 rolls/case (Speed Pack)
5171	3 in x 2.5 yds	Cream	16 rolls/case (Speed Pack)
5172	1 in x 2.5 yds	Cream	144 rolls/case (12 rolls/box, 12 box/case)
5174	2 in x 2.5 yds	Cream	72 rolls/case (6 rolls/box, 12 box/case)
5175	3 in x 2.5 yds	Cream	48 rolls/case (4 rolls/box, 12 box/case)
5177	4 in x 2.5 yds	Cream	36 rolls/case (6 rolls/box, 6 box/case)

This Product Contains Dry Natural Rubber.

TensoSport®

Soft cotton elastic tape, centerline for ease of application

TensoSport® is an elastic adhesive bandage suitable for medium support or compression. Excellent stretch and regain provides firm, controlled compression and support. Fluffy non-fray edges prevent the bandage from rucking or marking.

- High adhesive properties
- Provides strong compression and long-lasting support
- Water repellent

TensoSport Elastic Adhesive Bandage – Boxed Packs

REF/BNR	Measurements	Color	UOM
7205015	2 in x 5 yds	White	36 rolls/case
7205016	3 in x 5 yds	White	24 rolls/case

Caution: This product contains natural rubber latex which may cause allergic reactions



Tensoplast®

Soft cotton elastic tape

Tensoplast® is an elastic adhesive bandage suitable for medium support or compression. Excellent stretch and regain provides firm, controlled compression and support. Fluffy non-fray edges prevent the bandage from rucking or marking.

- High adhesive properties
- Provides strong compression and long-lasting support
- Water repellent



Tensoplast® Elastic Adhesive Bandage – Individually Boxed

REF/BNR	Measurements	Color	UOM
02598002	1 in x 5 yds	Beige	36 rolls/case
02599002	2 in x 5 yds	Beige	36 rolls/case
02600002	3 in x 5 yds	Beige	36 rolls/case
02601002	4 in x 5 yds	Beige	36 rolls/case
02602002	6 in x 5 yds	Beige	12 rolls/case
02593002	1 in x 5 yds	White	36 rolls/case
02594002	2 in x 5 yds	White	36 rolls/case
02595002	3 in x 5 yds	White	36 rolls/case
02596002	4 in x 5 yds	White	36 rolls/case
02597002	6 in x 5 yds	White	12 rolls/case

Tensoplast® Elastic Adhesive Bandage – Bulk Packs

REF/BNR	Measurements	Color	UOM
04410001	1 in x 5 yds	Beige	32 rolls/case
04412001	2 in x 5 yds	Beige	24 rolls/case
04413001	3 in x 5 yds	Beige	16 rolls/case
04414001	4 in x 5 yds	Beige	16 rolls/case

Caution: This product contains natural rubber latex which may cause allergic reactions

Tensoplast® LF

Tensoplast® LF is not made with natural rubber latex and assures the quality you trust. It provides flexible, strong support and compression thanks to its porously high-quality substrate. Its excellent stretch and adhesive molds and holds to any part of the body, keeping the bandages firmly in place even throughout lengthy treatment. Trusted patient comfort and compliance for those sensitive to latex make Tensoplast® LF the optimal choice for a wide range of indications.

- Permeable material
- High-quality elastic substrate
- Adhesive coating
- Yellow center line
- Not made with natural rubber latex



REF/BNR	Measurements	Color	UOM
7623200	2 in x 4.9 yds	Cream	12 rolls/case
7623201	3 in x 4.9 yds	Cream	12 rolls/case
7623202	4 in x 4.9 yds	Cream	12 rolls/case

1.2 LIGHTWEIGHT ELASTIC ADHESIVE TAPES

Lightplast® Pro

Tearable Adhesive Stretch Bandage

Lightplast® Pro is a tearable, elastic adhesive bandage, suitable for light support to compression. The unique fabric backing helps to conform to body movement. Ideal for all-purpose taping.

- Hand-tearable
- Excellent adhesion
- Light and close fitting

REF/BNR	Measurements	Color	UOM
76952	1 in x 5 yds	White	48 rolls/case
76953	1.5 in x 5 yds	White	32 rolls/case
76954	2 in x 5 yds	White	24 rolls/case
76955	3 in x 5 yds	White	16 rolls/case
76956	2 in x 7.5 yds	White	24 rolls/case
76957	3 in x 7.5 yds	White	16 rolls/case
72880	2 in x 5 yds	Black	24 rolls/case
72881	3 in x 5 yds	Black	16 rolls/case

This Product Contains Dry Natural Rubber.



1.3 NON-ELASTIC ADHESIVE STRAPPING TAPES

Leukotape® P

Rigid Strapping Tape

Leukotape® P is a high strength rigid strapping tape with a very strong zinc oxide adhesive. It is porous yet offers strong support for sprains and strains.

- Porous
- Easy unwinding tension
- High strength adhesive
- Ideal for patellofemoral taping

REF/BNR	Measurements	Color	UOM
76168	1.5 in x 15 yds	Beige	30 rolls/case

This Product Contains Dry Natural Rubber.



COACH® by Actimove®

Rigid Cloth Porous Adhesive Strapping Tape

- Support and conformity helps protect against sprains and injury
- High tensile strength
- Adheres securely yet tears easily by hand
- 100% cotton back cloth allows perspiration to escape
- Smooth and consistent unwind



REF/BNR	Measurements	Color	UOM
5188	1.5 in x 15 yds	White	32 rolls/case (Speed Pack)

This Product Contains Dry Natural Rubber.

ZONAS® by Actimove®

Rigid Cloth Porous Adhesive Strapping Tape

- Prophylactic use in active sport
- Ideal for conditions requiring immobilization and support
- Highest tensile strength
- Adheres securely yet tears by hand
- 100% cotton back cloth allows perspiration to escape
- Smooth and consistent unwind



REF/BNR	Measurements	Color	UOM
5104	1 in x 10 yds	White	144 rolls/case (12 rolls/box, 12 box/case))
5190	1.5 in x 15 yds	White	32 rolls/case (Speed Pack)

This Product Contains Dry Natural Rubber.

1.4 ADHESIVE PROTECTION SPRAY

Tensospray®

Adhesive Protection Spray

Tensospray® is an adhesive spray which leaves a thin protective film over the surface of the skin. This provides an ideal layer for adhesive bandages and tapes, promoting excellent traction and reducing skin trauma on removal. The product is not affected by disinfectant.

- Aerosol presentation makes for quick and easy application
- Provides a copolymer layer which improves the effectiveness of the strapping
- Protects skin and encloses hairs – aids comfort and reduces trauma of bandage removal



REF/BNR	Measurements	UOM
7160200	300ml/can	12 cans/case

Leukotape® K

Kinesiology Therapeutic Adhesive Tape

The wave pattern adhesive of **Leukotape® K** allows air to circulate. It is stretchable and works ideally with the skin's natural elasticity.

- High wearing comfort and skin friendly, with a formulation not made with natural rubber latex
- Stays in place up to 9 days, even after showering or sweating
- Residue free removal
- Dermatologically tested



REF/BNR	Measurements	Color	UOM / Pack
7297810	1 in x 5.5 yds	● Tan	5 rolls/pack
7297811	2 in x 5.5 yds		
7297812	3 in x 5.5 yds		
7297815	1 in x 5.5 yds	● Red	
7297816	2 in x 5.5 yds		
7297817	3 in x 5.5 yds		
7297820	1 in x 5.5 yds	● Blue	
7297821	2 in x 5.5 yds		
7297822	3 in x 5.5 yds		
7297823	2 in x 5.5 yds	● Black	
7297824	2 in x 5.5 yds	● Lt Blue	
7297825	2 in x 5.5 yds	● Pink	

1.6 NON-ADHESIVE SUPPORT BANDAGES

BSN Elastic Bandage

Non-Adhesive Bandage

- Adjustable and constructed of breathable elastic
- Hook and loop closure
- No metal clips
- Not made with natural rubber latex
- Individually wrapped



REF/BNR	Measurements	Color	UOM / Pack
1037023	2 in x 5 yds	● Beige	10 rolls/pack, 5 packs/case
1037033	3 in x 5 yds		
1037043	4 in x 5 yds		
1037053	6 in x 5 yds		

1.7 COHESIVE LIGHT SUPPORT BANDAGES

Co-Plus® LF

Cohesive Bandage

- Self-adhesive
- Non-slip support
- Hand-tearable in both directions
- End-of-bandage indicator for ease of use
- Not made with natural rubber latex
- Available in beige and color pack: neon pink, light blue, purple, neon green and red



REF/BNR	Measurements	Color	UOM
7210014	1 in x 5 yds	● Beige	60 rolls/box
7210016	2 in x 5 yds		36 rolls/pack
7210017	3 in x 5 yds		24 rolls/pack
7210018	4 in x 5 yds		18 rolls/pack
7210019	6 in x 5 yds		12 rolls/pack
7210023	2 in x 5 yds	Assorted	36 rolls/pack
7210020	3 in x 5 yds		24 rolls/pack
7210021	4 in x 5 yds		18 rolls/pack

1.8 TUBULAR SUPPORT BANDAGES

Tensogrip®

Easy to apply and help control edema

- Ideal for the light support of strains, sprains and soft tissue injuries
- Provides light support, light compression and dressing retention
- Constructed from a cotton/polyamide blend with natural rubber latex yarn



REF/BNR	Measurements	Color	UOM
7579	A - 1.75 in x 11 yds	White	1/roll
7580	B - 2.5 in x 11 yds	White	1/roll
7581	C - 3 in x 11 yds	White	1/roll
7582	D - 3.5 in x 11 yds	White	1/roll
7583	E - 3.5 in x 11 yds	White	1/roll
7584	F - 4 in x 11 yds	White	1/roll
7585	G - 4.5 in x 11 yds	White	1/roll
7586	J - 7 in x 11 yds	White	1/roll
7587	K - 8.5 in x 11 yds	White	1/roll



Caution: This product contains natural rubber latex which may cause allergic reactions



Fracture Management

As the global market leader in Fracture Management, drawing on a strong brand heritage, Delta-Cast® has a comprehensive and innovative portfolio to rely on.

Portfolio overview

CSAT TAPE



Delta-Lite®
Plus

Delta-Lite®
Conformable

PADDING



Delta-Rol®

Specialist®

Sof-Rol®

Delta-Dry®

AquaCast®
Liner

STOCKINETTE



Delta-Net®



Specialist®



Terry-Net™



Terry-Net™
Thumb Spica

CAST ACCESSORIES



Delta-Net®



Sifoam™



Delta Terry-Net™



Delta Terry-Net™



Zip Stick



American
Orthopaedic™



Essity® Pail



Saw Stop®
DE-FLEX™



Stainless Steel
Cast Spreader



FCT Scissors



Esmarch
Scissors

SPLINTING



ORTHO-GLASS®

Delta-Lite® Plus

Fiberglass Cast Tape

Delta-Lite® Plus provides extra rigid casts with smooth surfaces and edges. The advanced water activated resin system provides low tack properties allowing for ease of application.

- Smooth surface
- Rigid support
- Good conformability
- Excellent lamination
- Consistent and reliable working and set time
- Available in 11 colors and mixed packs
- Not made with natural rubber latex

REF/BNR	Measurements	Color	UOM
7345800	1 in × 4 yds	White	10 rolls/box
7345801	2 in × 4 yds	White	10 rolls/box
7345802	3 in × 4 yds	White	10 rolls/box
7345803	4 in × 4 yds	White	10 rolls/box
7345804	5 in × 4 yds	White	10 rolls/box



Delta-Cast® Conformable

Polyester Cast Tape

Delta-Cast® Conformable is an easy-to-use polyester cast tape that produces versatile, lightweight and strong casts with high conformability. Ideal for applications where patient fit and comfort are critical.

- High conformability
- Smooth finish and soft edges
- Rigid immobilization
- High radiolucency
- Excellent lamination between layers
- Excellent for FCT casting or removable casts
- colors and mixed packs
- Not made with natural rubber latex

REF/BNR	Measurements	Color	UOM
6601	1 in × 2 yds	White	10 rolls/box
6602	2 in × 4 yds	White	10 rolls/box
6603	3 in × 4 yds	White	10 rolls/box
6604	4 in × 4 yds	White	10 rolls/box
6605	5 in × 4 yds	White	10 rolls/box



Delta-Rol®

Synthetic Cast Padding

Delta-Rol® is made of breathable synthetic material. Its conformable stretch allows narrow widths fitting around small anatomies without cutting and tearing. Also available in sterile version.

- Conformable stretch
- Low linting
- Synthetic fibers
- Available in sterile version
- Not made with natural rubber latex

REF/BNR	Measurements	Color	UOM
6882	2 in × 4 yds	White	72 rolls/box
6883	3 in × 4 yds	White	36 rolls/box
6884	4 in × 4 yds	White	36 rolls/box
6886	6 in × 4 yds	White	18 rolls/box



Specialist®

Cotton Blend Cast Padding

Specialist® Cotton Blend Cast Padding has a unique micropleated design for excellent conformability, protection and ease of application without tucking or tearing.

- Micropleated construction provides excellent conformability
- Highly absorbent
- Cotton rayon blend
- Easy application

REF/BNR	Measurements	Color	UOM
9062	2 in × 4 yds	White	144 rolls/box
9043	3 in × 4 yds	White	72 rolls/box
9044	4 in × 4 yds	White	72 rolls/box
9046	6 in × 4 yds	White	36 rolls/box



Sof-Rol®

Rayon Absorbent Cast Padding

Sof-Rol® is a highly absorbent rayon material. This soft natural material creates a comfortable undercast surface for either synthetic or natural casting materials.

- Highly absorbent
- Easy application
- 100% viscose rayon
- Not made with natural rubber latex

REF/BNR	Measurements	Color	UOM
9052	2 in × 4 yds	White	144 rolls/box
9033	3 in × 4 yds	White	72 rolls/box
9034	4 in × 4 yds	White	72 rolls/box
9036	6 in × 4 yds	White	36 rolls/box



Delta-Dry®

Water Resistant Padding

The patented **Delta-Dry®** material offers water resistance. The open knit structure of the fibers accelerate drainage after submerging in water. The patient's body temperature will then evaporate residual moisture quickly through the airy mesh (high MVTR). In this way, cast and skin dry quickly and skin feels dry in about 90 minutes after daily hygiene.

- Allows patients to shower, swim, bathe and wash hands
- Helps reduce odor under the cast
- Excellent conformability to body contours
- Provides comfortable, fast-drying casts
- Not made with natural rubber latex

REF/BNR	Measurements	Color	UOM
7344300	2 in × 2.6 yds	White	12 rolls/box
7344301	3 in × 2.6 yds	White	12 rolls/box
7344302	4 in × 2.6 yds	White	12 rolls/box



AquaCast® Liner

Waterproof Cast Padding

AquaCast® Liner is a comfortable waterproof, breathable cast padding for use under waterproof casting materials. It is moisture vapor permeable and made with pillows and channels to enhance comfort and facilitate drainage.

- Quick drying of cast allows patients to bathe, swim and exercise
- Decrease odor and enhance hygiene with the ability to thoroughly wash and rinse skin
- Supports favorable skin condition following cast removal
- Less itching due to smooth material surfaces
- Less recasting due to wet or damaged casts
- Not made with natural rubber latex



REF/BNR	Measurements	Color	UOM
Aquacast® Liner Rolls			
ACL-2-S	2 in x 1.8yds	White	12 roll/case
ACL-3-S	3 in x 1.8yds	White	12 roll/case
ACL-4-S	4 in x 1.8yds	White	12 roll/case
ACL-6-S	6 in x 1.8yds	White	6 roll/case
Aquacast® Form Fit Liners			
FFS	Form Fit Liner - Small	White	4 units / bag
FFM	Form Fit Liner - Medium	White	4 units / bag
FFL	Form Fit Liner - Large	White	4 units / bag
Aquacast® Liner Starter Pack and Cast Kits			
ACLSP	Starter Kit	White	1 kit
2ACL-2-S	Cast Kit 2": 2 rolls 2", 1 roll SawStop	White	1 kit
3ACL-2-S	Cast Kit 2": 3 rolls 2", 1 roll SawStop	White	1 kit
2ACL-3-S	Cast Kit 3": 2 rolls 3", 1 roll SawStop	White	1 kit
3ACL-3-S	Cast Kit 3": 3 rolls 3", 1 roll SawStop	White	1 kit
2ACL-4-S	Cast Kit 4": 2 rolls 4", 1 roll SawStop	White	1 kit
3ACL-4-S	Cast Kit 4": 3 rolls 4", 1 roll SawStop	White	1 kit

The AquaCast® Liner Starter Pack includes 2 rolls of 2", 3", 4" AquaCast® Liner plus a full-length roll of SawStop® DE-FLEX® Protective Strip.

Delta-Net®

Orthopedic Synthetic Stockinette

- Soft and conformable
- Used for Casting and under bracing

REF/BNR	Measurements	Color	UOM
6861	1 in × 25yds	White	2 rolls/case
6862	2 in × 25yds	White	2 rolls/case
6863	3 in × 25yds	White	2 rolls/case
6864	4 in × 25yds	White	2 rolls/case
6866	6 in × 25yds	White	1 rolls/case
6868	8 in × 25yds	White	1 rolls/case
6810	10 in × 25yds	White	1 rolls/case



Specialist®

Orthopedic Cotton Stockinette

- Soft and conformable
- Used for Casting and under bracing

REF/BNR	Measurements	Color	UOM
9071	1 in × 25yds	White	1 rolls/case
9072	2 in × 25yds	White	9 rolls/case
9073	3 in × 25yds	White	8 rolls/case
9074	4 in × 25yds	White	8 rolls/case
9076	6 in × 25yds	White	6 rolls/case



Terry-Net™ S

Synthetic Terry Cloth Stockinette

- Soft and conformable
- Terry-net cushion
- High stretch
- Used for Casting and under bracing

REF/BNR	Measurements	Color	UOM
53012	2 in × 15yds	White	2 rolls/case
53013	3 in × 15yds	White	2 rolls/case
53014	4 in × 15yds	White	1 rolls/case
53015	5 in × 15yds	White	1 rolls/case
53016	6 in × 15yds	White	1 rolls/case



Terry-Net™ Thumb Spica

Precut Synthetic Terry Cloth Stockinette

- Soft and conformable
- Terry-net cushion
- High stretch
- Used for Casting and under bracing

REF/BNR	Measurements	Color	UOM
53120	2 in × 16yds (short arm)	White	10 box/case
53121	2 in × 19yds (long arm)	White	10 box/case
53101	3 in × 16yds (short arm)	White	10 box/case
53102	3 in × 19yds (long arm)	White	10 box/case



2.4 CAST ACCESSORIES

Delta Terry-Net™

Adhesive Felt Liner

Delta Terry-Net™ Adhesive Felt Padding is an adhesive orthopedic felt, used for added patient comfort. Its adhesive backing makes attaching the padding simple and easy. It may also be used as a cast edge padding or to provide protection for bony prominences

- Soft orthopedic felt
- Adhesive back

REF/BNR	Measurements	Color	UOM
54012	23 in x 39 in	White	1 roll/case



Sifoam™

Adhesive Back Padding

- Offloads pressure around bony prominences
- Helps prevent shearing
- Is compressable
- Unloads weight

REF/BNR	Measurements	Color	UOM
7279900	8 in x 9 in	Black	10 units/box



Delta Terry-Net™

Adhesive Fleece Edger

Delta Terry-Net™ Adhesive Fleece Edger features a high quality adhesive for secure placement. Its unique stretch allows the user to pad easily around corners and contours, ensuring smooth cast edges.

- Highly conformable
- Unique stretch
- Secure adhesive
- Works great with thermoplastics

REF/BNR	Measurements	Color	UOM
56012	1.25 in × 15 yds	White	1 roll/case
58012	1.25 in × 15 yds	Blue	1 roll/case



Delta Terry-Net™

Adhesive Hook & Stretch Loop Closures

This easy-to-use product simplifies the FCT technique. It provides closure without creating the discomfort caused by traditional loops' rough edges. Its high quality adhesive supports secure placement and the simple dispensing box provides convenient storage.

- Soft stretch loop
- High quality adhesive
- Dispensing box

REF/BNR	Measurements	Color	UOM
Adhesive Hook			
52100	1 in × 10 yds	White	1 roll/case
52200	2 in × 10 yds	White	1 roll/case
Stretch Loop			
52101	1 in × 15 yds	White	1 roll/case
52201	2 in × 15 yds	White	1 roll/case



Zip Stick

Cast Removal Aid

This convenient cast aid provides additional patient protection during cast removal. The stick slides under the cast and provides a smooth path for scissors or cast saw. Standard and metric rule adds functionality.

REF/BNR	Color	UOM
294304	White	1 ea



American Orthopaedic™

Cast Cutter

The American Orthopaedic™ Cast Cutter provides precisely what is needed when needed: power, speed, quiet operation, durability, ease of operation and a balance of weight and power. Includes one 2.5" diameter multi-purpose blade and wrench.

- Balance of weight and power helps reduce user fatigue
- Internally lubricated for longer service
- Designed for rugged use
- Switch is designed and placed so that it will not break if dropped
- Hex drive provides six (6) blade rotations extending blade life by 50%

REF/BNR	Color	UOM
0295-200	Blue	1/box



Casting Pail

This sturdy pail can be used for water supply for cast activation.

REF/BNR	Description	UOM
7204626	Essity Logo Casting Pail	1 / each



Saw Stop® DE-FLEX™

These protective strips are tough, waterproof, flexible ribbons that add an added layer of protection from a cast saw.

- Adhesive Backing for easy application

REF/BNR	Measurements	Color	UOM
SS10	1.25 in × 10 ft	Blue	1/box
SS30	1.25 in × 10 ft	Blue	3/box
SS100	1.25 in × 10 ft	Blue	10/box
SSW10	1.25 in × 10 ft	White	1/box
SSW30	1.25 in × 10 ft	White	3/box
SSW100	1.25 in × 10 ft	White	10/box



Casting Instruments

Casting instruments from Essity are effective for cutting casting material. They are handcrafted with high quality stainless steel to provide precision performance.

Stainless Steel Cast Spreader

- Spring loaded
- Excellent to open cut casts in cast removal

REF/BNR	Size	UOM
21-0136	9 in	1 unit



FCT Scissors

- Sharp, blunt tip scissors
- Clean cut
- Excellent for cast trimming

REF/BNR	Size	UOM
28230	8 in	1 unit



Esmarch Scissors

- Sharp, blunt tip scissors
- Excellent to cut bandages

REF/BNR	Size	UOM
7204407	8 in	1 unit



ORTHO-GLASS®

Fiberglass Splinting System

ORTHO-GLASS® is a cost-effective splinting system of choice for every phase of immobilization. Features Interlocking Performance Technology (ILP). ILP eliminates the risk of delamination, ensuring a strong, durable splint.

- Eliminates risk of delamination
- No wrong side application
- Easy stretch padding covers all edges
- Roll form eliminates waste
- Full educational support system – workshops, in-services, manuals, wall-charts, videos



ORTHO-GLASS® Splinting System

REF/BNR	Measurements	UOM
OG-1L2	1 in x 15 ft	2 rolls / case
OG-2L2	2 in x 15 ft	2 rolls / case
OG-3L2	3 in x 15 ft	2 rolls / case
OG-4L2	4 in x 15 ft	2 rolls / case
OG-5L2	5 in x 15 ft	2 rolls / case
OG-6L2	6 in x 15 ft	2 rolls / case
66000687	8 in x 7.5 ft	1 roll

ORTHO-GLASS® Pre-Cut Splints

REF/BNR	Measurements	UOM
OG-2PC	2 in x 12 in	10 splints / box
OG-3PC	3 in x 12 in	10 splints / box
OG-335PC	3 in x 35 in	5 splints / box
OG-4PC	4 in x 15 in	10 splints / box
OG-430PC	4 in x 30 in	5 splints / box
OG-5PC	5 in x 30 in	5 splints / box
OG-545PC	5 in x 45 in	5 splints / box
OG-610PC	6 in x 10 in	10 splints / box
OG-6PC	6 in x 30 in	5 splints / box
OG-645PC	6 in x 45 in	5 splints / box



Orthopedic Soft Goods

Actimove® is proud to provide medical devices in a broad range of excellent quality products that are reliable, thoroughly tested and tailored to meet the individual patient's needs.

Whether customers doing exercises, recovering from an injury, overcoming daily arthritis or movement-pain related challenges, they can rely on Actimove® to help them with pain relief and provide the support they need.

Portfolio overview

Actimove® UPPER EXTREMITY SUPPORTS



Actimove®
Shoulder
Support



Actimove®
Mitella Comfort
Arm Sling



Actimove®
Elbow
Strap



Actimove®
Elbow
Support



Actimove®
Episport®
Epicondylitis
Clasp



Actimove®
Wrist
Stabilizer



Actimove®
Thumb
Stabilizer



Actimove®
Alufoam



Actimove®
Rhizo Forte
Thumb

Actimove® LOWER EXTREMITY SUPPORTS



Actimove®
Knee
Brace



Actimove®
Knee
Brace



Actimove®
Knee
Brace



Actimove®
Stabilizer



Actimove®
Support



Actimove®
GenuMotion
Knee Support



Actimove®
Post Op Rom



Actimove®
Thigh Support



Actimove®
PowerMotion
Thigh Support



Actimove®
PowerMotion
Calf Support



Actimove®
Standard
Walker
Low / High



Actimove®
Standard
Air Walker
Low / High



Actimove®
Closed Shell
Air Walker
Low / High



Actimove®
Talocast
Ankle Brace



Actimove®
Ankle Stabilizer



Actimove®
Patella Strap

Actimove® Shoulder Support

Extra Pocket for Optional Hot/Cold Pack

- Pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Extra pain relief when optional hot/cold therapy pack is placed in build-in pocket
- Supports and protects the shoulder while allowing flexibility
- Adjustable chest and arm closure for a custom fit
- Can be worn left or right



	REF/BNR	Size	Measurements	Color	UOM
	7557310	S	31.5 in - 41.5 in	Black	1 ea
	7557311	M	41.5 in - 51 in		
	7557312	L	51 in - 61 in		

Actimove® Mitella Comfort

Comfort Arm Sling

- Easy to apply, adjust or remove with one hand thanks to the hook and loop Y-Tabs and open and-close pocket design
- For good patient wearing comfort, reduced neck pressure though the 5.5 cm wide and stable sling replacing a separate neck pad
- Helps to relieve pain by stabilizing the shoulder joint and keeping the flexed elbow in a pocket which limits the range of motion



	REF/BNR	Size	Measurements	Color	UOM
	7281961	S	< 37 cm	Blue	1 ea
	7281962	M	37 cm - 40 cm		
	7281923	L	40 cm - 43 cm		
	7281964	XL	> 43 cm		

Actimove® Elbow Strap

Hot/Cold Pack

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Extra pain relief supported by a removable therapeutic hot/cold pack
- Silicone dot lining for secure fit
- Can be worn left or right

	REF/BNR	Size	Measurements	Color	UOM
	7574230	Universal	11 in-15 - .75 in	Black	1 ea



Actimove® Elbow Support

Adjustable

- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Quick and easy to fit with slip-on design and individually adjustable straps
- Can be worn left or right

	REF/BNR	Size	Measurements	Color	UOM
	7561730	Universal	7.5 in - 12.25	Black	1 ea



Actimove® EpiSport®

- Provides support for Epicondylitis, Tennis Elbow, and Muscle Strain
- Dual action compression works simultaneously to treat lateral and medial epicondylitis (forehand and backhand tennis elbow)
- Comfortable shock absorbing liner helps wick moisture away

	REF/BNR	Size	Measurements	Color	UOM
	7347000	XS	7.75 in - 9 in	Black	1 ea
	7347001	S	9 in - 9.75 in		
	7347002	M	9.75 in - 11 in		
	7347003	L	11 in - 11.75 in		
	7347004	XL	11.75 in - 13.5 in		



Actimove® Wrist Stabilizer

Removable Metal Stay, Left or Right

- Anatomically shaped stay for secure stabilization, managing sprains and strains
- Pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Stabilizes wrist through an anatomically shaped removable stay for left and right hand application
- Quick and easy to fit with slip-on design and individually adjustable straps
- Special comfort in a premium soft material and seamless design
- Can be easily worn under clothing due to low profile
- Can be worn left or right



COOLMAX®
PRO

	REF/BNR	Size	Measurements	Color	UOM
	7572931	S	5 in - 6.25 in	Black	1 ea
	7572932	M	6.25 in - 7.5 in	Black	
	7572933	L	7.5 in - 9 in	Black	

Actimove® Thumb Stabilizer

Extra Stays

- Anatomically shaped stays, managing thumb injuries, keeping finger dexterity
- Pain relief and healing are supported through balanced warmth and medical compression
- Comfort and good breathability through neoprene-free performance material
- Stabilization of the thumb through two integrated anatomically shaped stays
- Comfortable to wear because of premium quality soft material
- Firm and adjustable wrist support from powerful elastic strap
- Full finger movement and hand dexterity from special open design



	REF/BNR	Size	Measurements	Color	UOM
	7563630	S/M	5.5 in - 7.5 in	Black	1 ea
	7563631	L/XL	7.5 in - 9 in		
	7563632	S/M	5.5 in - 7.5 in	Beige	
	7563633	L/XL	7.5 in - 9 in		

Actimove® Alufoam

Maleable Aluminum Cushioned Finger Strips

- Excellent finger protection
- Lightweight adjustable, trimmable finger splint

REF/BNR	Size	Measurements	UOM
7295600	S	.5 in - 18 in	10/box
7295601	M	.75 in - 7.5 in	
7295602	L	1 in - 9 in	



Actimove® Rhizo Forte

Thumb Orthosis

- Improved design for optimized fit and comfort around the thumb
- For everyday wear, as wrist and fingers are free to move while the thumb is effectively restricted
- Provides excellent fit and high comfort by customizing anatomical shape
- Individually adjustable to the palm by molding the aluminum core
- Maximum wearing comfort thanks to the open, low-profile design



	REF/BNR		Size	Measurements	Color	UOM
	Left	Right				
	7623801	7623800	S	2.75 in - 3 in	Gray	1 Thumb Orthosis
	7623803	7623802	M	3 in - 3.5 in		
	7623805	7623804	L	3.5 in - 4 in		
	7623819	7623818	S	2.75 in - 3 in	Black	1 Thumb Orthosis
	7623821	7623820	M	3 in - 3.5 in		
	7623823	7623822	L	3.5 in - 4 in		

Actimove® Knee Brace

Adjustable Horseshoe, Simple Hinges, Condyle Pads

- Hinged knee sleeve for extra stability, adjustable horseshoe for improved patellar tracking
- Pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage knee pain from soft tissue and ligament injuries through extra stability with medial and lateral hinges
- Optional patellar support with an adjustable horseshoe
- Stays in place with improved fit through elastic hook and loop support strap
- Easy application with convenient finger-loops
- Improved patellar tracking through increased pressure by adjusting the horseshoe either inferiorly, superiorly, medially or laterally
- Can be worn left or right


	REF/BNR	Size	Measurements	Color	UOM
	7568750	XS	12 in - 14 in	Black	1 ea
	7568751	S	14 in - 16 in		
	7568752	M	16 in - 18 in		
	7568753	L	18 in - 20 in		
	7568754	XL	20 in - 22 in		
	7568755	XXL	22 in - 24 in		
	7568756	XXXL	24 in - 26 in		



Actimove® Knee Brace

Wrap Around, Polycentric Hinges, Condyle Pads

- Supports easy application for sore knees through front wrap closure, with polycentric hinges designed to match natural movement
- Pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage knee pain from soft tissue, ligament injuries or osteoarthritis through extra stability from medial and lateral hinges
- Stays in place with improved fit through elastic hook and loop support straps
- Can be worn left or right

	REF/BNR	Size	Measurements	Color	UOM
	7311800	XS	12 in - 14 in	Black	1 ea
	7311801	S	14 in - 16 in		
	7311802	M	16 in - 18 in		
	7311803	L	18 in - 20 in		
	7311804	XL	20 in - 22 in		
	7311805	XXL	22 in - 24 in		
	7311806	XXXL	24 in - 26 in		



Actimove® Knee Brace

Wrap Around, Simple Hinges, Condyle Pads

- Front wrap closure system supports application, simple hinges for extra stability
- Pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage knee pain from soft tissue, ligament injuries or osteoarthritis through extra stability from medial and lateral hinges
- Stays in place with improved fit through elastic hook and loop support strap
- Can be worn left or right



COOLMAX®
PRO

	REF/BNR	Size	Measurements	Color	UOM
	7550010	XS	12 in - 14 in	Black	1 ea
	7550011	S	14 in - 16 in		
	7550012	M	16 in - 18 in		
	7550013	L	18 in - 20 in		
	7550014	XL	20 in - 22 in		
	7550015	XXL	22 in - 24 in		
	7550016	XXXL	24 in - 26 in		

Actimove® Knee Stabilizer

Adjustable Horseshoe and Stay

- Sleeve design with adjustable horseshoe for improved patellar tracking
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Manage patellar pain syndrome with an adjustable horseshoe
- Added knee support from medial and lateral spiral stays
- Stays in place with improved fit through elastic hook and loop support strap
- Easy application with convenient finger-loops
- Improved patellar tracking through increased pressure by adjusting the horseshoe either inferiorly, superiorly, medially or laterally
- Can be worn left or right



COOLMAX®
PRO

	REF/BNR	Size	Measurements	Color	UOM
	7245300	XS	12 in - 14 in	Black	1 ea
	7245301	S	14 in - 16 in		
	7245302	M	16 in - 18 in		
	7245303	L	18 in - 20 in		
	7245304	XL	20 in - 22 in		
	7245305	XXL	22 in - 24 in		
	7245306	XXXL	24 in - 26 in		

Actimove® Knee Support

Closed Patella

- Fast, easy and secure protection for minor knee injuries
- Enjoy sports, as pain relief and healing from contusions, minor Meniscal damage and injuries, swelling from sprains and strains are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Feel secure as knee is protected by a snug fit that stays in place
- Fast & easy slip-on style



	REF/BNR	Size	Measurements	Color	UOM
	7558620	XS	12 in - 14 in	Navy	1 ea
	7558621	S	14 in - 16 in		
	7558622	M	16 in - 18 in		
	7558623	L	18 in - 20 in		
	7558624	XL	20 in - 22 in		
	7558625	XXL	22 in - 24 in		
	7558626	XXXL	24 in - 26 in	Black	1 ea
	7558627	XS	12 in - 14 in		
	7558628	S	14 in - 16 in		
	7558629	M	16 in - 18 in		
	7558630	L	18 in - 20 in		
	7558631	XL	20 in - 22 in		
	7558632	XXL	22 in - 24 in		
	7558633	XXXL	24 in - 26 in		



Actimove® Post Op Rom

Post Operative Knee Brace

- Telescopic hinge to fit most patients
- Range of motion brace. Set from -10 degrees to 120 degrees

	REF/BNR	Size	Color	UOM
	7627858	Universal	Black	1 ea



Actimove® GenuMotion Knee Support

Knee Support

- Helps pain relief, as the support offers stability through the anatomical shape and the 3D flat-knitted structure, providing gradient medical compression and a secure fit
- Contributes to pain relief and healing by reducing edema, through targeted pressure and an intermittent massage effect thanks to the anatomically shaped circular pressure pad
- For more wearing comfort especially designed soft, smooth Motion-Comfort-Zone that helps to prevent wrinkling in the popliteal area
- Anti-slip for pressure-reduced distal and proximal cuffs without an elastic band preventing constriction
- Added stabilization with great mobility thanks to the lateral spiral stays
- Skin-friendly and breathable 3D knit using OEKO-TEX® certified combination of yarn materials



	REF/BNR	Size	Measurements		Color	UOM
			Above Knee	Below Knee		
	7346810	XXS	13.75 in - 15 in	9.875 in - 11 in	White	1 ea
	7346811	XS	15 in - 16.125 in	11 in - 12.25 in		
	7346812	S	16.125 in - 17.375 in	12.25 in - 13.375 in		
	7346813	M	17.375 in - 18.5 in	13.375 in - 14.625 in		
	7346814	L	18.5 in - 19.75 in	14.625 in - 15.75 in		
	7346815	XL	19.75 in - 20.875 in	15.75 in - 17 in		
	7346816	XXL	20.875 in - 22 in	17 in - 18.125 in		
	7346850	3XL	22 in - 23.75 in	18.125 in - 19.25 in		
	7346851	4XL	23.25 in - 24.375 in	19.25 in - 20.5 in		
	7346817	M+	18.5 in - 19.75 in	13.375 in - 14.625 in		
	7346809	L+	19.75 in - 20.875 in	14.625 in - 15.75 in		
	7346808	XL+	20.875 in - 22 in	15.75 in - 17 in		
	7346807	XXL+	22 in - 23.25 in	17" - 18.125 in		
	7346852	3XL+	23.25 in - 24.375 in	18.125 in - 19.25 in		
	7346853	4XL+	24.375 in - 25.625 in	19.25 in - 20.5 in		
	7346843	XXS	13.75 in - 15 in	9.875 in - 11 in	Charcoal	1 ea
	7346844	XS	15 in - 16.125 in	11 in - 12.25 in		
	7346845	S	16.125 in - 17.375 in	12.25 in - 13.375 in		
	7346846	M	17.375 in - 18.5 in	13.375 in - 14.625 in		
	7346847	L	18.5 in - 19.75 in	14.625 in - 15.75 in		
	7346848	XL	19.75 in - 20.875 in	15.75 in - 17 in		
	7346849	XXL	20.875 in - 22 in	17 in - 18.125 in		



Actimove® Thigh Support

Adjustable

- Adjustable, accommodating fit that stays in place for strains and contusions
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Feel secure as thigh is protected by an accommodating fit that stays in place
- Adjust the level of support for comfort, no matter the activity
- Can be worn left or right



COOLMAX®
PRO

	REF/BNR	Size	Measurements	Color	UOM
	7560202	Universal	17.75 in - 31.5 in	Black	1 ea

Actimove® PowerMotion Thigh Support

- Aids pain relief and enhances healing by stimulating blood circulation through the unique Y-shaped strap design providing additional adjustable compression
- Helps to reduce edema and enhances pain reduction by simulating taping techniques normally used for strains and fiber ruptures of the calf muscle
- Promotes injured thigh muscle regeneration by controlling its elongation through targeted compression and Y-shaped straps
- Comfortable to wear over the highly sensitive injured area thanks to the specially designed 3D flat-knitted Pressure-Relief-Zone
- Secure fit due its anatomical shape and adjustable straps. Fast and easy to apply for the patient when compared to professional taping techniques, thanks to the individually adjustable closure strapping system



	REF/BNR	Size	Measurements	UOM
	7286706	S	16.125 in - 17.75 in	1 ea
	7286707	M	17.75 in - 19.75 in	
	7286708	L	19.75 in - 22.125 in	
	7286709	XL	22.125 in - 24.875 in	
	7286710	XXL	24.875 in - 27.5 in	



Actimove® PowerMotion Calf Support

- Aids pain relief and enhances healing by stimulating blood circulation through the unique Y-shaped strap design providing additional adjustable compression
- Helps to reduce edema and enhances pain reduction by simulating taping techniques normally used for strains and fiber ruptures of the calf muscle
- Promotes injured calf muscle regeneration by controlling its elongation through targeted compression and Y-shaped straps
- Comfortable to wear over the highly sensitive injured area thanks to the specially designed 3D flat-knitted Pressure-Relief-Zone
- Secure fit due its anatomical shape and adjustable straps
- Fast and easy to apply for the patient when compared to professional taping techniques, thanks to the individually adjustable closure strapping system



	REF/BNR	Size	Measurements	UOM
	7286604	S	11.5 in - 13.375 in	1 ea
	7286605	M	13.375 in - 15.75 in	
	7286606	L	15.75 in - 18.5 in	



Actimove® Talocast Ankle Brace

- Engineered for functional performance to help manage pain while providing a comfortable fit and ease of use
- Anatomically shaped stabilizers that support the ankle joint
- Comes with a long lasting and easy to adjust heel strap
- Offers optional cold therapy for extra pain relief

Actimove Professional Line TaloCast AirGel

	REF/BNR	Size	Measurements	UOM
	7311900	S/M Trainer	< 5'4 in / < 165 cm	1 Ankle Brace
	7311901	L/XL Standard	< 5'4 in / < 165 cm	1 Ankle Brace



- Engineered for functional performance to help manage pain while providing a comfortable fit and ease of use
- Anatomically shaped stabilizers that support the ankle joint
- Comes with a long lasting and easy to adjust heel strap

Actimove Professional Line TaloCast Air

	REF/BNR	Size	Measurements	UOM
	7311904	S/M Trainer	< 5'4 in / < 165 cm	1 Ankle Brace
	7311905	L/XL Standard	< 5'4 in / < 165 cm	1 Ankle Brace



Actimove® Standard Walker Low/High

- The Actimove® Standard Walker offer impressive features
- The high-quality rigid outer skeleton for protection of the foot, ankle and lower leg
- Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient
- Walking with Actimove Walker products feels natural through the rocker bottom which simulates a physiological gait



REF/BNR	Size	Shoe Sizes		Color	UOM
		Male	Female		
7627050	XS - Low	< 4.5	< 6	Black	1 Ankle Walker
7627051	S - Low	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627052	M - Low	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627053	L - Low	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627054	XL - Low	> 12.5	> 13.5	Black	1 Ankle Walker
7627055	XS - High	< 4.5	< 6	Black	1 Ankle Walker
7627056	S - High	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627057	M - High	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627058	L - High	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627059	XL - High	> 12.5	> 13.5	Black	1 Ankle Walker

Actimove® Standard Air Walker Low/High

- The Actimove® Standard Air Walker offer impressive features
- High-quality rigid outer skeleton for protection of the foot, ankle and lower leg
- Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient
- Walking with Actimove Walker products feels natural through the rocker bottom which simulates a physiological gait



REF/BNR	Size	Shoe Sizes		Color	UOM
		Male	Female		
7627250	XS - Low	< 4.5	< 6	Black	1 Ankle Walker
7627251	S - Low	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627252	M - Low	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627253	L - Low	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627254	XL - Low	> 12.5	> 13.5	Black	1 Ankle Walker
7627255	XS - High	< 4.5	< 6	Black	1 Ankle Walker
7627256	S - High	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627257	M - High	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627258	L - High	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627259	XL - High	> 12.5	> 13.5	Black	1 Ankle Walker

Actimove® Closed Shell Air Walker Low/High



- The Actimove® Closed Shell Walker offer impressive features
- High-quality rigid outer skeleton for protection of the foot, ankle and lower leg
- Circumferential strapping, a shock absorbing sole and breathable materials provide great comfort for the patient
- Walking with Actimove Walker products feels natural through the rocker bottom which simulates a physiological gait

REF/BNR	Size	Shoe Sizes		Color	UOM
		Male	Female		
7627630	XS - Low	< 4.5	< 6	Black	1 Ankle Walker
7627631	S - Low	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627632	M - Low	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627633	L - Low	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627634	XL - Low	> 12.5	> 12	Black	1 Ankle Walker
7627635	XS - High	< 4.5	< 6	Black	1 Ankle Walker
7627636	S - High	4.5 - 7	6 - 8	Black	1 Ankle Walker
7627637	M - High	7.5 - 10	8.5 - 11	Black	1 Ankle Walker
7627638	L - High	10.5 - 12.5	11.5 - 13.5	Black	1 Ankle Walker
7627639	XL - High	> 12.5	> 12	Black	1 Ankle Walker
7627796	One Size	Actimove Heel Wedges For Closed Shell Walker		Black	1 Heel Walker

Actimove® Ankle Stabilizer

Criss-Cross Straps

- Simulating sports taping, increased stability for ankle and arch of foot
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Great comfort with exceptional breathability through quick drying performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Increased stability around the ankle and arch of foot with criss-cross straps simulating sports taping
- Stays in place during athletic activity with extra support from circular strap
- Can be worn left or right



	REF/BNR	Size	Measurements	Color	UOM
	7561130	Universal	8 in - 10.5 in	Black	1 ea

Actimove® Patella Strap

Adjustable

- Targeted patella pain relief that stays in place
- Enjoy sports, as pain relief and healing are supported through balanced warmth and medical compression
- Comfort and breathability through quick drying perforated performance material
- Neoprene-free, made with COOLMAX® AIR technology
- Targeted pain relief by increased pressure on patella tendon from the high density EVA shaped insert
- Stays in place under activity due to silicone dot lining
- Adjust the level of support, for comfort no matter the activity
- Can be worn left or right



	REF/BNR	Size	Measurements	Color	UOM
	7558920	Universal	11.5 in - 16.75 in	Black	1 ea



Wound Care

Essity has a full line of wound care products to support your athletes on and off the field.

Portfolio overview

STRETCHABLE LIGHT-
WEIGHT FABRIC DRESSING



Leukoplast® elastic

TRANSPARENT DRESSING
WITH ABSORBENT PAD



Leukomed® T Plus

WIDE AREA FIXATION



Leukoplast® Hypafix®

ELASTIC GAUZE BANDAGE



Elastomull®

HYDROPHOBIC BACTERIA
BINDING POST OP DRESSING



Leukomed® Sorbact®

CONFORMING, ADHESIVE
FIXATION BANDAGE



Cover-Roll® Stretch

WOUND & EYE
IRRIGATION



SplashCap®



EyeCap®

Leukoplast® elastic

Fabric Adhesive Bandages

- Absorptive, non-adherent pad
- Adhesive surrounds the island pad to seal the affected area
- Contours to different areas on the body
- Variety of shapes and sizes



formerly known as Coverlet®

REF/BNR	Size	Description	UOM
		Strips	
7645617	¾ in x 3 in		100/box
7645606	1 in x 3 in		100/box
7645604	1 in x 3 in		300/box
		Spots	
7645616	7/8 in in round		100/box
7645613	1¼ in in oval		100/box
		Patches	
7645612	1½ in x 2 in		100/box
7645610	2 in x 3 in		50/box
7645605	4 in x 2¾ in		50/box
		Knuckles	
7645607	1½ in x 3 in		100/box
		Toe Shield (non-sterile)	
7645618	2½ in x 1⅜ in		100/box
		Digits	
7645608	Sm. finger tip		100/box
7645609	Lg. finger tip		50/box
		4-Wing	
7645611	3 in x 3 in		50/box
7646600	—	Eye Occluder REG (non-sterile)	30/box
7646700	—	Eye Occluder JR (non-sterile)	30/box
7253702	—	Cabinet (empty)	1 unit
7253703	—	Cabinet (filled)	1 unit

Leukomed® T Plus

Transparent, Total-Island Dressing with Absorbent Non-Adherent Pad

- Waterproof and barrier to bacteria
- Easy to apply with the unique red stripe system
- Should clinical conditions allow, the dressing may remain in place for up to 7 days




REF/BNR	Dressing Size	Wound Pad Size	UOM
7238200	2 in x 3 in	1 1/8 in x 1.5 in	50/box
7238201	3 1/8 in x 4 in	1.5 in x 2 3/8 in	50/box
7238202	3 1/8 in x 6 in	1.5 in x 4 3/8 in	50/box
7238214	4 in x 7 7/8 in	2 in x 6 3/8 in	50/box
7238203	4 in x 10 in	2 in x 8 1/8 in	50/box
7238204	4 in x 11 7/8 in	2 in x 9 3/4 in	50/box
7238205	4 in x 13 3/4 in	2 in x 11 3/4 in	50/box

4.3 WIDE AREA FIXATION

Leukoplast® Hypafix®

Wide Area Fixation

- Secure, rapid fixation that is excellent on contoured parts of the body, can be used as **Turf-Tape**
- Entire wound dressing is protected against contamination
- Convenient cut-to-size roll eliminates unnecessary waste



REF/BNR	Size	UOM
7193921	2 in x 2 yds (5 cm x 2 m)	36/box
7193919	4 in x 2 yds (10 cm x 2 m)	12/box
7193920	6 in x 2 yds (15 cm x 2 m)	12/box
7144301	2 in x 11 yds (5 cm x 10 m)	24/box
7144302	4 in x 11 yds (10 cm x 10 m)	24/box
7144303	6 in x 11 yds (15 cm x 10 m)	24/box

Elastomull®

Elastic Gauze Bandage

- High density fabric for a soft, rich texture
- Performance and stretchability to help control bleeding and edema
- Precision fit in areas with constant movement
- Helps prevent slippage once bandage is in place



REF/BNR	Size	UOM / Case
Non-sterile		
0209400	1 in x 4.4 yds	20 rolls
0209500	2 in x 4.4 yds	20 rolls
0209600	3 in x 4.4 yds	20 rolls

4.5 HYDROPHOBIC BACTERIA BINDING POST OP DRESSING

Leukomed® Sorbact®

Post-Operative Dressing with Sorbact® Technology

- Transparent and waterproof film with its wound contact layer coated with Sorbact® Technology
- Sorbact® Technology binds, inhibits, and removes hydrophobic bacteria utilizing a physical mode of action, without releasing active substances into the wound
- Reliable adhesion that is breathable and shower-proof
- Absorbent wound pad for low to moderate wound exudate
- Can be worn for up to 7 days depending on the clinical conditions of the wound
- For low to moderate exudate
- Sterile and not made with natural rubber latex



REF/BNR	Size	UOM
7619935	2 in x 3 in	50/box
7619936	3 in x 4 in	20/box
7619937	3 in x 6 in	20/box
7619938	4 in x 8 in	20/box
7619939	4 in x 10 in	20/box
7619940	4 in x 12 in	20/box
7619941	4 in x 14 in	20/box

Cover-Roll® Stretch

Conforming, Adhesive Fixation Bandage

- Stretchable non-woven polyester material adapts well to body contours, can be used as **Turf-Tape**
- Highly conformable to contour to all parts of the body requiring flexibility
- Secure fixation

REF/BNR	Size	UOM
4554700	2 in x 2 yds	36 rolls/case
4554800	4 in x 2 yds	12 rolls/case
4554900	6 in x 2 yds	12 rolls/case
4555000	8 in x 2 yds	12 rolls/case
4555100	12 in x 2 yds	12 rolls/case
4555200	2 in x 10 yds	12 rolls/case
4555300	4 in x 10 yds	12 rolls/case
4555400	6 in x 10 yds	12 rolls/case
4555500	8 in x 10 yds	12 rolls/case
4555600	12 in x 10 yds	12 rolls/case



4.7 WOUND & EYE IRRIGATION

SplashCap®

Irrigate and flush wounds rapidly

- Rapidly flushes contaminants from wounds
- Acute wounds
- Use with standard irrigation bottle not included

REF/BNR	Color	UOM
7446800	Clear	25/case
7446801	Clear	100/case



EyeCap®

Irrigate and flush eyes quickly

- Rapidly flushes contaminants from eyes
- Eye contamination
- Use with standard irrigation bottle not included

REF/BNR	Color	UOM
20025	Clear	25/case



