

Moisture management solutions

Powerful technology to fight skin breakdown



Medline makes it easy for you to protect against MASD

Moisture-associated skin damage (MASD) encompasses distinct skin conditions caused by excessive and continued exposure to moisture: wound exudate, urinary and/or fecal incontinence, or perspiration. Identifying the cause of skin damage helps ensure appropriate management and prevention interventions.¹


Negative outcomes related to poor moisture management

- + Incontinence-Associated Dermatitis (IAD)**
IAD occurs when urine and/or feces come into sustained contact with skin. This results in painful, weepy partial-thickness skin erosions.¹
- + Pressure injuries**
Incontinence, moisture, friction and shear make skin vulnerable to pressure injuries.²
- + Intertriginous Dermatitis (ITD)**
ITD occurs when perspiration becomes trapped in skin folds and can't evaporate, inflaming and denuding skin, making it more prone to infection.¹

Maintaining a favorable microclimate is a vital factor in helping control MASD and associated skin erosion

However, the likelihood that your staff has it top of mind all day every day isn't likely. Especially when a patient presents with multiple comorbidities.

That's where we come in.



How does your facility manage moisture?

We have your comprehensive solution

Pairing the right products with the right guidance gives your team the tools and assurance they need to manage the microclimate intuitively, with confidence. Medline's moisture management solutions empower your staff to consistently deliver best practices in skin care, positively impacting your clinical and financial outcomes.

Our comprehensive moisture management portfolio and selection guide helps your caregivers make the right choice for every patient, every time.

We can help your facility reduce cost, improve outcomes and increase patient comfort



Realize cost savings

Our products help reduce the risk and associated cost of skin damage, which allows staff to be more efficient with their care



Improve clinical outcomes

Properly absorbed and channeled moisture may help reduce moisture-associated skin damage in healthcare settings



Increase patient satisfaction

When moisture is controlled and risks are reduced, patient satisfaction levels increase

How it works

See what's inside our moisture-wicking products



Nonwoven

Provides a soft feel to product and adds strength/durability



Super-absorbent polymer

Locks in moisture and provides a dry feeling



Wood pulp fluff

Wicks moisture into product and throughout the core



Air-permeable poly

Allows air flow through the product, but prevents moisture from seeping through

Medline has industry-leading moisture management solutions

	Ultrisorbs Page 6	QuickChange Page 8	Underpads Page 10	Dry sheets Page 12
				
Clinical use	Open airing, incontinence and moisture protection	Light-to-moderate urinary incontinence for men	Procedures and bed protection	Skinfold moisture management
Area of hospital	ICU, MedSurg, O.R., Outpatient Surgery	ICU, MedSurg, ER, Neuro	Patient care, ER, Labor and Delivery	Patient care
Benefits	Draws in moisture, locks it away from the skin and feels dry to the touch	Contains urinary output before it contacts the skin to support skin health	Absorbs moisture and protects bedding, equipment and surfaces	Wicks and locks away moisture from skinfolds

Briefs

Page 14



Heavy urinary
and bowel
incontinence

Patient care

Dual core wicks
away moisture
and provides
extended
protection and
containment

Underwear

Page 17



Light-to-
moderate urinary
incontinence

Patient care

Exceptional
comfort and
discreet protection;
easy to pull
on and off

Maternity and pads

Page 19



Absorbency
and containment
designed for
postpartum flow

Labor and
Delivery,
Post-Op

Comfort, security
and protection
for new moms

Ultrisorbs®



To help prevent and manage MASD, it's important to use products that absorb and hold moisture away from the skin.³ Unlike conventional underpads, Ultrisorbs feature advanced polymers that help ensure ongoing dryness of skin⁴, help reduce friction and shear, allow air circulation and support low air-loss mattresses.

Recommended for:

- Facilities that practice open air or briefless incontinence management
- Patients on a low air loss, alternating therapy bed or other specialty support surface
- Non-ambulatory and incontinent patients where moisture management and skin breakdown is a concern
- Caregivers who want an all-in-one solution for moisture management and patient turns, boosts and repositioning



High-strength backsheet

Provides exceptional resistance to tearing



Ultra-soft, air permeable materials

For use with low air-loss mattresses and reducing heat buildup



Advanced super-absorbent core

Lies flat when soiled without bunching, swelling or disintegrating



AquaShield film

Prevents strikethrough

Ultrisorbs®

- Superabsorbent polymer core wicks away and absorbs fluid and lies flat without bunching, swelling or disintegrating
- Ultrasoft, breathable materials provide protection against friction and shear to offer patient comfort
- Air-permeable back sheet is for use on low air loss mattresses and specialty support surfaces
- Weight capacity options available for boosting/repositioning patients

	Ultrisorbs	Ultrisorbs Advanced	Ultrisorbs Advanced+	Ultrisorbs Complete
Weight Capacity	N/A	300–375 lbs.	300–375 lbs.	400 lbs.
Core Style	Standard	Standard	Extended*	Rapid wicking
Air Permeability	Yes	Yes	Yes	Yes

*20% larger core

Ultrisorbs

Item No.	Description	Weight Cap.	Pkg.
ULTRASORB1016	Ultrisorbs 10" x 16"	n/a	10/bg, 100/cs
EXTRASRB2336A	Ultrisorbs 23" x 36"	n/a	10/bg, 70/cs
EXTRASRB3036A	Ultrisorbs 30" x 36"	n/a	5/bg, 70/cs



EXTRASRB3036A

Ultrisorbs Advanced

Item No.	Description	Weight Cap.	Pkg.
MSC8003003	Ultrisorbs Advanced 30" x 36"	300 lbs.	5/bg, 70/cs
EXTSB3036A350	Ultrisorbs Advanced 30" x 36"	350 lbs.	5/bg, 70/cs
USAP3136ESLP	Ultrisorbs Advanced 30" x 36"	375 lbs.	10/bg, 50/cs
MSC383500*	Ultrisorbs Advanced with SlipResist 30" x 36"	300 lbs.	5/bg, 70/cs

*Patent-pending



MSC8003003

Ultrisorbs Advanced+

Item No.	Description	Weight Cap.	Pkg.
ULTRASORB2436	Ultrisorbs Advanced+ 23" x 36"	300 lbs.	10/bg, 70/cs
ULTRASORB3136	Ultrisorbs Advanced+ 30" x 36"	300 lbs.	10/bg, 40/cs
USAP3136ES	Ultrisorbs Advanced+ 30" x 36"	375 lbs.	10/bg, 40/cs
MSC282900	Ultrisorbs Advanced+ 40" x 57"	375 lbs.	5/bg, 60/cs



ULTRASORB3136

Ultrisorbs Complete

Item No.	Description	Weight Cap.	Pkg.
MDP733036	Ultrisorbs Complete 30" x 36"	400 lbs.	5/bg, 60/cs



MDP733036

QuickChange[®] male incontinence wrap

The one-size-fits-all QuickChange wrap applies directly around the male anatomy to absorb urinary output before it makes contact with the skin, reducing the risk of incontinence-associated dermatitis or pressure injury

Recommended for:

- “Open to air” facility/brief-free
- Sleeping and non/semi ambulatory patients



Innovative, one-size-fits-all design

Ensures application regardless of differing male anatomy



Superabsorbent core

Draws moisture in and away from skin to help prevent liquid from pooling onto skin






Skin safe hook-and-loop closure

Removes the need to apply adhesives to the skin

How to use

1. Place QuickChange onto body, absorbent/white side facing up
2. Fold pointed flap around the penis
3. Wrap hook-and-loop side around to secure QuickChange closed



QuickChange® wrap

- Absorbent core wicks urinary output in to keep skin dry
- Apply directly to the male anatomy
- One-size-fits-all for non/semi ambulatory male patients

Item No.	Description	Pkg.
UIM1025	Male incontinence wrap	25/bg, 100/cs
UIM1010	Male incontinence wrap	10/bg, 100/cs



UIM1025

How QuickChange helps improve outcomes



Skin health

Absorbs urine before it pools onto skin or soils linens when open airing



Staff and patient safety

No need to log roll to change product. One staff member job



Patient dignity

Keeps patient dry, improving patient dignity and comfort

Underpads

Medline underpads provide outstanding absorbency to keep moisture away, while protecting bedding, furniture, equipment and surfaces

Recommended for:

- Bed and furniture protection
- Clinical procedures
- Wheelchair seats
- IV starts
- Dressing changes



Waterproof backing

Helps protect bedding and holds pad in place



Advanced super-absorbent core

Polymer and fluff core locks in moisture



Moisture wicking topsheet

Helps keep wetness away from skin

Disposable underpads

- Polypropylene backing helps keep underpad in place
- Provide moisture protection for a variety of clinical surfaces
- Disposable underpad features a combination fluff and polymer core for high absorbency and odor control

17" x 24" (24 x 61 cm)

Item No.	Color	Absorbency	Quilted	Pkg.
MSC281224C	Blue	Light	✓	25/bg, 300/cs
MSC281230	Blue	Moderate	—	300/cs

For IV starts
and blood draws

23" x 24" (58 x 61 cm)

Item No.	Color	Absorbency	Quilted	Pkg.
MSC281225	Blue	Light	—	200/cs
MSC281241LB	Blue	Moderate	—	5/bg, 200/cs

For seat covers,
wheelchair covers, IV starts,
dressing changes

23" x 36" (58 x 91 cm)

Item No.	Color	Absorbency	Quilted	Pkg.
MSC281229C	Blue	Light	✓	25/bg, 150/cs
MSC281236P	Blue	Moderate	—	10/bg, 150/cs
MSC281264	Blue	Moderate	✓	25/bg, 150/cs
MUP2060P	Green	Heavy	—	10/bg, 80/cs

For bedding
protection

30" x 30" (76 x 76 cm)

Item No.	Color	Absorbency	Quilted	Pkg.
MSCB281245LB	Blue	Light	✓	5/bg, 150/cs
MUP1266P	Blue	Moderate	—	10/bg, 90/cs
MSC282045	Teal	Heavy	✓	10/bg, 100/cs
MUP0305P	Green	Heavy	—	10/bg, 90/cs

For clinical
procedures

30" x 36" (76 x 91 cm)

Item No.	Color	Absorbency	Quilted	Pkg.
MSC282033P	Green	Moderate	—	10/bg, 100/cs
MUP0365P	Green	Heavy	—	10/bg, 80/cs
MSC282030LBC	Green	Heavy	✓	5/bg, 75/cs

For bedding
protection,
labor and
delivery
and O.R.

36" x 36" (91 x 91 cm)

Item No.	Color	Absorbency	Quilted	Pkg.
MSC281271P	Green	Moderate	✓	5/bg, 50/cs
MSC282070LB	Green	Heavy	✓	5/bg, 50/cs

For bedding
protection



MSC281236P



MUP0365P

Specialty underpads

Item No.	Description	Color	Absorbency	Quilted	Pkg.
MSC280914	9" x 14" (23 x 35 cm), crib liner	White	Light	✓	50/bg, 600/cs
MSC282110	21" x 36" (53 x 91 cm), breathable	Blue	Heavy	✓	5/bg, 50/cs

Skinfold dry sheets



Medline Skinfold Dry Sheets absorb moisture, lock it away from skinfolds and keep patients dry reducing the risk of Intertriginous Dermatitis

Recommended for:

- Skinfold moisture
- High risk of ITD



Superabsorbent polymer core

Draws moisture in
away from skin



Pre-cut


Each sheet is pre-cut
to the perfect size



Soft and non- chafing materials


For cooling
and comfort

How to use




Assess and cleanse

Apply product to clean and dry skinfold*



Apply dry sheet

Place as deeply in skinfold as possible; fold sheet if necessary



Change as needed

Allow skinfold to close evenly over sheet.
Change product daily or when sheet becomes wet

*Always follow facility protocol when applying this product.

Dry sheets

- Designed to absorb moisture and help keep skinfolds dry
- Change once a day or when sheet becomes wet
- Each sheet measures 6" x 14" and should not be cut
- One sheet holds 8-oz. of liquid

Item No.	Description	Pkg.
MSC10614	6" x 14" (15 x 35 cm)	10/bg, 250/cs



MSC10614

How dry sheets help improve outcomes



Skin health

Absorbs moisture so it does not become trapped within skinfolds



Patient dignity

Keeps patient dry, improving patient dignity and comfort



Infection prevention

Reduces skin-on-skin friction that often leads to inflamed skin increasing the risk of infection

Briefs



FitRight® briefs address two important concerns—keeping skin dry and maintaining skin integrity with a unique dual core that wicks away moisture and provides extended protection

Recommended for:

- Bedridden patients
- Severe urinary and bowel incontinence
- Containment and maximum protection against moisture



State-of-the-art design

Allows improved fit around the legs and superior dryness



Breathable side panels

Allows airflow to the skin to increase user comfort



Soft backsheet

Cloth-like for extra softness and discreet feel

FitRight® briefs

- 4D-core with odor protection helps wick away fluid to help promote dryness
- Skin-safe closures provide secure, repeated refastening
- Ultrasoft cloth-like back sheet provides a discreet, garment-like feel and breathable side panels allow airflow for comfort
- Moisture alert indicator changes color when the garment has been soiled

FitRight briefs

Item No.	Description	Color	Pkg.
FITBASICMD	Medium, 32"-44" (81-112 cm)	White	25/bg, 100/cs
FITBASICLG	Large, 44"-56" (112-142 cm)	Blue	25/bg, 100/cs
FITBASICXLG	X-Large, 56"-64" (142-163 cm)	Beige	25/bg, 100/cs
FITBASICXXL	XX-Large, 60"-70" (152-178 cm)	Green	25/bg, 100/cs



FITBASICLG

FitRight Extra briefs

Item No.	Description	Color	Pkg.
FITEXTRASM	Small, 20"-32" (51-81 cm)	Peach	20/bg, 80/cs
FITEXTRAMD	Medium, 32"-44" (81-112 cm)	White	20/bg, 80/cs
FITEXTRALG	Large, 44"-56" (112-142 cm)	Blue	20/bg, 80/cs
FITEXTRAXLG	X-Large, 56"-64" (142-163 cm)	Beige	20/bg, 80/cs
FITEXTRAXXL	XX-Large, 60"-70" (152-178 cm)	Green	20/bg, 80/cs



FITEXTRALG

FitRight Plus briefs

Item No.	Description	Color	Pkg.
FITPLUSMD	Medium, 32"-44" (81-112 cm)	White	20/bg, 80/cs
FITPLUSLG	Large, 44"-56" (112-142 cm)	Blue	20/bg, 80/cs
FITPLUSXLG	X-Large, 56"-64" (142-163 cm)	Beige	20/bg, 80/cs
FITPLUSXXL	XX-Large, 60"-70" (152-178 cm)	Green	20/bg, 80/cs



FITPLUSLG

Anti-leak
guards

FitRight® stretch briefs

- 4D-core with odor protection helps wick away fluid to promote dryness and help maintain skin integrity
- Large back panel for increased backside coverage
- Ultrasoft cloth-like back sheet provides a discreet, garment-like feel and breathable side panels allow airflow for comfort
- Moisture alert wetness indicator changes color when garment has been soiled

FitRight Stretch Extra briefs

Item No.	Description	Color	Pkg.
FRSE1	Medium/Regular, 30"-52" (76-132 cm)	Lavender	20/bg, 80/cs
FRSE2	Large/X-Large/XX-Large, 51"-70" (130-178 cm)	Green	20/bg, 80/cs



FitRight stretch ultra briefs

Item No.	Description	Color	Pkg.
FRSU1	Medium/Regular, 30"-52" (76-132 cm)	Lavender	20/bg, 80/cs
FRSU2	Large/X-Large/XX-Large, 51"-70" (130-178 cm)	Green	20/bg, 80/cs



FRSU2



Reduced size system
Easier to size and helps drive compliance



Two zones of stretch
Improved stretch for better strength and fit



Larger tabs
One large, centered tab for increased ease of application

FitRight bariatric briefs

- High-capacity core with odor protection helps wick away fluid to promote dryness and help maintain skin integrity for long-lasting protection
- Flexible stretch panels offer a secure and comfortable fit and offer ease of use for staff
- Cloth-like outer cover provides a soft, natural feeling on the skin and discreet fit
- Skin-safe closures offer secure, repeated refastening and the moisture alert wetness indicator changes color when the garment has been soiled

Item No.	Description	Color	Pkg.
MTB084	3XL, up to 84" (213 cm),	White	12/bg, 48/cs
BARIBRIEFC	4XL, 65"-94" (165-239 cm)	White	8/bg, 32/cs
BARIBRIEF	4XL, 65"-94" (165-239 cm)	Beige	8/bg, 32/cs
BARIBRIEFPLUS	5XL, 90"-106" (229-269 cm)	Beige	8/bg, 32/cs



BARIBRIEF

Underwear

FitRight underwear provides exceptional comfort and worry-free containment. Discreet and disposable, they can be pulled on and off just like normal underwear.

Recommended for:

- Light to heavy urinary incontinence and/or occasional bowel incontinence episodes
- Patients able to walk with or without assistance
- Discreet, garment-like protection
- Those who are lightly incontinent during the day and may need a little extra protection at night



Body-contoured design

Designed to match the feel and comfort of everyday underwear



SensiSoft® fabric

Feels soft and soothing on skin



Superabsorbent core

Draws moisture into core to keep skin dry and maintain skin integrity

FitRight® underwear

- Features a body-contoured anatomical design, soft nonwoven fabric and multiple absorbency options
- Soft anti-leak guards help reduce leakage and improve containment (available on FitRight Ultra)
- The acquisition and dryness layer wicks moisture away from the skin and into the core for ideal dryness and leakage protection
- Linear embossing technology distributes fluid throughout the core for exceptional dryness and containment

FitRight ultra underwear

Item No.	Description	Pkg.
FIT23005A	Medium, 28"-40" (71-102 cm)	20/bg, 80/cs
FIT23505A	Large, 40"-56" (102-142 cm)	20/bg, 80/cs
FIT23600A	X-Large, 56"-68" (142-173 cm)	20/bg, 80/cs



FIT23505A

FitRight extra underwear

Item No.	Description	Pkg.
MSC13005A	Medium, 28"-40" (71-102 cm)	20/bg, 80/cs
MSC13505A	Large, 40"-56" (102-142 cm)	20/bg, 80/cs
MSC13600A	X-Large, 56"-68" (142-173 cm)	20/bg, 80/cs



MSC13505A

FitRight bariatric underwear

- Features a body-contoured anatomical design, SensiSoft fabric, soft anti-leak guards and a super-absorbency core
- The acquisition and dryness layer wicks moisture away from the skin and into the core for ideal dryness and leakage protection
- Linear embossing technology distributes fluid throughout the core for exceptional dryness and containment

Item No.	Description	Pkg.
FIT700A	XXL, 68"-80" (173-203 cm)	20/bg, 80/cs
FIT800	XXXL, 75"-94" (190-238 cm)	20/bg, 80/cs



FIT700A

Maternity and pads



Medline's maternity and incontinence pads are designed to keep users comfortable and dry. These pads can be paired with the user's own underwear or with Medline's knit or mesh pants.

Recommended for:

- Maternity
- Bladder control pads for light bladder leakage
- FitRight liners for moderate to heavy incontinence



Soft outer layer

Wicks moisture away from the skin and ensures user comfort



Superabsorbent core

Draws moisture in and away from skin to help prevent liquid from pooling onto skin



Multiple sizes

Available for versatility

FitRight® bladder control pads

- Discreet contoured design features absorbent core to lock in moisture and odor and keep skin dry
- Waterproof backsheet is soft, comfortable and fabric-like; won't make crinkle noises like plastic backsheets
- Moisture-wicking topsheet features a tasteful design
- Adhesive strip gently holds the pad in place in any type of underwear; individually wrapped
- Compares to national brands*

Item No.	Description	Pkg.
FBC1010	Light, 3½" x 9" (9 x 23 cm)	20/bg, 240/cs
FBC1020	Moderate, 4" x 11" (10 x 28 cm)	16/bg, 192/cs
FBC1030	Maximum, 5" x 14" (13 x 36 cm)	10/bg, 120/cs
FBC1040	Ultimate, 6" x 16" (15 x 40 cm)	10/bg, 120/cs
FBC1100	Overnight, 8.5" x 16" (21.6 x 40 cm)	24/bg, 144/cs

*Data on file.



FBC1030

Medline ContourPlus bladder control pads

- Cloth-like outer cover provides a soft feeling on the skin and discreet fit
- Elastic leg gathers help prevent leaks
- Super absorbent core keeps moisture away from skin and controls odors
- Contoured design offers a comfortable fit

Item No.	Description	Pkg.
BCPE01	Moderate, 5½" x 10½" (14 x 26 cm)	28/bg, 336/cs
BCPE02	Maximum, 6½" x 13½" (16 x 34 cm)	28/bg, 168/cs
BCPE03	Ultimate, 8" x 17" (20 x 43 cm)	28/bg, 168/cs



BCPE02

FitRight liners


- 4D-core with odor protection for super absorbency
- Cloth-like backsheet for increased comfort and discretion
- Low-profile core for improved comfort and fit

FitRight Plus liners

Item No.	Description	Color	Pkg.
FITLINER300	Plus, 13" x 30" (33 x 76 cm)	Green	20/bg, 80/cs
FITLINER300Z	Plus, 13" x 30" (33 x 76 cm)	Green	20/bg

FitRight Extra liners

Item No.	Description	Color	Pkg.
FITLINER100	Extra, 13" x 30" (33 x 76 cm)	Yellow	20/bg, 80/cs
FITLINER100Z	Extra, 13" x 30" (33 x 76 cm)	Yellow	20/bg



FITLINER300

Medium maternity pads

- Maximum absorbency for protection from leaks
- Variety of sizes, features and packaging options

With adhesive

Item No.	Description	Pkg.
NON241278	Maxi, non-sterile, 9" (20 cm), indiv. boxed	1/bx, 250/cs
NON241279	Maxi, non-sterile, 9" (20 cm), indiv. wrapped	12/bg, 600/cs
NON241286	Maxi, non-sterile, 11" (28 cm), indiv. wrapped	12/bg, 288/cs



NON241286

With adhesive and wings

Item No.	Description	Pkg.
NON241289	Maxi, non-sterile, 11" (28 cm), with wings, indiv. wrapped	12/bg, 288/cs



NON241281

Adhesive-free

Item No.	Description	Pkg.
NON241280	Non-sterile, 11" (28 cm), with tails	12/bg, 288/cs
NON241281	Sterile, 11" (28 cm), with tails, indiv. wrapped	12/bg, 288/cs
NON241283	Sterilizable, 11" (28 cm), with tails, indiv. wrapped	12/bg, 288/cs

Large maternity pads

- Super absorbent core locks fluids away from skin
- Adhesive strip attaches securely to underwear

Item No.	Description	Pkg.
MSC323950C	Contoured peach pad, 7" x 14"	20/bg, 240/cs
MSC323000C	Embossed liner, 7" x 17"	12/bg, 240/cs
MSC322505	Embossed liner, 9" x 24"	20/bg, 80/cs



MSC323950C

Medline Black knit pants

- Full-coverage black pant
- Perfect for postpartum wear, with pads and ice
- Can be used for C-section deliveries due to the high waist and won't irritate incisions
- Comfortable polyester and spandex blend
- Machine washable

Item No.	Description	Color	Pkg.
KNITPANTBL	Medium/Large, 20"-60" (51-152 cm) waist	Black	2/bg, 50/cs
KNITPANTBXL	X-Large, 45"-70" (114-178 cm) waist	Black	2/bg, 50/cs
KNITPANTBXXL	XX-Large, 50"--75" (127-191 cm) waist	Black	2/bg, 50/cs



KNITPANTBL

Medline Knit pants

- Versatile design works with any of Medline's incontinence liners and pads
- Can be washed multiple times
- Seamless design helps eliminate seams pressing into skin
- Knit weave is soft against the skin
- Color-coded waistbands make size identification simple

Premium

Item No.	Description	Color	Pkg.
MSC86300	Medium/Large, 20"-60" (51-152 cm) waist	Blue/Brown	5/bg, 100/cs
MSC86400	X-Large, 45"-70" (114-178 cm) waist	Green	5/bg, 100/cs
MSC86500	XX-Large, 50"-75" (127-191 cm) waist	Purple	5/bg, 100/cs
MSC86700	XXX-Large, 65"-85" (165-216 cm) waist	Red	2/bg, 50/cs



MSC86300



MSC76400

Deluxe

Item No.	Description	Color	Pkg.
MSC76400	Large/X-Large, 22"-55" (60-140 cm) waist	Green	100/cs
MSC76500	XX-Large, 24"-64" (61-163 cm) waist	Purple	100/cs

Protection Plus mesh pants

- Designed to hold a liner in place when pants are sized correctly
- Soft, comfortable fit
- Disposable, intended for one-time use

Item No.	Description	Color	Pkg.
MBP3202	Large, 30"-45" (76-114 cm)	White	400/cs
MBP3703	X-Large, 32"-48" (81-122 cm)	Green	300/cs
MBP3704	XX-Large, 38"-58" (97-147 cm)	Purple	300/cs



MBP3703

Wipes



Medline offers a full line of dry wipes to clean and care for delicate skin

Recommended for:

- Moisture cleanup
- Applying ointments, creams and protective barriers
- Pericare and bathing
- Use on fragile skin, where rough terry cloth may cause irritation
- General cleaning purposes



A softer, gentler touch

Ultra-soft, disposable cleansing wipes are gentle on skin



Convenient and easy to use

Ready-to-use pre-moistened cloths increase caregiver efficiency. Just open the pack, cleanse, and toss



Versatile

Different options for different needs



FitRight® UltraSoft™ dry wipes

- Reduce waste with single-wipe dispensing
- Save time and hassle with new pop-top dispenser
- Resealable container helps reduce the risk of cross-contamination

Item No.	Description	Pkg.
FRDW263854	Pop top, soft, cleansing dry wipes, 8" x 8"	1200/cs



FRDW263854

Absorbent Wipers+

- Disposable dry wipes with more strength and absorbency
- Can be used for general EVS, surface clean up, and incontinence care

Item No.	Description	Pkg.
AWP0500	Medline Absorbent Wipers+, 12" x 13"	90/bg, 1080/cs



AWP0500

NatureSoft dry wipes

- Dispersible material allows it to be flushed without clogging pipes
- Strong wipe is gentle to the touch yet tough enough for virtually all patient cleaning needs

Item No.	Description	Pkg.
NATURESOFT913	Flushable dry cleansing wipes, 9" x 13"	500/cs



NATURESOFT913

UltraSoft™ dry wipes

- Made of Soft-Spun fabric that is softer than reusable cloths
- Perfect for incontinence care, patient bathing, baby care, surface cleaning and more
- Disposable, non-flushable cloths won't dry out like traditional wet wipes
- No additives: Alcohol free, fragrance free, paraben free

Item No.	Description	Pkg.
ULTRASOFT1013	Soft, cleansing dry wipes, 10" x 13"	500/cs
ULTRASOFT1214	Soft, cleansing dry wipes, 12" x 14"	1080/cs
ULTRASOFT713	Soft, cleansing dry wipes, 7" x 13"	1200/cs
NON1013	Soft, cleansing dry wipes, 10" x 13"	1000/cs



ULTRASOFT1013

Clinical and product support

We understand that introducing new protocols or products can be a significant investment in time and effort. We can help. Our team of experts offers you clinical and product support to help ensure successful trials and implementations. We provide:



**Guided process
for trials and
implementations**



**Utilization
review and SKU
consolidation**



**30+ product
specialists to support
in onsite education**



**Business reviews
and follow ups**



**Communication
support and tools**



**Product usage and
best practice education**

Medline University

**Be a moisture management scholar.
Earn free CE credits from Medline University.**

No tuition. No classrooms. No kidding. When you enroll in Medline University you have unlimited access to hundreds of free courses about a variety of topics across the continuum of care, including moisture management.

Simply log on and learn from courses developed just for us by clinical leaders. Plus, you can earn and save your continuing education credits in one convenient place. Genius.



medlineuniversity.com

**Sign up
today**



**We make
healthcare
run better™**

Medline Industries, Inc.

Three Lakes Drive, Northfield, IL 60093
Medline United States | 1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada

1-800-268-2848 | medline.ca | canada@medline.com

Medline México

01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us



References: **1.** Dowsett D, Allen L. Moisture-Associated Skin Damage Made Easy. Wounds UK. 2013;9(4): 1-4. Available at http://www.wounds-uk.com/pdf/content_10961.pdf Accessed November 9, 2015. **2.** Cooper K. Evidence-Based Prevention of Pressure Ulcers in the Intensive Care Unit. Critical Care Nurse. 2013;33(6):57-66. Available at <http://ccn.aacnjournals.org/content/33/6/57> Accessed November 9, 2015. **3.** Dowsett D, Allen L. Moisture-Associated Skin Damage Made Easy. Wounds UK. 2013;9(4). Available at http://www.wounds-uk.com/pdf/content_10961.pdf Accessed November 9, 2015. **4.** Data on file. Available upon request. **5.** Data on file. Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

© 2021 Medline Industries, Inc. Medline, FitRight, Comfort Glide, 4D-Core, FitRestore, Medline Remedy, Protection Plus, Baribrief are registered trademarks and Comfort-Aire, Naturesoft, and Advancing the Health of Healthcare are trademarks of Medline Industries, Inc. MKT19W1679010 / ne19634 / 55