

Ostomy Product Guide

Julian
SenSura® Mito &
Brava® user



Welcome to O

Supporting

2-Pi

2-P


Fistula/Wound Man

SenSura® Mio Click

All the benefits of BodyFit
Technology in a 2-piece
mechanical coupling system



- Reassurance and security with a tactile lock ring
- Mechanical coupling can be repositioned while it remains on the patient's body




2-Piece SenSura® Mio Click

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click Barrier

- Extended wear
- Flat
- Belt tabs for use with 4237, 4247
- Not made with natural rubber latex



Code	Stoma Size	Color	Units: 5	HPCS*
Cut-to-fit				
10502	3/8-1 3/8" (10-35 mm)	Green		A4409
10512	3/8-1 3/4" (10-45 mm)	Red		A4409
10522	3/8-2 1/8" (10-55 mm)	Blue		A4410
10531	3/8-2 9/16" (10-65 mm)	Yellow		A4410
Pre-cut				
10503	3/4" (20 mm)	Green		A4409
10506	7/8" (22 mm)	Green		A4409
10504	1" (25 mm)	Green		A4409
10505	1 1/8" (28 mm)	Green		A4409
10513	1 3/16" (30 mm)	Red		A4409
10514	1 1/4" (32 mm)	Red		A4409
10515	1 3/8" (35 mm)	Red		A4409
10516	1 1/2" (38 mm)	Red		A4409
10523	1 9/16" (40 mm)	Blue		A4409
10524	1 3/4" (45 mm)	Blue		A4409
10525	2" (50 mm)	Blue		A4409

SenSura® Mio Click Convex Light Barrier

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex



Code	Stoma Size	Color	Units: 5	HPCS*
Cut-to-fit				
16901	3/8 - 3/4" (10-20 mm)	Green		A4407
16911	5/8 - 1 3/16" (15-30 mm)	Red		A4407
16921	5/8 - 1 9/16" (15-40 mm)	Blue		A4407
16931	5/8 - 2 1/16" (15-53 mm)	Yellow		A4408
Pre-cut				
16903	1 1/16" (18 mm)	Green		A4407
16913	1 3/16" (21 mm)	Red		A4407
16914	1" (25 mm)	Red		A4407
16915	1 1/8" (28 mm)	Red		A4407
16923	1 1/4" (31 mm)	Blue		A4407
16924	1 3/8" (35 mm)	Blue		A4407
16933	1 1/2" (38 mm)	Yellow		A4408
16934	1 5/8" (41 mm)	Yellow		A4408

**SenSura® Mio Click Convex Deep Barrier**

- Extended wear
- Belt tabs for use with 4237, 4247
- FlexShell is designed for stability and flexibility
- Not made with natural rubber latex



Code	Stoma Size	Color	Units: 5	HCPCS*
Cut-to-fit				
16941	3/8 - 3/4" (10-20 mm)	Green	A4407	
16951	5/8 - 1 3/16" (15-30 mm)	Red	A4407	
16961	5/8 - 1 3/16" (15-40 mm)	Blue	A4407	
16971	5/8 - 2 1/16" (15-53 mm)	Yellow	A4408	
Pre-cut				
16942	1 1/16" (18 mm)	Green	A4407	
16953	1 3/16" (21 mm)	Red	A4407	
16954	1" (25 mm)	Red	A4407	
16955	1 1/8" (28 mm)	Red	A4407	
16963	1 1/4" (31 mm)	Blue	A4407	
16964	1 3/8" (35 mm)	Blue	A4407	
16973	1 1/2" (38 mm)	Yellow	A4408	
16974	1 5/8" (41 mm)	Yellow	A4408	

SenSura® Mio Click Convex Flip Barrier

- Extended wear
- Designed for curves, bulges and hernias
- Curved, star-shaped barrier conforms to outward areas with less creasing and folding
- "Flip to fit" design allows for easy application
- Belt tabs for use with 4237, 4247
- Not made with natural rubber latex



Code	Stoma Size	Color	Units: 5	HCPCS*
Cut-to-fit				
18500	3/8 - 1 3/8" (10-35 mm)	Green	A4407	
18510	3/8 - 1 3/4" (10-45 mm)	Red	A4408	
18520	3/8 - 2 1/8" (10-55 mm)	Blue	A4408	
18530	3/8 - 2 9/16" (10-66 mm)	Yellow	A4408	
Max Cut				
Code	Stoma Size	Color	Units: 5	HCPCS*
Pre-cut				
18501	3/4" (20 mm)	Green	A4407	
18511	1" (25 mm)	Red	A4407	
18512	1 3/16" (30 mm)	Red	A4407	
18513	1 3/8" (35 mm)	Red	A4407	
18521	1 3/16" (30 mm)	Blue	A4407	
18522	1 3/8" (35 mm)	Blue	A4407	
18523	1 9/16" (40 mm)	Blue	A4407	
18524	1 3/4" (45 mm)	Blue	A4407	
18531	1 3/4" (45 mm)	Yellow	A4408	

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems.

SenSura® Mio Click MAXI Drainable Pouch

- EasiClose™ WIDE outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- Not made with natural rubber latex



11" (28 cm) / 655 mL Units: 10 HCPCS*

Code	Code	Color	HCPCS*
Transparent			
11463	11462	Green	A4427
11473	11472	Red	A4427
11483	11482	Blue	A4427
11487	11486	Yellow	A4427

SenSura® Mio Click MAXI Drainable Pouch

- EasiClose™ WIDE outlet
- **Without filter**
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



11" (28 cm) / 655 mL Units: 10 HCPCS*

Code	Color	HCPCS*
Transparent		
11474	Red	A4426
11484	Blue	A4426
11488	Yellow	A4426


SenSura® Mio Click MIDI Drainable Pouch

- EasiClose™ WIDE outlet
- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- Not made with natural rubber latex



10 1/2" (26 cm) / 485 mL Units: 10 HCPCS*


Code	Color	HCPCS*
Opaque gray		
11461	Green	A4427
11471	Red	A4427

**SenSura® Mio Click High Output Pouch with Tap Outlet**

- **Full-circle filter designed for ileostomy output**
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque includes inspection window
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- Not made with natural rubber latex

12½" (32 cm) / 1080 mL Units: 10

Code	Color	HCPCS*
Transparent		
18610		A4413
18611		A4413
18612		A4413
18613		A4413

**SenSura® Mio Click High Output Pouch with Soft Outlet**

- **Full-circle filter designed for ileostomy output**
- Neutral gray, water-resistant textile
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque includes inspection window
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- Not made with natural rubber latex

12½" (32 cm) / 1080 mL Units: 10

Code	Color	HCPCS*
Transparent		
18640		A4413
18641		A4413
18642		A4413

**SenSura® Mio Click Drainable Pouch with Tap Outlet**

- **Full-circle filter designed for ileostomy output**
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- Wider-lumen outlet provides effective drainage for patients with thicker output; connects to Coloplast fistula drain bag #14010
- Not made with natural rubber latex

9" (23 cm) / 650 mL Units: 10

Code	Color	HCPCS*
Opaque gray		
18620		A4427
18621		A4427
18622		A4427

**SenSura® Mio Click MAXI Closed Pouch**

- **Full-circle filter**
- Neutral gray, water-resistant textile
- Opaque includes inspection window
- Not made with natural rubber latex

8¼" (20½ cm) / 570 mL Units: 30

Code	Color	HCPCS*
Transparent		
N/A		A4423
N/A		A4423
N/A		A4423
11426		A4423

**SenSura® Mio Click MIDI Closed Pouch**

- **Full-circle filter**
- Neutral gray, water-resistant textile
- Peek-a-boo inspection window for easy application
- Not made with natural rubber latex

7¼" (18½ cm) / 275 mL Units: 30

Code	Color	HCPCS*
Opaque gray		
11401		A4423
11412		A4423
11421		A4423


Urostomy

SenSura® Mio Uro

Wide anti-reflux valves support varying ways of managing stents:

- Stents threaded into a pouch connected to a night drainage bag
- Rolling the stents in the top of the pouch, making a loop
- Stents that have been cut

Wide anti-reflux valves



Welded "multi-chambers" minimize ballooning

SenSura® Mio Click MAXI Urostomy Pouch

- Multi-chambers minimize sloshing
- Wider anti-reflux valves prevent urine from washing over the stoma
- Soft outlet can be connected to urostomy night bag #21365 or any standard night or leg bag
- Neutral gray, water-resistant textile
- Not made with natural rubber latex



10½" (26½ cm) / 480 mL

Units: 10

Code	Code	Color	HCPCS*
Transparent	Opaque gray		
11492	11491	Green	A4433
11498	11497	Red	A4433
11500	11499	Blue	A4433

SenSura® Mio Click Barrier + Drainable Pouch Single Use Kit

- Extended wear
- Barrier #10522 Flat with belt tabs
- Pouch #11483 MAXI / EasiClose™ WIDE outlet
- Instructions for use
- Stoma measuring guide
- Not made with natural rubber latex



Code	Stoma Size	Color	HCPCS*
Cut-to-fit / Transparent			
21008	¾ - 2½" (10-55 mm)	Blue	A4410, A4427

NOTE: SenSura® Mio Click is NOT compatible with SenSura® Click or Assura® systems. SenSura® Click and Assura® are compatible.

SenSura® Click Flat Barrier

- Standard wear
- Belt tabs for use with 4215, 4220
- Not made with natural rubber latex



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10011	¾ - 1⅜" (10-35 mm)	Green	A4414
10021	¾ - 1¾" (10-45 mm)	Red	A4414
10031	¾ - 2⅛" (10-55 mm)	Blue	A4414
10041	¾ - 2½" (10-65 mm)	Yellow	A4415*

Pre-cut

10012	1⅜" (20 mm)	Green	A4414
10013	1" (25 mm)	Green	A4414
10022	1⅜" (30 mm)	Red	A4414
10023	1⅜" (35 mm)	Red	A4414
10032	1⅞" (40 mm)	Blue	A4414

SenSura® Click Xpro Flat Barrier

- Extended wear
- Belt tabs
- Not made with natural rubber latex



Units: 5

Code	Stoma Size	Color	HCPCS*
Cut-to-fit			
10015	¾ - 1⅜" (10-35 mm)	Green	A4409
10025	¾ - 1¾" (10-45 mm)	Red	A4409
10035	¾ - 2¼" (10-55 mm)	Blue	A4409
10045	¾ - 2½" (10-65 mm)	Yellow	A4410

Pre-cut

10016	1⅜" (20 mm)	Green	A4409
10017	1" (25 mm)	Green	A4409
10026	1⅜" (30 mm)	Red	A4409
10027	1⅜" (35 mm)	Red	A4409
10036	1⅞" (40 mm)	Blue	A4409