

SenSura® Mio Convex Light MAXI Drainable Pouch with Soft Outlet

- Extended wear
- Belt tabs for use with 4237, 4247
- **Special full-circle filter designed specifically for ileostomy output**
- Transparent pouches feature a printed scale for visual indication of output volume
- Opaque pouches feature a peek-a-boo inspection window and tuck-away option to keep the outlet from rubbing on the wearer's leg
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or any standard night bag
- **Not made with natural rubber latex**



11" (28 cm) / 630 mL				Units: 10
Code	Code	Cut	Stoma Size	HPCPS*
Transparent Opaque				
18676	18675	Cut-to-fit	$\frac{3}{8}$ - 1 $\frac{11}{16}$ " (10-43 mm)	A5057

SenSura® Flat MAXI Drainable Pouch

- Standard wear
- Soft outlet provides effective drainage; connects to a standard drain bag or Coloplast fistula drain bag #14010
- **Filter**
- Soft cloth on front and back
- Peek-a-boo inspection window allows easy application
- **Not made with natural rubber latex**



Opaque 12 $\frac{1}{2}$ " (32 cm) / 630 mL				Units: 30
Code	Cut	Stoma Size	HPCPS*	
19004	Cut-to-fit	$\frac{3}{8}$ - 3" (10-76 mm)	A4424	

For bed drainage bag, see page 72.

SenSura® Mio 1-piece Post Op Pouch with Window

- Neutral gray, water-resistant textile
- Printed scale provides visual indication of output volume
- Easy-to-open and close post-op window attached to pouch
- Soft outlet can connect to Coloplast's fistula drain bag #14010, or a drain bag with a standard connection
- **Not made with natural rubber latex**



Transparent				Units: 5
Code	Code	Cut	Volume	Stoma Size
Non-sterile Sterile				
18690	18680	Cut-to-fit	680mL	$\frac{3}{8}$ - 2 $\frac{3}{4}$ " (10-70 mm)
18691	18681	Cut-to-fit	725mL	$\frac{3}{8}$ - 4" (10-100 mm)

 Indicates convexity