

Introducing SenSura® Mio Kids

So kids can move like kids

SenSura® Mio Kids

SenSura Mio Kids
has BodyFit Technology®
to fit the contours of kids'
small, round bellies.



There's a Mio for everybody



The first SenSura® Mio pouching system designed specifically for kids

SenSura Mio Kids features the same elastic adhesive as SenSura Mio pouching systems for adults, enabling it to adapt to individual body movements. It's available in both a 1-piece and 2-piece solution.

Additional kid-friendly features




Star-shaped elastic barrier

The adhesive is star-shaped with special petals designed to hug kids' small, round bellies while giving more flexibility in terms of placement on the abdomen.



Fit zones

Special fit zones in the rim and center of the barrier increase its flexibility around the stoma for optimal range of movement.



Unique, full-circle pre-filter

This specially developed feature is proven to postpone filter clogging. It also reduces incidents of ballooning by 61%.*

*Not included in - or relevant for - urostomy pouches



Discreet material

SenSura Mio Kids features a neutral-colored, water-repellent, comfortable textile for a non-medical, discreet appearance.

Product assortment list

	Item code	Product name	Coupling size	Units / box	HCPCS*
Kids 1-piece 	18710	1-piece drainable pouch w/ inspection window, midi, neutral gray, 0 - 1 3/8"	–	10	A4424
	18711	1-piece drainable pouch, midi, transparent, 0 - 1 3/8"	–	10	A4424
	18712	1-piece urostomy pouch, transparent, 3/8 - 1 3/8"	–	10	A5071
	18713	1-piece urostomy pouch, neutral gray, 3/8 - 1 3/8"	–	10	A5071
Kids 2-piece pouches 	18720	2-piece Flex drainable pouch w/ inspection window, mini, neutral gray		10	A4425
	18721	2-piece Flex drainable pouch, mini, transparent		10	A4425
	18722	2-piece Flex drainable pouch w/ inspection window, midi, neutral gray		10	A4425
	18723	2-piece Flex drainable pouch, midi, transparent		10	A4425
	18724	2-piece Flex drainable pouch w/ inspection window, midi, neutral gray		10	A4425
	18725	2-piece Flex drainable pouch, midi, transparent		10	A4425
	18730	2-piece Flex urostomy pouch, neutral gray		10	A5073
	18731	2-piece Flex urostomy pouch, neutral gray		10	A5073
	18735	2-piece Flex barrier, cut-to-fit 0 - 1 1/16"		5	A4414
Kids 2-piece barriers 	18747	2-piece Flex barrier, cut-to-fit 0 - 1 3/8"		5	A4414

1. Coloplast, Clinical Study, 2010, 61% for ileostomates, 53% for colostomates.

*Reimbursement Disclaimer: Coloplast Corp. provides this information for your general reference and related to the reimbursement of Coloplast products only. Reimbursement, coverage and payment policies can vary from one insurer and region to another, and may change over time. Coloplast does not guarantee coverage or payment of products.

Ostomy Care / Continence Care / Wound & Skin Care / Interventional Urology

Coloplast Corp. Minneapolis, MN 55411 / 1-800-533-0464

www.coloplast.us The Coloplast logo is a registered trademark of Coloplast A/S. © 2019 Coloplast Corp. All rights reserved.



PM-05794 01.19