

Upholstered Drop Arm Wheeled Commode

HCPCS: E0165




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G




FIGURE H

- Simple, tool-free assembly
- Fully knock-down for easy storage and transport
- Drop-arms allow for easier lateral transfers
- Swing-away detachable footrests
- Collection pan included secures underneath the seat and can be easily removed from the backside of the unit
- Comfortable, wide backrest
- Durable and easy-to-clean padded vinyl seat cover

PRODUCTS	Description
11120SV-1F	Wheeled Commode with Footrests, 1/CV
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	27.99 lb
Seat Dimensions	16" (W) x 16" (D)
Weight Capacity	350 lb
Carton Shipping Dimensions	22.25" (L) x 22.45" (W) x 24.02" (H)
Carton Shipping Weight	30.87 lb
HCPCS	E0165
Width Between Arms	18.5" Outside Arms: 22.25"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]