



MQ5700



FILTER MQ0054 10/BAG



MQ5700FFC FILTER CAP



SPECIFICATIONS

Power	AC 120V, 60Hz / 230v,50Hz
Current	Max. 1.5A / 0.8A
Medication Capacity	5ml
Particle Size	0.5 to 10 µm
MMAD	4 µm
Sound Level	≤ 62 dBA (1 m away)
Maximum Medication capacity	13ml (cc)
Average Nebulization Rate	≥ 0.18ml/min
Compressor Pressure Range	17.4 to 29 Psi (0.12 to 0.2 Mpa)
Operating Pressure Range	7.25 to 10 Psi (0.05 to 10.07 Mpa)
Liter Flow Range	5 - 8 lpm
Operating Temperature Range	10°C to 40°C (50°F to 104°F)
Operating Humidity Range	10 to 95% RH
Storage Temperature Range	-20°C to 60°C (-4°F to 140°F)
Storage Humidity Range	10 to 95 % RH
Dimension (L x W x H)	170*136*90mm
Weight	1.18(kg)
Standard Accessories	nebulizer, air tube, filters (5pcs), mouthpiece

NEB KIT- REUSEABLE=18254

DISPOSABLE= NEB KIT 500



01/25/2023