

# COMPRESSION BANDAGES

A complete range of bandage products for the management of lymphedema. Designed with the patient and therapist in mind for optimal Complete Decongestive Therapy (CDT) outcomes.

## Multi-Layer Bandaging



### COMPRILAN®

- Short stretch compression bandage
- 100% cotton
  - Gold standard in bandaging for lymphedema management
  - Provides high working pressure and low resting pressure
  - Easy care and durable



### JOBST® COMPRIFOAM

- Open cell foam bandage
- 100% polyurethane
  - Air and moisture permeable
  - Reusable
  - May be steam sterilized



### TRICOFIX®

- Lightweight absorbent stockinette
- Used underneath Comprilan to provide a smooth layer
  - Breathable



### ARTIFLEX®

- Non-woven padding bandage
- Used underneath Comprilan to help protect sensitive skin
  - Provides padding around bony or irregular shapes
  - Ensures even distribution of pressure
  - Air permeable, non-absorbent



### ELASTOMULL®

- Soft comfortable bandage
- Provides light compression
  - Ideal for fingers and toes
  - Available in sterile or non-sterile



### JOBST® FOAM PAD

- Kidney shaped foam pad
- 100% synthetic rubber
  - Free of natural rubber
  - Evenly distributes pressure
  - Hand washable/air dry
  - Reusable

## 2-Layer Bandaging JOBST® COMPRI2

The JOBST® Compri2 is a 2-Layer Bandage System specifically designed for easy application. It provides graduated and reliable compression for up to seven days to help achieve effective edema reduction and helps heal the wound while skin stays dry and comfortable.



### Easy Application

- 1st layer of padding is customizable and easy to apply
- 2nd compression layer can be re-applied or adjusted as needed
- Indicator ensures correct compression level is achieved

### High Patient Comfort

- Padding can be added in sensitive areas
- Excellent mobility increases compliance and satisfaction
- Patients can wear with most footwear



### 1st Layer - Non Woven Padding

- **Stretchable** - simplifies application
- **Active Moisture Management** - keeps skin dry and comfortable
- **Customized Padding** - allows for anatomical adjustments and additional protection to sensitive areas



### 2nd Layer - Short Stretch Bandage

- **Cohesive** - delivers consistent and reliable compression
- **Adjustable** - allows modification during application
- **Visual** - indicator offers simple and effective application of correct compression