

SPONGES

Ducare® Woven Gauze Sponges

(Comparable to Curity® and J&J and Medline's Woven Gauze Sponges)

These woven sponges are made from 100% cotton. They are the gold standard for general use or for debriding, covering and packing wounds. These dressings come in a variety of sizes and are available in sterile and non-sterile. **Does not contain natural rubber latex.**



Product	Description	Pkg Unit/Case
STERILE		
82208	2" x 2", 8 ply	2 ea/pouch, 50 pouches/box, 30 boxes/case, 3000 ea/case
82312	3" x 3", 12 ply	2 ea/pouch, 25 pouches/box, 48 boxes/case, 2400 ea/case
82408	4" x 4", 8 ply	2 ea/pouch, 25 pouches/box, 24 boxes/case, 1200/case
82410	4" x 4", 12 ply	10 ea/pouch, 32 pouches/tray, 4 trays/case, 1280 ea/case
82412	4" x 4", 12 ply	2 ea/pouch, 25 pouches/box, 24 boxes/case, 1200 ea/case
82812	8" x 4", 12 ply	2 ea/pouch, 25 pouches/box, 20 boxes/case, 1000 ea/case
NON-STERILE		
90208	2" x 2", 8 ply	200/sleeve, 25 sleeves/case, 5000/case
90212	2" x 2", 12 ply	200/sleeve, 40 sleeves/case, 8000/case
90312	3" x 3", 12 ply	200/sleeve, 20 sleeves/case, 4000/case
90408	4" x 4", 8 ply	200/sleeve, 20 sleeves/case, 4000/case
90412	4" x 4", 12 ply	200/sleeve, 10 sleeves/case, 2000/case
90812	8" x 4", 12 ply (H Fold)	200/sleeve, 10 sleeves/case, 2000/case

Also offered in USP Type VII Gauze — See your Derma Sciences Representative for more information

Dupaque® X-Ray Detectable Woven Gauze Sponges

(Comparable to Vistec®, Ray-Tec®, and Medline's X-Ray Detectable Woven Gauze Sponges)



These woven sponges are made from 100% cotton and are low-linting, fluffy, and highly absorbent. They feature x-ray detectable barium sulfate strands for easy detection. These sponges are available as sterile and non-sterile. **Does not contain natural rubber latex.**

Product	Description	Pkg Unit/Case
STERILE		
9464S	4" x 4", 16 ply	10 ea/tray, 128 trays/case, 1280/case
9467S	8" x 4", 12 ply	10 ea/tray, 80 trays/case, 800/case

Also offered in USP Type VII Gauze

Dusoft® Non-Woven Sponges

(Comparable to Versalon®, NuGauze® and Avant Gauze®)



These non-woven sponges are ideal for wound cleansing and also can be used as a secondary dressing for added absorbency. They produce less lint and are less adherent than woven gauze sponges, so they cause less disruption to newly formed granulation tissue. The rayon/polyester blend provides fast wicking properties while being highly absorbent, minimizing skin maceration. **Does not contain natural rubber latex.**

Product	Description	Pkg Unit/Case
STERILE		
84122	2" x 2", 4 ply	2 ea/pouch, 25 pouches/box, 60 boxes/case, 3000 ea/case
84133	3" x 3", 4 ply	2 ea/pouch, 25 pouches/box, 48 boxes/case, 2400 ea/case
84144	4" x 4", 4 ply	2 ea/pouch, 25 pouches/box, 24 boxes/case, 1200 ea/case
84644	4" x 4", 6 ply	2 ea/pouch, 25 pouches/box, 12 boxes/case, 600 ea/case
NON-STERILE		
94122	2" x 2" 4 ply	200 ea/sleeve, 40 sleeves/case, 8000 ea/case
94133	3" x 3" 4 ply	200 ea/sleeve, 20 sleeves/case, 4000 ea/case
94144	4" x 4" 4 ply	200 ea/sleeve, 10 sleeves/case, 2000 ea/case