

Bandages/Rolls



Kerlix™ AMD Rolls

- Contain PHMB (Polyhexamethylene Biguanide).
- Resists bacterial colonization within the dressing.
- Reduces bacterial penetration through the dressing.
- Broad-spectrum effectiveness provides protection against gram-, gram+ and fungi/yeast microorganisms including MRSA and VRE.
- Limits cross-contamination from and to the patient, clinician and environment.
- May be used as a primary or secondary dressing
- Hypoallergenic.

Code	Description	Ship Case
Sterile in Rigid Plastic Tray 6 ply, 1 roll/tray		
3331	4.5" x 4.1 yds (11.4 cm x 411.5 cm)	60 rolls
Sterile in Soft Pouch 6 ply, 1 roll/Package		
3332	4.5" x 4.1 yds (11.4 cm x 411.5 cm)	100 rolls
Non-sterile Bulk		
1892AMD	4.5" x 4.1 yds (11.4 cm x 411.5 cm)	100 rolls

Kerlix™ Rolls

- Made of prewashed, fluff-dried 100% woven gauze with a unique crinkle-weave pattern for loft and bulk to cushion and protect wound areas.
- Six-ply construction with finished edges reduces loose ends and lint.
- As primary dressing, open-weave design provides fast wicking action, superior aeration and excellent absorbency.
- As secondary dressing, it provides bulk, cushioning and superior conformability.
- Ideal for bandaging heads, limbs and difficult-to-dress wounds; i.e., burns, plastic and orthopedic wounds.
- Choice of three sizes for all bandaging needs.
- Available in convenient tray, soft pouch, and a wide array of non-sterile put-ups.
- Code 3324 comes packaged in a convenient dispensing carton.



Code	Description	Per Box	Ship Case
Sterile in Soft Pouch 6 ply, 1 roll/package			
6715	4.5" x 4.1 yds		100 rolls
6725	3.4" x 3.6 yds		96 rolls
6720	2.25" x 3.0 yds	12	96 rolls
Sterile Soft Pouch 8 ply, 1 roll/package			
6716	4.5" x 3.1 yds		100 rolls
Sterile in Rigid Plastic Tray 6 ply, 1 roll/tray			
6730	4.5" x 4.1 yds		60 rolls
Non-sterile 6 ply, Bulk and Packed 12s			
1892	4.5" x 4.1 yds		100 rolls
3324	4.5" x 4.1 yds	12	48 rolls
6735	3.4" x 3.6 yds		96 rolls
1801	2.25" x 3.0 yds	12	96 rolls
Burn Pack, Sterile in Soft Pouch			
6760	4.5" x 4.1 yds	5 per pouch;16 pouches per box	80 rolls