

Latex Free

Alternating Pressure Pad and Pump System

- Bubble style pad effectively reduces pressure
- Automatic 5-minute cycle alternately inflates and deflates
- Heavy gauge, medical grade vinyl is easy to clean, flame retardant and anti-bacterial
- Pump mounts to bed with built-in mounting bracket and features an illuminated on/off switch
- Available with standard fixed pressure pump or deluxe variable pressure pump
- 78" x 34" x 2 1/2"
- 275 lb maximum weight capacity
- 90-day limited warranty on pad
- 1 year warranty on pump
- Latex Free


AQ1000	Standard Pump and Bubble-style Pad, (120V)	E0180	1 ea
AQ1000-2	Standard Pump only	E0182	1 ea
AQ1000-220V	Standard Pump and Bubble-style Pad, (220V)	E0180	1 ea
AQ2000	Deluxe Variable Pressure Pump and Bubble-style Pad	E0181	1 ea
AQ2000-2	Deluxe Variable Pressure Pump only	E0182	1 ea
AQ2000-1	Replacement Pad for AQ1000 and AQ2000		1 ea



AQ1000



AQ2000



	AQ1000	AQ2000
Power	120 Volt/60Hz, 1 Amp	120 Volt/60Hz, 1 Amp
Pressure	1.9 psi	1.0 – 2.0 psi
Air Output	4.5 liters/min	4.5 liters/min
Cycle Time	5 min	5 min
Power Cord Length	10 ft	10 ft
Pump Size	8 1/2" x 4 1/4" x 3 5/8"	9 3/4" x 5 1/8" x 3 7/8"
Pump Weight	2.9 lb	2.5 lb