

Powdercrush®

Automated Medication Crusher

An ideal product for long-term care, pharmacies, hospitals and correctional facilities



The Powdercrush saves time, eliminates the problem of cross contamination and crushes meds to a superfine powder. The Powdercrush is battery operated: powered by 2 long-lasting rechargeable batteries (charger included).

The Powdercrush is used with disposable medication pouches available through Crocus Medical. Stop hand crushing your meds and start powdercrushing with the powerful, portable Powdercrush from Crocus Medical.

The Powdercrush will help reduce the chance of repetitive strain injuries for nurses.

Medication Crushing as Easy as 1, 2, 3



1. Place pills to be crushed in a Crocus Medical supplied disposable medication pouch. Fold over the top of the pouch to keep the powder inside the pouch.
2. With one hand, insert the disposable medication pouch filled with pills in the slot on top of the machine. Press down on the start button to crush meds.
3. Typically, 4-5 seconds will produce the desired results. For a finer powder operate for a longer period.

Specifications

The Powdercrush system is used with disposable medication pouches available through Crocus Medical. The Powdercrush is fully rechargeable using lithium ion batteries.

Width: 8.25 cm (3.25 in)
Length: 27.9 cm (10.98 in)
Height: 9.5 cm (3.74 in)
Weight: 2.7 kg (5.95 lbs.)

Our Mission

Our mission at Crocus Medical is to help health care providers save lives by eliminating medication errors.

Our Commitment

Crocus Medical is dedicated to improving the safety of medication preparation through the development and supply of specialized error prevention products and equipment to the acute, long term, and community care markets. We specialize in medication compliance and are backed by over 35 years of industry experience.

[Click here to view detailed product information and pricing.](#)

Call 800-397-5899 to order today!



Crocus Medical Inc. is a US Distributor for Manrex Limited.