

Venturi masks

AirLife offers a wide selection of both standard and Dual Dial™ Venturi-style masks that feature either six diluter jets or two oxygen concentration dials. All masks are under-the-chin style and supplied with a six-inch flextube, humidification hood and seven-foot crush-resistant tubing. We also offer a Venturi-style with a tracheostomy mask for use with a 50 psi oxygen source and flowmeter.


Venturi mask		
Catalog no.	Description	Qty. (cs.)
001250	Pediatric, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ tubing	50
001240	Adult, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ tubing	50
001248	Adult, tracheostomy, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ tubing	50
001363	Adult, 6" flextube, 6 diluter jets (24%, 28%, 31%, 35%, 40%, 50%), humidification hood, 7' crush-resistant O ₂ U/Connect-It tubing	50
001259	Pediatric, 6" flextube, Dual Dial™, humidification hood, 7' crush-resistant O ₂ tubing	50
001255	Adult, 6" flextube, Dual Dial, humidification hood, 7' crush-resistant O ₂ tubing	50



Aerosol masks

These soft, high-quality AirLife disposable aerosol masks are made from high-grade resins. The smooth, feathered edges are designed for patient comfort by reducing irritation points. All AirLife masks are designed for anatomical compatibility. A comprehensive line is available to meet requirements in pediatric and adult sizes.

Aerosol mask		
Catalog no.	Description	Qty. (cs.)
001263	Pediatric	50
001261	Pediatric, short style	50
001206	Adult	50
001208	Adult, 6' corrugated tubing	50



Tracheostomy masks

The AirLife disposable tracheostomy masks are available in pediatric and adult sizes, and incorporate several patient features. The neckband is made from a comfortable, nonbiting material, and the swivel tubing connector permits access from either side of the patient. The special strap clips allow the mask to be removed with minimal disturbance to the patient.

Tracheostomy mask		
Catalog no.	Description	Qty. (cs.)
001225	Pediatric	50
001226	Adult	50
001227	Adult, large	50