

## Deluxe Folding Walker, Two Button with 5" Wheels

HCPCS: E0143



FIGURE A



FIGURE B




FIGURE C

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces (Figure A)
- Vinyl-contoured hand grip
- Comes with wheels (Replacement wheels item # 10109)
- The 10210-4ASM Walker comes fully assembled out of the box
- Piggy back label is supplied on product hangtag

### PRODUCTS

### Description

10210-1	Walker Adlt 2but w/5 whl.1/ea
10210-4ASM	Walker,2BtnAdult,5"WhlsAttached,4/cs
10211-4ASM	Dlx.Walker,2BtnJunior,5"WhlAttached,4/cs

### SPECIFICATIONS

### Description

Carton Shipping Weight	30   29   34   9.5 lbs.
Overall Width	24"
Warranty	Limited Lifetime
Weight	7.5   7   7.5 lbs.
Weight Capacity	350 lbs
Width	17" inside hand grip   20" inside back legs