

## Basic Full- and Semi-Electric Beds

### Light

- Easy to transport with a foot section that weighs just 42 lbs (19 kgs)

### Easy-To-Use

- Convenient DC motor box snaps off quickly, making delivery, setup, and routine maintenance easier

### Safe

- A 9V battery back up allows movement to lowest bed position during power failures (batteries not included)
- Headspring supports recessed to resist pinching and injuries
- Foot assembly recessed to reduce chance of pinching

### High Quality

- 5-year warranty on electronics and 2-year warranty on mechanical components
- 100% system tested and certified

Comfort and reliability at an economical price

## Basic Full-Electric Beds

Provides high reliability at an economical price



## Basic Semi-Electric Bed

Features both a smooth high-low motor and a traditional hand crank



### Features

- Light weight – foot section weighs just 42 lbs. (19 kg)
- Convenient DC motor box snaps off quickly, making delivery, setup, and routine maintenance much easier
- 9V battery backup allows lowering during power failures (battery not included)
- Recessed headspring supports and foot assembly help reduce the chance of pinching and injuries


### Specifications

Height range	15"-20" (38 cm-51 cm)
Bed deck dimensions	80" (203 cm) x 36" (91 cm)
Assembled weight	Full 155 lbs. (70 kg) Semi 147 lbs. (66.7 kg)
Weight capacity	450 lbs. (204 kg)
Casters	Four (2 locking)
Warranty	5-years for electrical and 2-years for mechanical components; lifetime for the bed frame



### Transport your Basic bed with ease. Just fold, attach and roll!

- Makes beds easy to transport and facilitates efficient storage



Item Number	Description
MDR107003E	Full Electric Bed
MDR107003ELO	Full Basic Electric Low Bed
MDR107002E	Semi-Electric Bed
MDSBEDCART	Transport Bed Cart

For Mattresses and Accessories see pages 108-110