

Med-Aire Plus 10" Bariatric Alternating Pressure and Low Air Loss Mattress Replacement System



INSET A

- Active combination bariatric therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal all stages of pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, auto-recovery, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent “bottoming out” and provide up to 24 hours of power outage protection
- 10" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable quilted Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS	Description
14030	Bariatric Alternating Pressure Mattress System, E0277, 1/ea
14048	Bariatric Alternating Pressure Mattress System, E0277, 1/ea
14054	Bariatric Alternating Pressure Mattress System (Replaces Mason # LS9000B-54), E0277, 1/ea
14060	Bariatric Alternating Pressure Mattress System (Replaces Mason # LS9000B-60), E0277, 1/ea

SPECIFICATIONS	Description
Control Unit Airflow	12 LPM
Control Unit Cycle Time	10, 15, 20, 25 Minutes
Control Unit Type	Compressor
Control Unit Alarms	Visual / Audible
Mattress Dimensions	42", 48", 54", 60" (W) X 80" (L) X 10" (H)
Mattress Cover Material	Fluid-Resistant, 4-way Stretch
Weight Capacity	600 lbs. - 1000 lbs.
Warranty	18 Month Limited

ACCESSORIES	Description
14427	Electric Pump (Sold Separately) (Replaces Mason # 8151-EP), 1/ea
14030XP	Replacement pump for 14030, 14048, 14054 & 14060, 1/ea
14030MN	Replacement Mattress for 14030, 1/ea
14048MN	Replacement Mattress for 14048, 1/ea
14054MN	Replacement Mattress for 14054, 1/ea

ACCESSORIES	Description
14060MN	Replacement Mattress for 14060, 1/ea
14030C	Replacement Cover, 1/ea
14048C	Replacement Cover, 1/ea
14054C	Replacement Cover, 1/ea
14060C	Replacement Cover, 1/ea
14048-CELL	Replacement Bladder, standard, 1/ea
14030-CELL	Replacement Bladder, standard, 1/ea
14051-CELL	Standard Cell, 1/ea
14054-CELL	Replacement Bladder, standard, 1/ea
14060-CELL	Replacement Bladder, standard, 1/ea
14030-LASERCELL	Replacement Ventilated Bladder, 1/ea
14048-LASERCELL	Replacement Ventilated Bladder, 1/ea
14054-LASERCELL	Replacement Ventilated Bladder, 1/ea
14060-LASERCELL	Replacement Ventilated Bladder, 1/ea