



Titan P22 Standard Powerchair

HCPCS: K0822



INSET A



INSET B



INSET C



INSET D

[Click here to view detailed product information and pricing.](#)
Call 800-397-5899 to order today!

- Height-adjustable swivel seat features fold-down backrest
- Includes interchangeable color panels in red and blue
- Seat back reclines and features adjustable headrest
- Adjustable-length controller mount (Figure B)
- Armrests are padded, adjustable and removable
- Available with 18" and 20" wide Pan Seats (Figure A)
- Optional swing-away controller arm (Item #AA4100) sold separately
- Solid, flip-up footplate with adjustable height, depth and angle
- Positioning belt included
- Optional elevating leg rests (requires Item #AE2500 & Item #LK3JELR) sold separately
- Programmable controller
- Flat-free, non-marking tires (Figure C)
- Easy free wheel operation (Figure D)
- Anti-tip wheels

PRODUCTS	Description
TITAN18P22	Titan P22, 1/cs
TITAN20P22	Titan P22 K0822 FWD 20 x 18 Pan Seat, 1/cs

SPECIFICATIONS	Description
HCPCS	K0822
Weight Capacity	300 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	25"
Climbing Angle	6°
Ground Clearance	2.5"
Dimensions	40" (L) x 24" (W)
Floor -to-Seat Height	20.5" - 23"
Controller	Dynamic Shark 40A
Motors (2)	24V x 200W x 3900 rpm
Batteries (pair)	12V x 35AH (U1)
Charger	3A Offboard
Brakes	Electromagnetic
Caster Wheels	6" x 2" Solid
Drive Wheels	10" x 3" Flat-free
Base Weight	68 lbs.
Battery Weight (pair)	44 lbs.

SPECIFICATIONS	Description
Seat Weight	43 lbs.
Total Weight	158 lbs.
Warranty on Frame	Limited Lifetime
Warranty on Electronic Controller	24 Months
Warranty on Drive Train Components	24 Months
Warranty on Batteries	6 months

ACCESSORIES	Description
AA4100	Swing-Away Joystick, 1/ea
AE2500	ELR Bracket, 1/ea
AB1000	Arm Bag, 1/ea
AZ0060	Cup Holder, 1/ea

PARTS	Description
TITAN-01	Footplate and Attaching Bracket, 1/ea
TITAN-14	Footrest Depth Adj. Bolt, 1/ea
TITAN-02	Drive Wheel Complete, 1/ea
TITAN-11	Drive Wheel Complete, 1/ea
TITAN-03	Caster with Bearings, 1/ea
TITAN-04	Caster Axle, 1/ea
TITAN-05	Frame, 1/ea
TITAN-06	Seat Post Bolt, 1/ea
TITAN-38	Seat Post, 1/ea
TITAN-39	Seat Post, 1/ea
TITAN-07	Battery Tray Complete, 1/ea
TITAN-08	Battery Hold Down Strap, 1/ea
TITAN-09	Motor and Gear Box Right, 1/ea
TITAN-10	Motor and Gear Box Left, 1/ea
TITAN-12	Battery Harness, 1/ea
LRA402307	Circuit Breaker, 1/ea
TITAN-13	Controller, 1/ea
TITAN-16	Shroud, 1/ea
TITAN-17	Decoration Set Red, 1/ea
TITAN-18	Decoration Set Blue, 1/ea
TITAN-29	Captain Seat Assy. 16 x 16, 1/ea
TITAN-31	Captain Seat Assy 18x16, 1/ea
TITAN-19	Captain Seat Assy 18x18x16.5, 1/ea
TITAN-20	Captain Seat Assy 18x18x18, 1/ea
TITAN-32	Captain Seat Assy 22x20, 1/ea
TITAN-33	Captain Seat Assy 20x20, 1/ea
STS01-181818	Pan Seat 18", 1/ea
ST502-202020	Pan Seat 20", 1/ea
LDM400001-07	Joy Stick Bracket Knob, 1/ea
TITAN-24	Joystick Bracket, 1/ea
LDM400001-05	Joystick Bracket sub, 1/ea
LPW182001	Joystick, 1/ea
DK-REMD01	Joystick Overlay, 1/ea
LPW182001-02	Joystick Cable, 1/ea
LPW182001-03	Joystick Extension Cable, 1/ea
LDM400001-04	Joystick Arm, 1/ea
LDM400001-01	Joystick Screw, 1/pr
LDM400013	Caster Fork, 1/ea
LDR4100006-01	Fork Stem Bearing, 1/ea
LDR410011L	Armrest Assembly Left, 1/ea
LDR410011R	Armrest Assembly Right, 1/ea