

# CHAIR SCALES

## MECHANICAL CHAIR SCALES



Detecto 475

Frequent weight checks of frail, handicapped, and elderly patients can be a tremendous, time-consuming task without the proper equipment. Detecto simplifies that job with a selection of innovatively-designed chair scales with wheels.

MODEL	CAPACITY
475	400 lb x 4 oz
4751	400 lb x 4 oz / 175 kg x 100 g
475-CFR	Covered foot rest

- Plastic-molded seat for patient comfort
- Lift-away arms and footrest for easy access
- Die-cast beam readable from both sides

## DIGITAL CHAIR SCALES



Detecto 6475



Developed to ensure patient safety and comfort as well as to reduce weighing time, these digital chair scales with wheels are conveniently portable and offer easy side or front access for your special patients. Features crisp waist-high LCD readouts.

MODEL	CAPACITY
6475	400 lb x 0.2 lb / 180 kg x 0.1 kg
6475K	180 kg x 0.1 kg
475-CFR	Covered foot rest

- Large 0.75 in / 19 mm high LCD digits
- Lift-away arms and footrest for easy access
- Fast readouts - accuracy to 0.2 lb or 0.1 kg

